ENZYWELL

TREPONEMA IgM

REF 91051 (96 tests)

Prodotto da/Manufactured by/Fabricado por:
DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX / INDICE / CONTEÚDO

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17. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S.A. Lukehart et al. Characterization of monoclonal antibodies to Treponema Pallidum. J. Immunology 134: 585 (1985).
3. S.A. Lukehart: Identification and characterization of Treponema Pallidum antigens by monoclonal antibodies. in: Monoclonal Antibodies, Academic Press 1986, p. 1.
4. L. Lewis et al.: Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
5. M. Norgard et al.: Sensitivity and specificity of monoclonal antibodies directed against antigenic determinants of Treponema Pallidum Nichols in the diagnosis of syphilis. J. Clin. Microbiol. 20: 711 (1984).
6. M.J. Bailey et al.: Production of murine monoclonal antibodies to the major axial filament polypeptide of Treponema pallidum. J. Gen. Microbiol. 133: 1805 (1987).

INSTRUCTIONS FOR USE

ENZYWELL TREPONEMA IgM

REF 91051

(English)

1. INTENDED USE

IMMUNOENZYMATIC CAPTURE METHOD FOR THE QUALITATIVE DETERMINATION OF IgM-CLASS ANTIBODIES TO TREPONEMA PALLIDUM IN HUMAN SERUM.

2. SUMMARY AND EXPLANATION OF THE TEST

The serological diagnosis of Syphilis is performed by demonstrating the presence of significant levels of specific antibodies for *Treponema pallidum* in the patient's serum.

The reference method used is the FTA-ABS technique which allows the detection of both IgG and IgM. However, its execution is laborious and the interpretation of the results is not simple; alternative methods have therefore been found to simplify the procedure. ELISA gives completely automated, reproducible and specific results, and for this reason is proposed by many operators. The specific assay of IgM as opposed to IgG is of great importance in the diagnosis of congenital Syphilis.

3. PRINCIPLE OF THE TEST

The test for the assay of *Treponema pallidum* IgM is based on the principle of the capture of these immunoglobulins and the subsequent identification of those which are specific, making use of their ability to bind an antigen conjugated to peroxidase. The capture is performed using monoclonal antibodies bound to the solid phase (microtiter wells). The antigen is composed of purified recombinant proteins of *Treponema Pallidum*, labelled with peroxidase (1-6).

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.
- **Bring to room temperature before use.**

MT PLATE MICROPLATE 12x8 wells coated with anti-human IgM monoclonal antibodies.

Use: open the package at the opposite end from the code (MM followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + POSITIVE CONTROL 1 x 1.6 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT-OFF CUT OFF CONTROL 1 x 2.5 mL

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONJ CONJUGATE 1x16 mL

Contents: Purified recombinant proteins of *Treponema pallidum* conjugated to peroxidase, in phosphate buffer with phenol 0.05% and Bronidox 0.02%.

Preparation: ready for use.

CONTROL IgM - IgM NEGATIVE CONTROL (PF93900) 1 x 1.6 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Diluted human serum in Phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF 10x WASH BUFFER 10X (PF93603) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL 10 DILUENT 10 (PF93621) 1 x 100 mL **INTERCHANGEABLE BETWEEN LOTS**

To be used to dilute samples.

Contents: PBS containing proteins 10% w/v and sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619) 1x12 mL Ready for use **INTERCHANGEABLE BETWEEN LOTS**

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602) 1x16 mL **INTERCHANGEABLE BETWEEN LOTS**

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED

- Incubator at 37°C
- Microplate reader, wavelength 450 or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µL solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Controls	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY.

Caution:

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin. Waste disposal: the used devices containing materials of human origin must be handled as infectious material and disposed of according to the full requirements of this issue.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.09% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides.

If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Solution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Avoid to use self-defrosting freezers for the storage of the samples.
8. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
9. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
10. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
11. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
12. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
13. Use of the kit with automatic instruments must be validated by the user.
14. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
 - installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C, and can be thawed a maximum of 3 times. Do not keep the samples in auto-defrosting freezers. Defrosted

samples must be carefully shaken before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results. Strongly lipemic, icteric or contaminated samples should be avoided.

The test cannot be performed on human plasma.

8. TEST PROCEDURE

Manual Technique

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x (100 mL + 900 mL H₂O).

Dilute samples 1:51 distributing 20 µL of serum into 1 mL of diluent. Dispense 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control. Leave one well for the blank, performed using 100 µL of the substrate mixture.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of the conjugate (purified recombinant proteins of *Treponema pallidum* conjugated to peroxidase) to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film.

The plate is washed again 4 times, as described above. Finally, the substrate is distributed, 100 µL/well.

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop Solution.

The adsorbance (O.D.) is read at 450 nm or at 450/620 nm within 30 min.

9. SCHEME OF TEST PROCEDURE

Manual Technique

STEP 1 Place 100 µL of diluted sample//controls in the wells of the strips.

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µL)

STEP 2 Add 100 µL of tracer to each well

-

Incubate for 45 min. at 37°C

-

Wash 4 times (300 µL)

STEP 3 Add 100 µL of Substrate to each well

-

Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop Solution

-

Read absorbance at 450 nm within 30 min

10. VALIDATION OF THE TEST

Subtract the value of the blank (≤ 0.150) from all the other readings. The O.D. values of the control Cut-off serum must be within 25% of the mean value. Disregard any abnormal value and recalculate the mean. The Positive control must have an O.D. at least 1.5 times that of the Cut-off serum. The ratio between Negative Control and Cut-off must be ≤ 0.6 .

The O.D. of the cut-off must be ≥ 0.2 at 450 nm and ≥ 0.16 at 450/620 nm.

11. INTERPRETATION OF THE RESULTS

Qualitative results

If the adsorbance of the sample is higher than that of the Cut-off, the sample is positive for the presence of specific IgM.

Calculate the ratio between the average O.D. value of the sample and that of the Cut-off.

The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE PROCEDURE

The results of this test should be interpreted in conjunction with other clinical and laboratory data.

As in the IgM response in other infectious diseases, a small number of Syphilis patients may present persistent low levels of IgM to *T. pallidum* even after the normal maximum period of 1-2 years.

In cases of suspected very early Syphilis in which the test may result negative, these patients should be retested about 1 week later.

13. ANALYTICAL SPECIFICITY

The following samples were tested containing potentially interfering substances:

- Highly hemolyzed samples (n = 3)
- Rheumatoid Factor (up to 1048 IU/ml) (n = 5)
- Bilirubin (up to 11 mg/dl) (n = 5)
- Triglycerides (up to 1281 mg/dl) (n = 4).

In no case was the result of the Enzywell Treponema IgM test influenced by the type of sample.

14. DIAGNOSTIC SENSITIVITY AND SPECIFICITY

In a clinical trial performed in a hospital laboratory, 275 samples were analysed. The samples were analysed with another commercial immunoenzymatic method. The results are summarized in the following table:

		Reference	
		+	-
Diesse	+	37	0
	-	3	235

The ENZYWELL kit offers a sensitivity of 92% and a specificity of 100% when compared with another commercial kit used as reference.

15. PRECISION

In-run Precision

	TPM1	TPM2	TPM3	Cut off	Positive Control
N =	24	24	24	12	12
O.D.	0.22	0.50	1.45	0.29	2.49
CV%	19	5	5	2	4

Between-run Precision

INDEX

Sample	Average	CV%
TPM3	4.1	11
TPM2	1.5	16
TPM1	0.7	16
Positive Control	7.5	8

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing

		the plate (see instructions for use point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessiccant must be pale yellow). Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S.A. Lukehart et al. Characterization of monoclonal antibodies to Treponema Pallidum. J. Immunology 134: 585 (1985).
3. S.A. Lukehart: Identification and characterization of Treponema Pallidum antigens by monoclonal antibodies. in: Monoclonal Antibodies, Academic Press 1986, p. 1.
4. L. Lewis et al.: Evaluation of immunoglobulin M Western blot analysis in the diagnosis of congenital syphilis. J. Clin. Microbiol. 28: 296 (1990).
5. M. Norgard et al.: Sensitivity and specificity of monoclonal antibodies directed against antigenic determinants of Treponema Pallidum Nichols in the diagnosis of syphilis. J. Clin. Microbiol. 20: 711 (1984).
6. M.J. Bailey et al.: Production of murine monoclonal antibodies to the major axial filament polypeptide of Treponema pallidum. J. Gen. Microbiol. 133: 1805 (1987).

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