MASTAFLUOR™ SKMA

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Introduction

Auto-immune diseases are caused by a disorder of the cellular and/or humoral immunological reaction. These reactions which normally occur against external influences may under certain circumstances turn against the body itself and thereby cause various diseases.

Anti-skeletal muscle antibodies are directed against a variety of cytoplasmic components of the cross stripes of skeletal and heart muscle, also of thymic cells. They react predominantly with the I- but also with the A-band and the Z-stripes of the sarcomere and with the sarcoplasmic reticulum. The biochemical antigens are myosin, heavy-meromyosin, actin and others.

Anti-skeletal muscle antibodies with titers higher than 1:60 mainly occur with myasthenia combined with thymom, more rarely with polyomyositis and other muscular diseases. High titers indicate a thymom, a server form or an exacerbation. Permanent high titers may be marker of a poor prognosis, while the decrease of titer indicates an improvement.

These antibodies are observed more often in younger patients, at the beginning of the illness and in carriers of certain HLA antigens. Following thymectomy titres may not necessarily decrease.

Intended Use

The determination is performed using tissue sections of rat thigh muscle. A sample is considered positive if a green fluorescence of specific cell structures of the tissue sections is visible.

Test Principle

IFA slides are coated with combined sections of rat thigh muscle. Diluted patient serum or plasma is incubated on the wells of the slide. In the presence of specific antibodies a stable antigen-antibody complex is formed. Unspecific or unbound antibodies are removed in a washing step. Specific immune complexes are then detected by FITC-conjugated anti-human IgG antibodies. Again unspecific or non-bound antibodies are removed in a washing step. Results could be read visually using an immunofluorescence microscope. Positive reactions show specific fluorescent reactions with rat thigh muscle tissue.

Kit Contents

Content	100 Tests	50 Tests	Specifications
Slides	10 x 10 Wells	10 x 5 Wells	Coated with tissue section of rat thigh muscle
Positive control	0.5 mL	0.5 mL	Positive serum, ready to use, contains < 0.1 % NaN ₃
Negative control	0.5 mL	0.5 mL	Negative serum, ready to use, contains < 0.1 % NaN ₃
FITC-conjugate	3 mL	2 mL	FITC-conjugated anti-human IgG, ready to use, contains < 0.1 % NaN ₃
Evans Blue	3 mL	3 mL	Ready to use, contains < 0.1 % NaN ₃
Mounting Medium	3 mL	3 mL	Ready to use, contains < 0.1 % NaN ₃
PBS	2 Sachets	2 Sachets	1 sachet to be diluted in 1 L of distilled water, pH 7.4 ± 0.2

Materials required but not provided

1. Sterile test tubes
2. Micropipettes and tips
3. Staining dish or Coplin jar
4. Moist chamber
5. Volumetric flask for PBS
6. Distilled water or water of higher quality
7. Forceps
8. Transmitted or epi-fluorescent microscope with a 490 nm excitation filter and a 510 nm barrier filter
9. Cover slips
10. Wash bottle

Stability and Storage

MASTAFLUOR™ SKMA is stable until end of shelf life as indicated on kit label. The kit should be stored at 2–8 °C.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Serum and plasma specimens may be stored at 2–8 °C for up to 48 hours prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Warnings and Precautions

1. The reagents supplied in this kit are for *in vitro* diagnostic use only.
2. Read instructions carefully before conducting the assay. Do not modify the procedure.
3. Do not use kit beyond the expiry date.
4. Wear appropriate protective clothing and use appropriate facilities.
5. Do not mouth pipette.
6. Use disposable plastic ware where ever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.
7. Sera provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting disease. No guarantee is given that the sera are free of infection or microbial contamination.
8. Use distilled water or water of higher quality.
9. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
10. Do not cross-contaminate reagents or interchange caps on bottles. Use a separate pipette or pipette tips for each sample and reagent.
11. Do not allow wells to dry out during the assay procedures.
12. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
13. Contaminated plastic ware should be disposed of and incinerated. Contaminated glassware should either be autoclaved at 121 °C for 1 hour or decontaminated with a solution of 2.5 % sodium hypochlorite. All liquid waste should be decontaminated appropriately e.g. using sodium hypochlorite.
14. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.
15. Microbial contaminated serum samples should not be used.
16. Sodium azide is used as a preservative as marked. It may be toxic if ingested. Sodium azide may react with lead and copper plumbing to form highly explosive salts. Always dispose of by flushing to drain with plenty of water.

Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or water of higher quality.
3. Dilute patient samples 1:20 with PBS.
Example: 10 µL serum + 190 µL of PBS.
4. Controls are supplied ready to use and must not be diluted.
5. Carefully remove the required number of slides from their sachets and mark accordingly. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
6. Apply 20–25 µL of diluted specimens and ready to use controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered and that serum does not escape from the wells.

Direct contact of a pipette tip with the slide surface may result in damage to the antigen substrate and should be avoided.
7. Cover the moist chamber and incubate at room temperature for 30 min.
8. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the centre of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
9. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min. A change of the PBS buffer after 5 min increases washing strength.
10. Remove slides one at a time from the staining jar and drain off any excess PBS. Using a blotter, dry the area outside wells. Wipe off any excess PBS from the back of the slide. Do not touch the well surface.
11. Immediately transfer the slide to a moist chamber and add 20–25 µL of FITC-conjugate to each well. Ensure that all the wells are covered by the conjugate.
12. Incubate slides at room temperature in a covered moist chamber in the dark for 30 min.
13. Switch on microscope 15 min prior to use.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Optional counter stain: add 3–5 drops of Evans Blue to the PBS wash buffer in the last 5 min washing step. The counter stain increases the contrast of the immunofluorescent pattern.

Briefly wash the slides in PBS to remove all Evans Blue.
16. Remove the slides from the PBS after washing and blot dry as in step 10.
17. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.

Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering. Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip.

18. Examine under a fluorescence microscope at a total magnification of 200x to 400x.

Interpretation and diagnostic relevance

The evaluation should always be performed with the positive and negative controls.

The presence of anti-skeletal muscle antibodies is indicated by fine and continuous fluorescent cross lines vertical to the myofibrils.

This is caused by the binding of specific antibodies to the various regions of the sarcomere mainly in the regions of the I- but also of the A-band and Z-stripes.

To distinguish between the fluorescence of the I- and the A-band, a phase-contrast picture is advised.

Limitations / Cross reactions

Cross-reactivity

Cross-reactivities are unknown.

Precision and reproducibility

With this immunofluorescence assay, no difference in the Interlot and Interassay variability by using the controls could be detected.

Internal quality control for user

All controls should be carried with every test run. The desired values of the positive and negative controls should be taken from the Lot specific data given on the Quality Control Certificate. If controls give invalid levels then results from test samples are invalid, requiring retesting. For every new lot, an internal quality control should be performed by testing one or more internal control samples which should be clearly specified.

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