ENZYWELL

HERPES SIMPLEX 1 IgG

REF 91027 (96 tests)

Manufactured by: DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) - Italy

INDICE / INDEX

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17. BIBLIOGRAFIA

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) Italy
Tel. 0577-587111

INSTRUCTIONS FOR USE

ENZYWELL HERPES SIMPLEX 1 IgG

REF 91027

(English)

1. INTENDED USE

IMMUNOENZYMATIC METHOD FOR THE QUALITATIVE DETERMINATION OF IgG-CLASS ANTIBODIES TO HERPES SIMPLEX VIRUS (TYPE 1) IN HUMAN SERUM, AS AN AID IN THE DIAGNOSIS OF HERPES SIMPLEX VIRUS INFECTION.

2. SUMMARY AND EXPLANATION OF THE TEST

The Herpes simplex virus (HSV) is a member of the Herpesviridae family, of which two types are known: type 1 (HSV-1) and type 2 (HSV-2) which present slight antigenic differences. HSV-1 causes chiefly oral-facial lesions, while HSV-2 is mainly responsible for genital lesions, but this distinction is not binding, both types occasionally causing infection in either anatomical site. HSV may also cause a form of ocular cheratitis, and lesions of the central nervous system.

HSV can affect practically the whole population. The primary infection is often in a subclinical form and is rarely diagnosed. After a latency period of variable duration, reactivation may occur and viral replication may or may not give rise to clinical lesions. Infection contracted during birth is of particular interest, this being an important cause of morbidity and mortality. It is therefore important to determine the immunitary state of women during pregnancy in order to detect serum conversion. The assay of specific IgG is important for the diagnosis of the serological state of the patient.

3. PRINCIPLE OF THE TEST

The test is based on the ELISA technique (Enzyme linked Immunosorbent Assay).

The antigen, composed of purified, inactivated Herpes Virus type 1, is bound to the solid phase (strips of 1x8 wells). Through incubation with human serum diluted in a solution containing HSV-2 (Sorbent 2), the specific immunoglobulins are bound to the antigen.

After washings to eliminate the proteins which have not reacted, incubation is performed with the conjugate composed of anti-IgG monoclonal antibodies labelled with peroxidase. The unbound conjugate is eliminated and the peroxidase substrate added. The colour which develops is proportional to the concentration of specific antibodies present in the serum being tested.

4. KIT CONTENTS AND REAGENT PREPARATION

- Reagents are sufficient for 96 determinations.

- Bring to room temperature before use.

MT PLATE MICROPLATE. 12x8 wells coated with herpes simplex virus type 1.

Use: open the package at the opposite end from the code (H followed by the lot number) which is useful for identification purposes, remove the support and strips to be used from the foil package, and place the unused strips in the polythene bag with the silica gel, expel the air and seal by pressing the closure.

CONTROL + POSITIVE CONTROL (1 x 1.6 mL)

Contents: Diluted human serum, at known concentration of anti-HSV 1 IgG, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONTROL CUT OFF CUT-OFF CONTROL (1 x 2.0 mL)

Contents: Diluted human serum, at known concentration of anti-HSV 1 IgG, in Phosphate buffer 0.01 mol/L containing BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

Colour: the colour is proportional to the relative antibody titer.

CONJ CONJUGATE. 1x16 mL.

Contents: monoclonal antibodies labelled with Peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%. Ready for use without further dilution.

SORBENT 2. 6 vials of freeze-dried powder.

Contents: solution containing HSV-2.

Preparation: reconstitute with the volume of diluent reported on the label, mixing by inversion.

CONTROL - IgG NEGATIVE CONTROL (PF93910) 1 x 1.6 mL INTERCHANGEABLE BETWEEN LOTS

Contents: Human serum in Phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%, liquid, ready for use without further dilution.

WASH BUF WASH BUFFER 10X (PF93603). 1 x 100 mL. INTERCHANGEABLE BETWEEN LOTS

Contents: Phosphate buffered saline, concentrated 10 times; contains Brij 0.5% .

Preparation: dilute the required volume 1:10 with distilled water in order to obtain the washing buffer ready for use. If crystals are present, they should be dissolved at 37°C before dilution.

SAMP DIL DILUENT 2 (PF93611). 1 x 100 mL. INTERCHANGEABLE BETWEEN LOTS

To be used to dilute samples.

Contents: Proteic solution in phosphate buffer with sodium azide 0,09% containing methyl orange as dye.

SUBS TMB SUBSTRATE (PF93619). 12 mL. Ready for use. INTERCHANGEABLE BETWEEN LOTS

Contents: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilised in citrate buffer 0.05 mol/L (pH 3.8).

H₂SO₄ 0.3 M STOP SOLUTION (PF93602). 1x16 mL. INTERCHANGEABLE BETWEEN LOTS

H₂SO₄ 0.3 mol/L, in solution ready for use.

ADHESIVE FILMS (2).

POLYTHENE BAG (1).

MATERIALS REQUIRED BUT NOT PROVIDED.

- Incubator at 37°C
- Microplate reader, wavelength 450 or 450/620 nm, with OD linearity up to 2,000 (at least).
- Microplate washer (preferable) able to dispense volumes in the range 225-375 µL
- Distilled or deionized water
- Normal laboratory glassware: cylinders, test-tubes etc.
- Micropipettes for the accurate collection of 10, 100, 1000 µl solution
- Disposable gloves
- Timer
- Sodium Hypochlorite solution (5%)
- Containers for collection of potentially infectious materials
- Absorbent tissue.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
Microplate	5 weeks at 2/8°C, polythene bag
Control sera	5 weeks at 2/8°C
Conjugate	5 weeks at 2/8°C
Sorbent 2	the solution cannot be stored and must be used after reconstitution
Substrate	up to the expiry date at 2/8°C, 1 week at 15-30°C; store in the dark
Sample Diluent	up to the expiry date at 2/8°C
Wash Buffer	2 weeks at 2/8°C, 5 days at 15/30°C.
Stop Solution	up to the expiry date at 2/8°C.

6. PRECAUTIONS

For in vitro diagnostic use only.

This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Health and Safety Information

1. Do not pipette by mouth. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. The following reagents contain low concentrations of harmful or irritant substances:
 - a) The Wash Buffer contains detergents
 - b) The conjugate contains phenol
 - c) The substrate is acid
 - d) The controls contain 0.9% Sodium Azide which can react with lead and copper in plumbing forming highly explosive deposits of metal azides; dilute with large amounts of water to eliminate.
 If any of the reagents come into contact with the skin or eyes, wash the area extensively with water.
3. Non-disposable apparatus should be sterilized after use. The preferred method is to autoclave for 1 h at 121°C; disposables should be autoclaved or incinerated.
4. Sulphuric acid required for the Stop Soution and hydrochloric acid used for washing glassware are corrosive and should be handled with appropriate care. If they come into contact with the skin or eyes, wash thoroughly with water.
5. Neutralized acids and other liquid waste should be decontaminated by adding a sufficient volume of sodium hypochlorite to obtain a final concentration of at least 1.0%. A 30 minute exposure to 1% sodium hypochlorite may be necessary to ensure effective decontamination.
6. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature. **It is important to work at the correct temperature. Check that the thermostat does not go below 35°C or over 39°C.** Open the envelope containing the strips after at least ½ hr at room temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
6. Do not allow wells to become dry during the assay procedure.
7. Care must be taken not to cross-contaminate reagents. It is important that pipettes are dedicated for exclusive use with the various reagents.
8. Care should be taken to avoid touching or splashing the rim of the well with conjugate. Do not "blow-out" from microplates.
9. Enzyme immunoassays can occasionally exhibit an "edge effect" which must be minimized by increasing the humidity during incubation steps. Plates must be covered with their covers and incubated at 37°C either in a water bath with a rack or float to support the plates if necessary, or in an incubator. Alternatively, plates can be incubated in an approved analyzer. See the appropriate operating manual for further details. CO₂ incubators must not be used.
10. Ensure that the bottom of the plate is clean and dry, and that no bubbles are present on the surface of the liquid before reading the plate.
11. Use of highly hemolyzed samples, incompletely clotted sera, or samples with microbial contamination may give rise to erroneous results.
12. Use of the kit with automatic instruments must be validated by the user.

13. For each instrument used, read the manufacturer's instructions manual carefully to obtain additional information on the following points:
- installation and particular requisites
 - operating principles, instructions, precautions and risks
 - manufacturer's specifications and instrument performance
 - servicing and maintenance.

7. TYPE AND STORAGE OF SAMPLE

The sample is composed of serum collected in the normal manner from the vein and handled with all precautions dictated by good laboratory practice. The fresh serum may be stored for 4 days at 2/8°C, or frozen for longer periods at -20°C; avoid repeated freezing and thawing. Defrosted samples must be carefully shaken before performing the test. Heat inactivation can lead to erroneous results. The quality of the sample can be seriously affected by microbial contamination which leads to erroneous results.

Strongly lipemic, icteric or contaminated samples should be avoided. If a new sample cannot be obtained, such samples should be clarified by filtration (0.45 µm) or centrifugation (3000 rpm x 10').

The test cannot be performed on human plasma.

8. TEST PROCEDURE

- Prepare the required number of strips.
- Prepare the washing buffer by diluting the Wash Buffer 10x.
- Prepare Sorbent 2 by reconstituting with diluent (volume reported on the label).

Dilute samples 1:101 distributing 10 µL of serum into 1 mL of diluent. Leave one well for the blank, performed using 100 µL of the substrate mixture.

Dispense 30 µL of Sorbent 2 and 100 µL of each diluted sample per well (duplicate testing is recommended). Place UNDILUTED controls in a strip (100 µL in each well **without the Sorbent 2**). The minimum requisite is 1 negative control, 2 cut-off and 1 positive control.

Wells are covered with protective film and incubated for 45 minutes at 37°C. After washing four times for 30 seconds (300 µL), add 100 µL of conjugate to each well and incubate again for 45 minutes at 37°C, covering the wells with the protective film. The plate is washed again 4 times, as described above. Finally, the substrate is distributed (100 µL/well).

After 15 minutes at room temperature the enzymatic reaction is stopped with 100 µL of Stop solution.

The adsorbance (O.D.) is read at 450 nm or 450/620 nm within 30 min.

9. TEST PROCEDURE FOR HERPES 1 IgG

STEP 1 Place 100 µL of controls (positive, negative and cut off) in the wells of the strip (without the SORBENT).

Place 30 µL of Sorbent 2 and 100 µL of diluted Serum in the wells of the strip

Incubate for 45 min. at 37°C

Wash 4 times (300 µL)

STEP 2 Add 100 µL of conjugate to each well

Incubate for 45 min. at 37°C

Wash 4 times (300 µL)

STEP 3 Add 100 µL of Substrate to each well

Incubate for 15 min. at R.T.

STEP 4 Add 100 µL of Stop solution

Read absorbance at 450 nm within 30 min

10. VALIDATION OF THE TEST

Subtract the value of the blank (≤ 0.150) from all the other values. The value in OD of the Cut-off control must be within 25% of the average value if tested in triplicate. Discard any abnormal value and recalculate the average. The positive control must have an OD at least 1.5 times that of the Cut-Off. The ratio between Negative and Cut-Off must be ≤ 0.6 . The O.D. of the cut-off must be ≥ 0.2 at 450 nm and ≥ 0.16 at 450/620 nm.

11. QUALITATIVE RESULTS

If the adsorbance of the sample is higher than that of the Cut-off, the sample is positive for the presence of specific IgG. Calculate the ratio (Index) between the O.D. value of the sample and that of the Cut-off. The sample is considered:

Positive: if the ratio is > 1.2 .

Doubtful: $\pm 20\%$ of the Cut-off.

Negative: if the ratio is < 0.8 .

If the result is doubtful, repeat the test. If it remains doubtful, collect a new serum sample.

12. LIMITATIONS OF THE TEST PROCEDURE

A serum sample obtained during the acute phase of infection, when only IgM antibodies are present, may be negative by this procedure.

The Herpes Simplex Virus IgM level should be determined using Enzywell Herpes Simplex Virus IgM kit. Alternatively, a second serum sample obtained 8-14 days later, should be tested in parallel to determine an increase in the IgG antibody level.

The test result should be used in conjunction with information available from the evaluation of other diagnostic procedures.

13. ANALYTICAL SPECIFICITY

17 samples which were non reactive for Herpes Simplex and contained IgG antibodies against other viruses such as Rubella, Epstein Barr, Cytomegalovirus, Mumps, were tested. In no case did the presence of these antibodies interfere in the test.

14. PERFORMANCE CHARACTERISTICS

In a clinical trial performed in a hospital laboratory, 174 samples were analysed, in comparison with another commercial immunoenzymatic method. The results can be summarized as follows:

Other method		
	P	N
Diesse P	158	2
N	1	13

Diagnostic sensitivity: 99.4% (IC at 95% 96.5 - 99.9)

Diagnostic Specificity: 86.7% (IC at 95% 62.1 - 96.2%)

15. PRECISION

Table 1 "In run" Precision

Sample	Average O.D.	CV%
1 (n=24)	0.107	8
2 (n=24)	0.396	3
3 (n=24)	0.977	4
Cut Off (n=16)	0.375	8
C. Pos. (n=16)	1.838	6

Table 2 "Between run" Precision (n=3)

Index	1	2	3	average	S.D.	CV%
C Pos	3,70	4,05	4,90	4,2	0,6	15
1	0,15	0,14	0,18	0,2	0,02	13
2	1,13	1,31	1,41	1,3	0,14	11
3	3,30	4,36	3,54	3,7	0,56	15

16. TROUBLE SHOOTING GUIDE

PROBLEM	POSSIBLE SOURCE	TEST OR ACTION
Invalid run (all negative)	One or more reagents not added or added in wrong sequence	Recheck procedure Check for unused solutions. Repeat test.
	Unreactive plate	Check the code on the package containing the plate (see package insert point 4 for correct code).
		Check for moisture in unused plate. (Silica gel dessicant must be pale yellow).Repeat test
Invalid run (all positive)	Contamination of substrate	Take new aliquot of substrate.
	Inadequate washing	Ensure that wash apparatus works well
Poor precision	Incomplete washing of wells	Ensure that wash apparatus works well
	Inadequate aspiration of wells	Ensure that wash apparatus works well
	Pipetting error	Check pipette function
	Reagent addition too slow	Avoid drying of the plate after washing step. Add reagents immediately
	Presence of bubbles	Avoid air bubbles during pipetting.
	Optical pathway not clean	Check instrument light source and detector for dirt. Wipe bottom of plate with soft tissue.
Inadequate Color development	Incorrect incubation times or temperature	Check for temperature control and time monitoring
		Adhere to recommended instruction for use.
	Inadequate volume of substrate added to the plate	Check pipette function.

17. REFERENCES

1. G.B. Wisdom: Enzyme-Immunoassay. Clin. Chem. 22: 1243 (1976).
2. S. Land et al.: Rapid diagnosis of herpes simplex virus infections by enzyme-linked immunosorbent assay. J. Clin. Microbiol. 19: 865 (1984).
3. B. Gonik et al.: Comparison of two enzyme-linked immunosorbent assays for detection of herpes simplex virus antigen. J. Clin. Microbiol. 29: 436 (1991).
4. C. Gleaves et al.: Evaluation of an enzyme immunoassay for the detection of herpes simplex virus (HSV) antigen from clinical specimens in viral transport media. J. Virological Meth. 28: 133 (1990).
5. M. Morgan and T. Smith: Evaluation of an enzyme-linked immunosorbent assay for the detection of herpes simplex virus antigen. J. Clin. Microbiol. 19: 730 (1984).
6. D. Ho et al.: Indirect ELISA for the detection of HSV-2 specific IgG and IgM antibodies with glycoprotein G (gG-2). J. Virological Meth. 36: 249 (1992).
7. R. Eberle et al.: The immune response to herpes simplex virus: comparison of the specificity and relative titers of serum antibodies directed against viral polypeptides following primary herpes simplex virus type 1 infections. J. Med. Virology 16: 1247 (1985).
8. J.E. Kuhn et al.: Analysis of the IgM and IgG antibody response against herpes simplex virus type 1 (HSV-1) structural and nonstructural proteins. J. Medical Virology 23: 135 (1987).

DIESSE Diagnostica Senese
Via delle Rose 10
53035 Monteriggioni (Siena) Italy
Tel. 0577-587111