DIESSE Diagnostica Senese SpA

Via delle Rose 10, 53035 Monteriggioni SI
 Tel. ++39/0577/ 587111 - Fax ++39/0577/318690

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SYPHILIS FAST (TPHA MODIFICATO)

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(English)

1. INTENDED USE

LATEX AGGLUTINATION TEST FOR THE RAPID ASSAY OF TOTAL ANTI-TREPONEMA PALLIDUM ANTIBODIES IN HUMAN SERUM USING THE SLIDE METHOD

2. SUMMARY AND EXPLANATION OF THE TEST

The serological diagnosis of syphilis is performed by revealing the presence of significant levels of specific Treponema pallidum (TP) antibodies in the serum. The methods at present most widely used are: VDRL and TPHA, ELISA IgG and IgM, total Ig with "competitive" methods.

VDRL is used to assay non-specific antibodies; it detects acute-phase positive reactions, but poses problems of false-positive results and poor sensitivity. The hemoagglutination test (TPHA) assays antibodies specific for Treponema pallidum, but the use of erythrocytes of animal origin leads to non specific reactions (hemoagglutinating power) with some human sera.

The ELISA method is the most advanced technique in use at present, thanks to the possibility of its complete automation, its specificity and sensitivity.

3. PRINCIPLE OF THE TEST

The system is composed of polystyrene latex particles coated with three of the immunodominant proteins of Treponema pallidum (17K, 15K and 42K) obtained in E. coli by recombinant technology. When human serum containing anti-Treponema pallidum antibodies is placed in contact with the sensitised latex, ***agglutination takes place within 8 minutes, and is easily visible***; in the absence of antibodies (negative serum) the latex remains uniform, and there is no trace of agglutination.

4. KIT CONTENTS

Bring reagents to room temperature before use.

LATEX ADSORBED LATEX (R1) Freeze-dried, 1 x 1.1 mL

Contents: Blue-coloured polystyrene latex particles, coated with specific stabilised and freeze-dried Treponema pallidum proteins.

Preparation: To be reconstituted with 1.1 mL of distilled water. Bring the latex to room temperature before use.

CONTROL + POSITIVE CONTROL SERUM 1 x 0.5 mL

Contents: Positive human serum. The serum is calibrated to be considered highly positive.

CONTROL - NEGATIVE CONTROL SERUM 1 x 0.5 mL

Contents: Serum of animal origin, totally non-reactive with the adsorbed latex.

OTHER MATERIALS

- 10 Slides to perform up to 11 reactions simultaneously
- 100 Mixing sticks. Use a new stick for each test and throw away after use.

OTHER MATERIAL REQUIRED, BUT NOT SUPPLIED

- Distilled water.
- Micropipettes for distribution of latex and sera
- Physiological saline
- Source of light.

5. STORAGE AND STABILITY OF REAGENTS

Reagents must be stored at 2/8°C.

The expiry date is printed on each component and on the box label.

Reagents have a limited stability after opening and/or preparation

REAGENT	CONDITIONS
LATEX	2 months after reconstitution
POSITIVE CONTROL	2 months.
NEGATIVE CONTROL	2 months

6. PRECAUTIONS

FOR IN VITRO DIAGNOSTIC USE ONLY. STORE AT 2-8°C

Caution: This kit contains materials of human origin which have been tested and gave a negative response by FDA-approved methods for the presence of HbsAg and for anti-HIV-1, anti-HIV-2 and anti-HCV antibodies. As no diagnostic test can offer a complete guarantee regarding the absence of infective agents, all material of human origin must be handled as potentially infectious. All precautions normally adopted in laboratory practice should be followed when handling material of human origin.

Waste disposal: serum samples and reagents once used must be treated as infectious residuals and eliminated according to law.

Health and Safety Information

1. Wear disposable gloves and eye protection while handling specimens and performing the assay. Wash hands thoroughly when finished.
2. Spillage of potentially infectious materials should be removed immediately with adsorbent paper tissue and the contaminated area swabbed with, for example, 1.0% sodium hypochlorite before work is continued. Sodium hypochlorite should not be used on acid-containing spills unless the spill area is first wiped dry. Materials used to clean spills, including gloves, should be disposed of as potentially biohazardous waste. Do not autoclave materials containing sodium hypochlorite.
3. Control sera contain Sodium Azide 0.09% as preservative which may react with copper and lead in plumbing to form potentially explosive metal azides.

Analytical precautions

1. Allow all reagents and samples to come to room temperature (18-30°C) before use. Immediately after use return reagents to the recommended storage temperature.
2. Do not use the reagents beyond the stated expiry date. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
3. Do not modify the Test Procedure or substitute reagents from other manufacturers or other lots unless the reagent is stipulated as interchangeable. Do not reduce any of the recommended incubation times.
4. Any glassware to be used with the reagents should be thoroughly washed with 2M hydrochloric acid and then rinsed with distilled water or high quality deionized water.
5. Avoid the use of self-defrosting freezers for the storage of samples.
6. Do not expose reagents to strong light or hypochlorite fumes during storage or during incubation steps.
7. Use of highly hemolyzed samples, incompletely clotted sera, plasma samples containing fibrin or samples with microbial contamination may give rise to erroneous results.
8. Avoid contaminating the slide with the dust from disposable gloves.

7. SAMPLES AND STORAGE

Serum obtained from blood collected without any particular measures.

No special preparation of the patient is necessary. The use of anticoagulants such as citrate, heparin or EDTA does not interfere in the test. Samples can be stored for 7 days at 2-8°C. For longer storage, freeze at -20°C. Freezing/thawing of the samples or heat inactivation for 30 min at 56°C does not influence the results. Do not keep the samples in auto-defrosting freezers.

Strongly lipemic, icteric or contaminated samples should be avoided.

8. PROCEDURE

*** MANUAL TEST PROCEDURE**

Pipet 40 µL of whole human serum and 20 µL of coated latex (R1) into a division on the slide; **mix the latex with the sample over the whole surface of the circle on the slide** using a stick, and rotate for 8 minutes.

Examine under a lamp for the presence of agglutinate.

Include the positive and negative controls in every run.

*** SEMIAUTOMATIC PROCEDURE**

Use the volumes reported for the manual method, and a mechanical rotating shaker, with a radius of 3 cm, at a speed of 100-130 rpm.

*** SEMIQUANTITATIVE METHOD IN SLIDE**

Dilute the sample by doubling dilutions in saline and perform the test in slide on the different dilutions.

9. TEST VALIDATION

After 8 minutes incubation time, the positive control must present evident, total agglutination; the negative control must present no agglutination.

10. INTERPRETATION OF THE TEST*** RESULTS MANUAL/SEMIAUTOMATIC METHOD**

POSITIVE: Evident agglutination against the white background of the slide.

NEGATIVE: No agglutination; homogeneous appearance.

*** RESULTS SEMIQUANTITATIVE METHOD**

The titre of the serum sample is give by multiplying the dilution by the sensitivity limit 1/40.

Example:

Dilution	1:2	1:4	1:8
Slide Test	++	+	-
Titer	$1/2 \times 1/40 = 1/80$		

Samples with titers of 1/40 should be considered borderline positive.

11. INTERFERING SUBSTANCES

The following substances do not interfere in the test :

- Rheumatoid factor from 60 to 500 IU/mL
- Autoantibodies from 16 to 12.000 EU/mL
- Triglycerides <= 870 mg/dl
- Bilirubin <= 2,5 UI
- Heterophile antibodies <= 1:448
- Hemoglobin <= 30 mg/dl
- Anti-Borrelia burgdorfi antibodies

12. LIMITATIONS OF THE TEST

False positive results can be found in women during pregnancy.

13. PERFORMANCE

200 serum samples from healthy donors were analysed with Syphilis Fast and with an immunoenzymatic screening method taken as reference. Seven false positive results were found. The diagnostic specificity was 96.5%.

105 positive sera were tested, and 3 false negatives were found. The diagnostic sensitivity was 97.1%.

No prozone effect was observed with titers >= 81,920 TPHA.

14. ANALYTICAL SENSITIVITY

The sensitivity of the method is equal to 0.1 IU/mL, with reference to the 1st International Standard Preparation, WHO 1958.

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- 1.Luger A.F.H. In: Serological diagnosis of sexually transmitted diseases. New York: M. Dekker Inc. 1988, 249-674.