Vitassay qPCR

Zika + Dengue + Chikungunya

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- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96 TM Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)

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Termocicladores con bloque de bajo perfil	Termocicladores con bloque de alto perfil
Agilent Technologies	Abbott
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ / CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler ®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler ®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System*
	DTprime Real-time Detection Thermal Cycler*
	Eppendorf
	Masterycler™ ep realplex
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

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Intended use

Vitassay qPCR Zika+Dengue+Chikungunya test allows the detection and differentiation of Zika, Dengue and/or Chikungunya viruses by real-time RT-PCR in clinical samples. The product is intended for use in the diagnosis of Zika virus, Dengue virus and/or Chikungunya virus infections alongside patient clinical data and other laboratory tests outcomes.

References

Vitassay qPCR Zika+Dengue+Chikungunya 4x8 -well strip, low profile 7041005

Vitassay qPCR Zika+Dengue+Chikungunya 4x8-well strip, high profile 7042005

Materials/reagents provided

Reference	Reagent/Material	Colour	Amount
7041S005/ 7042S005	Zika+Dengue+Chikungunya Virus strips low/high profile	-	4 x 8-well strip
7C005	Zika+Dengue+Chikungunya Positive Control	red	1 vial
7001A	PCR grade water	white	1 vial x 1 mL
7002B	Resuspension Buffer	green	1 vial x 1,8 mL
7003N	Negative control	yellow	1 vial x 1 mL
7004O	Optical caps	-	4 x 8 cap strip

Transport and storage

- The reagents and the test can be shipped and stored at 2-40°C until expiration date stated in the label.
- The resuspended positive control should be stored at -20°C. In order to avoid repeated freeze/thaw cycles, it is recommended to distribute the content in different aliquots.
- Keep all reagents in the darkness.

Additional equipment and material required

- RNA extraction kit
- Real-time PCR instrument (thermocycler) (Attached I)
- Centrifuge for 1.5 mL tubes
- Vortexer

- Micropipettes (1-20 µL, 20-200 µL)
- Filter tips
- Powder-free disposal gloves

Summary

Zika (ZIKV), Dengue (DENV) and Chikungunya (CHIKV) viruses are transmitted by the sting of the same *Aedes* genus mosquito.

Besides, Zika virus was isolated in 1947 from a rhesus monkey in the Zika forest in Uganda. Typically, clinical manifestations of human ZIKV disease are similar to many others derived of arbovirus infections that occur without serious complications and may include a self-limiting febrile illness, arthralgia, myalgia, headache and maculopapular rash. Conjunctivitis, retroorbital eye pain, lymphadenopathy and diarrhoea have also been reported.

Dengue is also an acute viral illness caused by RNA virus of the family Flaviviridae. Presenting features may range from asymptomatic fever to dreaded complications such as haemorrhagic fever and shock. Acute onset high fever, muscle and joint pain, myalgia, cutaneous rash, haemorrhagic episodes, and circulatory shock are the commonly seen symptoms.

Chikungunya virus was isolated in Tanzania in 1953. Chikungunya virus causes severe, debilitating and often chronic arthralgia with high attack rates, resulting in severe morbidity and economic costs to affected communities.

These three arboviruses in the initial stages of the infection cause similar clinical presentations, although the disease management is different, for that, an early detection and differentiation is crucial. Confirmation and identification of these infections is based mostly on viral RNA detection in serum and urine clinical samples by using reverse transcription PCR (RT-PCR). In fact, serology is not advised since this technique can cause cross reactivity between the antibodies of these arbovirus.

Principle of the test

Vitassay qPCR Zika+Dengue+Chikungunya test is based on the real-time amplification of specific conserved fragments of the *envelope* gene (Zika virus), of the 3' Non-coding region (Dengue virus) and of the *NSP1* gene (Chikungunya virus). The viral RNA extracted is transcribed into cDNA using a specific primer by reverse transcription step followed immediately by polymerase chain reaction. The presence of the virus is detected by a fluorescence increment during the reaction upon hydrolysis of the fluorescent probe.

Vitassay qPCR Zika+Dengue+Chikungunya test is a ready-to-use test which contains in each well all the necessary reagents for real-time PCR assay in a stabilized format. In addition, an internal control allows the detection of a possible reaction inhibition. The amplification of the Zika virus RNA target sequence is detected through the Cy5 channel, Dengue virus RNA target in FAM channel, Chikungunya virus RNA target in ROX channel whereas the internal control (IC) in HEX, VIC or JOE channel (depending on the equipment used select the proper detection channel, see Annex 2).

Precautions

- For *in vitro* diagnostic use.
- Do not use after expiration date.
- Do not mix components from different kits and/or lots.
- Do not use if package is open or damaged.
- Do not use the kit if the desiccant from the different reagents is not present or broken or if the foil has been broken or damaged.
- Protect the kit against humidity. Prolonged exposure to humidity may affect product performance.
- Follow Good Laboratory Practices. Wear protective clothing, use disposal gloves, goggles and mask.
- Do not eat, drink or smoke in the working area.
- The laboratory process must be one-directional, it should begin in the Extraction Area and then move to the Amplification and Detection Areas. Do not return samples, equipment and reagents to the area in which the previous step was performed.
- Specimens and reagents/materials that have been exposed to them must be treated as potentially infectious. Take necessary precautions during the collection, storage, treatment and disposal of samples.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.

Procedures

Specimen collection, processing and RNA extraction

For pre-treatment and RNA isolation, it is recommended to use your optimized manual or automatic system compatible with real-time PCR technology. The assay has been validated with the following extraction kits:

Maxwell® 16 Viral Total Nucleic Acid Purification Kit (Promega).

QIAamp Viral RNA Mini Kit, QIAcube instrument (QIAGEN).

High Pure Viral Nucleic Acid Kit (Roche).

Invisorb® Spin Universal Kit (Stratec).

Positive control preparation

Reconstitute the lyophilized Zika+Dengue+Chikungunya Positive Control (red tube) in 100 µL of PCR grade water (white tube). To ensure a complete resuspension, vortex the tube thoroughly. After the first use, dispense into aliquots in order to avoid multiple freeze-thaw cycles and store them at -20°C.

This component contains high copies number template and is a very significant contamination risk. Therefore, we recommend to open and to manipulate it in a separate laboratory area away from the other components and clinical samples.

Reaction setup

- Separate the number of required reactions including samples and controls. Remember that one positive and one negative control must be included in each run.
- Peel off protective aluminium seal from the strips.
- Pipette 15 µL of Resuspension buffer (green tube) and add it into each well.
- Pipette 5 µL of RNA sample, negative and positive controls and add it into each well.
- Cover the wells with the caps provided. Spin down briefly (optional).
- Place the strips in the Real-time PCR instrument.

Programme your thermocycler

Set your thermocycler following the conditions indicated below:

Step	Temperature	Time	Cycles
Reverse transcription	45°C	15 min	1
Initial denaturation	95°C	2 min	1
Denaturation	95°C	10 sec	45
Annealing/Extension (Data collection*)	60°C	50 sec	

Set the fluorescence data collection during the extension step (*) through the Cy5 (Zika virus), FAM (Dengue virus), ROX (Chikungunya) and HEX, JOE or VIC channels

(Internal Control (IC)). If you use the Applied Biosystems 7500 Fast or the Stratagene Mx3005P™ check that passive reference option ROX is none. (Attached 2).

Analysis and interpretation of results

The analysis of the results is done by the software itself of the used real-time PCR system following manufacturer's instructions.

For a valid diagnostic test run, the following control conditions must be met:

Positive control

The positive control used in each run, must show an amplification in channels: Cy5 (Zika virus), FAM (Dengue virus) and ROX (Chikungunya), which validates the reaction.

Negative control

The negative control included in each run, must show the absence of signal in FAM, ROX and Cy5 channels, which validates the reaction.

The experiment seems to be failed if there is signal of amplification in negative control or absence of signal in the positive well. The assay should be repeated.

The result interpretation is summarized in the following table:

Zika Virus	Dengue Virus	Chikungunya Virus	Internal control	Negative control	Positive control	Interpretation
+	+	+	+/-	-	+	Zika, Dengue and Chikungunya viruses Positive
-	-	-	+	-	+	Zika, Dengue and Chikungunya viruses Negative
+	-	-	+/-	-	+	Zika Virus Positive, Dengue and Chikungunya Viruses Negative
+	+	-	+/-	-	+	Zika and Dengue Viruses Positive, and Chikungunya Virus Negative
+	-	+	+/-	-	+	Zika and Chikungunya Viruses Positive, and Dengue Virus Negative
-	+	-	+/-	-	+	Dengue Virus Positive, Zika and Chikungunya Viruses Negative
-	+	+	+/-	-	+	Dengue and Chikungunya Viruses Positive, Zika Virus Negative
-	-	+	+/-	-	+	Chikungunya Virus Positive, Zika and Dengue Viruses Negative
+	+	+	+	+	+	Experiment fail
-	-	-	-	-	-	Experiment fail

(+) Positive: Amplification signal

(-) Negative: No amplification signal

If the negative samples do not show a positive result for the internal control, they should be retested using a dilution 1:10 from the original sample or the nucleic acid extraction has to be repeated due to possible problems caused by PCR inhibitors.

Quality Control

An Internal Control (IC) is included in each reaction for confirming the appropriate performance of the molecular diagnostic technique. Besides, a positive and a negative control must be included in each assay to interpret the results correctly.

Performance evaluation

Clinical sensitivity and specificity

A total of 102 samples from EQA panels were analyzed using Vitassay qPCR Zika + Dengue + Chikungunya. The results were compared with the final reports of the EQAs programs. Vitassay qPCR Zika + Dengue + Chikungunya detected 25 positive samples for the Zika virus, 27 positive samples for the Dengue virus and 29 positive samples for the Chikungunya virus.

The results show a high sensitivity and specificity to detect the Zika, Dengue and Chikungunya viruses using Vitassay qPCR Zika + Dengue + Chikungunya kit.

Analytical sensitivity

The analytical sensitivity was determined by analysis of 10-fold dilution series of Zika, Dengue and Chikungunya virus templates ranging from 10^7 to 10^1 copies/rxn. This assay has a detection limit of ≥ 10 viral RNA copies per reaction.

Analytical specificity

The analytical specificity for Zika, Dengue and Chikungunya viruses was tested within the panel of the following viruses, where no cross-reactivity was seen between any of the species:

Zika	Dengue	Chikungunya
Chikungunya S27 Petersfield	Chikungunya S27 Petersfield	Zika MR 766 (Uganda)
Dengue 1 Hawaii	Zika MR 766 (Uganda)	Zika 11474/16 (French Polynesia)
Dengue 2 New Guinea C	Zika 11474/16 (French Polynesia)	Zika 11468/16(French Polynesia)
Dengue 3 H87	Zika 11468/16(French Polynesia)	Zika Virus (African)
Dengue 4 H241	Zika Virus (African)	Zika PF13/251013-18 (Asian)
Plasmodium falciparum	Zika PF13/251013-18 (Asian)	Dengue 1 Hawaii
Japanese encephalitis	Japanese encephalitis	Dengue 2 New Guinea C
St Louis Encephalitis 17D	St Louis Encephalitis 17D	Plasmodium falciparum
West Nile H160/99	West Nile H160/99	Trypanosoma cruzi
West Nile Heja	West Nile Heja	Dengue 3 H87
West Nile Ug37	West Nile Ug37	Dengue 4 H241
Yellow Fever 17D	Yellow Fever 17D	St Louis Encephalitis 17D
Trypanosoma cruzi	Plasmodium falciparum	West Nile H160/99
-	Trypanosoma cruzi	West Nile Heja
-	-	West Nile Ug37
-	-	Yellow Fever 17D
-	-	Japanese encephalitis

Analytical reactivity

The reactivity of Vitassay qPCR Zika+Dengue+Chikungunya Virus Real Time PCR assay for Zika virus was confirmed by the real time amplification using Zika virus strain MR 766 (Uganda, 1947), Zika Virus strain 11474/16 (French Polynesia), Zika Virus strain 11468/16(French Polynesia), Zika Virus (African) and Zika virus strain PF13/251013-18 (Asian) as template.

The reactivity of Vitassay qPCR Zika+Dengue+Chikungunya Virus Real Time PCR for Dengue virus was confirmed by the real time amplification using the following strains as templates: Dengue 1 virus strain Hawaii, Dengue 2 virus strain New Guinea C, Dengue 3 virus strain H87 and Dengue 4 virus strain H241.

The reactivity of Vitassay qPCR Zika+Dengue+Chikungunya Virus Real Time PCR for Chikungunya virus was confirmed by the real time amplification using Chikungunya virus S27 Petersfield (African genotype), and Chikungunya virus Martinique isolate (Asian genotype) as template.

Compatibles real-time PCR equipment

Vitassay qPCR Zika+Dengue+Chikungunya has been validated on the following equipments:

- Cobas Z480 (Roche)
- 7500 Fast Real-Time PCR System (Applied Biosystems) ^{II}
- CFX96 TM Real-Time PCR Detection System (Bio-Rad)
- AriaMx Real-Time PCR System (Agilent Technologies)
- DTlite Real-Time PCR System (DNA-Technology)
- DTPrime Real Time Detection Thermal Cycler (DNA-Technology)
- Rotor-Gene® Q (Qiagen)^I
- SmartCycler® (Cepheid)^I

I: For Rotor-Gene® Q and SmartCycler® thermocyclers the product should be reconstituted following the procedure and transferred into specific Rotor-Gene® Q and/or SmartCycler tubes.

II: When using the Applied Biosystems 7500 Fast with strips it is recommend placing a plate holder for the tubes (Ref. PN 4388506).

Limitations

- This test provides a presumptive diagnosis of Zika, Dengue and Chikungunya virus infections. All the obtained results must be interpreted together with complementary clinical information and laboratory findings available to the physician.
- This assay was tried with serum and urine samples. The use of other samples has not been established.
- The quality of the test depends on the quality of the sample; proper RNA from clinical specimens must be extracted. Unsuitable collection, storage and/or transport of specimens may give false negative results.
- Extremely low levels of target below the limit of detection may be detected, but results may not be reproducible.
- There is a possibility of false positive results due to cross-contamination by Zika, Dengue and/or Chikungunya viruses, either samples containing high concentrations of target RNA or contamination due to PCR products from previous reactions.

Attached I: Compatibility of the most common real-time PCR equipment

The most common thermocyclers are listed in the following table separated by tube type. If you do not find your thermocycler, please contact with your supplier.

Low profile Block Thermocyclers	High profile Block Thermocyclers
Agilent Technologies	Abbott
AriaMx/AriaDx Real-Time PCR System	Abbott m2000 RealTime System
Applied Biosystems	Applied Biosystems
7500 Fast Real-Time PCR System	7300 Real-Time PCR System
7500 Fast Dx Real-Time PCR System	7500 Real-Time PCR System
QuantStudio™ 12K Flex 96-well Fast	7900 HT Real-Time PCR System
QuantStudio™ 6 Flex 96-well Fast	ABI PRISM 7000
QuantStudio™ 7 Flex 96-well Fast	ABI PRISM 7700
QuantStudio™ 5 Fast Real-Time PCR System	QuantStudio™ 12K Flex 96-well
QuantStudio™ 3 Fast Real-Time PCR System	QuantStudio™ 6 Flex 96-well
StepOne Plus™ Real-Time PCR System	QuantStudio™ 7 Flex 96-well
StepOne™ Real-Time PCR System	QuantStudio™ 5 Real-Time PCR System
ViiA™ 7 Fast Real-Time PCR System	QuantStudio™ 3 Real-Time PCR System
BIONEER	ViiA™ 7 Real-Time PCR System
Exicycler™ 96	Analytik Jena Biometra
Bio-Rad	TOptical
CFX96™ / CFX96™ IVD Real-Time PCR Detection System	qTOWER 2.0
Mini Opticon™ Real-Time PCR Detection System	BIONEER
Cepheid	Exicycler™ 96
SmartCycler®	Bio-Rad
Qiagen	CFX96™ / CFX96™ IVD Real-Time PCR Detection
Rotor-Gene® Q	iCycler iQ™ Real-Time PCR
Roche	iCycler iQ™ 5 Real-Time PCR
LightCycler ®480 Real-Time PCR System	MyiQ™ Real-Time PCR Detection System
LightCycler ®96 Real-Time PCR System	MyiQ™ 2Real-Time PCR Detection System
Cobas z480 Analyzer	Cepheid
	SmartCycler®
	DNA-Technology
	DTlite Real-Time PCR System
	DTprime Real-time Detection Thermal Cycler
	Eppendorf
	Masterycler™ ep realplex
	Qiagen
	Rotor-Gene® Q
	Stratagene / Agilent Technologies
	Mx3000P™ Real Time PCR System
	Mx3005P™ Real Time PCR System

* See Attached III to configure exposure settings.

Attached II: Detection channels of most common real time PCR equipment

The fluorescence detection channels of some of most common Real Time PCR Thermocyclers are specified in the following table:

REAL-TIME PCR THERMOCYCLER	VITASSAY CHANNEL	DETECTION CHANNEL	OBSERVATIONS
Bio-Rad CFX96™	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
ABI 7500 Applied Biosystems	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Roche Lightcycler®480II	FAM	465/510	Colour Compensation required
	HEX	533/580	
	ROX	533/610	
	Cy5	618/660	
Smartcycler® Cepheid	FAM	Channel 1	
	HEX	Channel 2	
	ROX	Channel 3	
	Cy5	Channel 4	
Abbott m2000rt	FAM	FAM	
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
Mx3000P™ Mx 3005P™ Stratagene	FAM	FAM	Passive reference option ROX is none
	HEX	VIC	
	ROX	ROX	
	Cy5	Cy5	
AriaMx Agilent	FAM	FAM	
	HEX	HEX	
	ROX	ROX	
	Cy5	Cy5	
Rotor-Gene® Q Qiagen	FAM	Green	In the Channel Setup, click on the "Gain Optimisation" button and then go to "Optimise Acquiring". The fluorescence Target Sample Range has to be between 5 and 10 FI for each channel. Also select the option "Perform Optimisation Before 1st Acquisition".
	HEX	Yellow	
	ROX	Orange	
	Cy5	Red	
Exicycler™ 96 BIONEER	FAM	FAM	
	HEX	JOE	
	ROX	ROX	
	Cy5	Cy5	

Attached III: Optical measurement exposure setting

The exposure values of some thermocyclers must be adjusted for proper operation. Set exposure values as follows:

Thermocycler	Vitassay channel	Exposure values
DTlite Real-Time PCR System (DNA-Technology)	FAM	500
	HEX	500
	ROX	500
	Cy5	500
DTprime Real-time Detection Thermal Cycler (DNA-Technology)	FAM	500
	HEX	1000
	ROX	1000
	Cy5	1000

Bibliography/Bibliografía

1. Genetic and serologic properties of Zika virus associated with an epidemic, Yap State, Micronesia, 2007 Lanciotti RS, Kosoy OL, Laven JJ, Velez JO, Lambert AJ, Johnson AJ, Stanfield SM, Duffy MR. *Emerg Infect Dis.* 2008 Aug;14(8):1232-9.
2. Quantitative real-time PCR detection of Zika virus and evaluation with field-caught mosquitoes. Faye O, Faye O, Diallo D, Diallo M, Weidmann M, Sall AA. *Virol J.* 2013 Oct 22;10:311.
3. Dengue virus: A global human threat: Review of literature. Hasan S, Jamdar SF, Alalowi M, Al Ageel Al Beaiji SM. *J Int Soc Prev Community Dent.* 2016 Jan-Feb 6(1):1-6.
4. One-step real-time RT-PCR assays for serotyping dengue virus in clinical samples. Alm E, Lindegren G, Falk KI, Lagerqvist N. *BMC Infect Dis.* 2015 Nov2;15:493.
5. Pan-dengue virus detection by PCR for travelers returning from the tropics. Dumoulin A, Marti H, Panning M, Hatz C, Hirsch HH. *J Clin Microbiol* 2008 Sep;46(9):3104-6.
6. External quality assessment studies for laboratory performance of molecular and serological diagnosis of Chikungunya virus infection. Jacobsen S, Patel P, Schmidt-Chanasit J, Leparc-Goffart I, Teichmann A, Zeller H, Niedrig M. *J Clin Virol* 2016 Mar;76:55-65.

Trademarks

CFX™ and IQ5™ are registered trademarks of Bio-Rad Laboratories.

ABI®, QuantStudio™ and ViiA™ are registered trademarks of Thermo Fisher Scientific Inc.

LightCycler® is a registered trademark of Roche.

Mx3000P™ and Mx3005™ are registered trademarks of Agilent Technologies.

Mastercycler™ is a registered trademark of Eppendorf.

Rotor-Gene® Q is a registered trademark of Qiagen.

SmartCycler® is a registered trademark of Cepheid.

Symbols for IVD components and reagents/ Símbolos utilizados para componentes y reactivos IVD

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