MAST ISOPLEX® *Malaria PfPv*

MAST ISOPLEX® *Malaria PfPv*

REF 67MALA-2

96 Tests

UDI-DI 4250729702501

**Gebrauchsanweisung /
Instructions for Use /
Instructions d'utilisation**

**Nur zur Verwendung durch Fachpersonal /
For professional use only /
Pour un usage professionnel uniquement**

Deutsch

Seiten

02–07

English

Pages

08–12

Français

Pages

13–17

MAST ISOPLEX® Malaria PfPv

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MAST ISOPLEX® Malaria PfPv

FOR IN VITRO DIAGNOSTIC USE

For professional use only**Intended purpose**

MAST ISOPLEX® Malaria PfPv is an in vitro diagnostic kit for qualitative detection of *Plasmodium falciparum* and *Plasmodium vivax* parasites in human blood samples.

It is intended to be used as an aid in the diagnosis of *P. falciparum* and *P. vivax* in symptomatic patients or individuals suffering from febrile episodes.

Background

Malaria infection is caused by the *Plasmodium* spp. pathogen. The pathogen can persist without symptoms appearing. Since there are many fever-causing pathogens in endemic areas worldwide, it is important to make an accurate diagnosis of which pathogen is responsible for the symptom of fever. The MAST ISOPLEX® Malaria PfPv assay has been adjusted to detect only those malaria infections that are responsible for the symptomatic fever episodes.

Test principle

MAST ISOPLEX® Malaria PfPv is based on the principle of loop-mediated isothermal amplification (LAMP, 1). The assay is used for the direct, qualitative detection of DNA from the malaria pathogens *P. falciparum* and *P. vivax* in human blood samples.

In the MAST ISOPLEX® Malaria PfPv, the DNA of the protozoa *P. falciparum* or *P. vivax* is amplified by a LAMP reaction. The use of novel "Mediator Displacement Probes" in combination with universal reporters (2) allows real-time multiplexing of the specific target sequences.

The kit contains all the components necessary for a LAMP reaction (specific primers, probes, universal reporters, dNTPs, buffer and DNA polymerase with strand displacement activity) as well as a Positive Control (PC) and an Internal Inhibition Control DNA (IC), respectively, to monitor the sample preparation and/or identify a possible inhibition of the amplification reaction.

The *P. falciparum* specific target sequence of the mitochondrial region is amplified and detected in the red channel (Cy5). The *P. vivax* specific target of the 18S ribosomal gene is amplified measured in the green channel (FAM). The IC target sequence is amplified and detected in the orange channel (ROX).

For isothermal amplification and detection standard real-time PCR instruments can be used or specific programmable real-time heat blocks capable to read fluorescent signals.

The amplified product can be detected within 45 minutes.

Kit Contents

The reagents of one kit are sufficient for 96 tests.

Table 1:

Kit code	Content	Volume	Lid colour	Additional packaging information
CONTROL +	Positive Control	1x 250 µL	red	Lyophilized, to be reconstituted before use
CONTROL IC	Internal Control	1x 50 µL	blue	Lyophilized, to be reconstituted before use
PEL	Lyo pellet PfPv	8x pellets for 12 reactions each	yellow	Lyophilized, to be reconstituted before use
RB-PEL	Reconstitution Buffer Pellets	2x 1.1 mL	black	RTU
RB-CONTROL	Reconstitution Buffer Controls	2x 1.1 mL	transparent	RTU

Further abbreviations

RTU ready to use

Materials required but not Provided

- Appropriate real-time PCR instrument or equipment capable of isothermal nucleic acid amplification and detection of amplified products using fluorescence (RotorGene 4 plex or higher from Qiagen has been validated, other instruments should be validated before use)
- Appropriate nucleic acid extraction system or kit (MAST ISOPLEX® DNA / RNA Extraction Kit has been validated, other extraction system should be validated before use)
- Standard DNase free supplies such as reaction tubes, pipette tips with filters
- Molecular grade water (DNase / RNase free)
- Appropriate reaction tubes for cyclor model or 96 well reaction plates with corresponding closing material
- Vortex mixer
- Tabletop centrifuge with rotor to hold 1.5 mL assay reaction tubes
- Centrifuge with rotor for microtiter plates, if using 96 well reaction plates
- Ice bucket with ice or cold block
- Disposable powder-free gloves
- Adjustable calibrated pipettes
- Personal protective equipment

Table 2: Recommended Real-time PCR equipment, consumables / extraction kits

Real-time PCR manufacturer	Model
Qiagen	RotorGene 4-plex or higher
Isothermal nucleic acid extraction systems	Type
MAST	MAST ISOPLEX® DNA / RNA Extraction Kit

Warnings

- Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
- The reagents of the MAST ISOPLEX® *Malaria PfPv* Kit are very sensitive, therefore take precautions to protect reagents and samples from contamination. Do not use contaminated reagents. Do not use reagents with damaged bottles or packaging.
- The kit and its components do not contain any infectious materials. The user, however, must take care when handling patient samples, because they may carry infectious organisms and viruses.
- The kit reagents may contain bovine components. To protect against infections with prions or zoonotic pathogens, standard occupational safety regulations (GLP) must be observed.

Precautions

- Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
- Do not use reagents after expiry date.
- Do not mix reagents between different lots, as reagents have been calibrated for each batch.
- Before pipetting, all reagents should be mixed thoroughly by gentle agitation.
- In accordance with Good Laboratory Practice (GLP) all laboratory equipment should be checked regularly regarding accuracy and precision.
- It is recommended to wear protective clothing e.g. lab coat, protective gloves and goggles while working with the kit, kit components and sample materials. Regional and workplace specific regulations must be observed.
- Reaction tubes should be kept closed at all times following addition of reagents and should be discarded, without opening following use, as per local guidelines.
- To avoid any contamination with the amplified product never open a vial after amplification in the environment where sample preparation or master mixes are prepared!

Disposal

- The tubes should be disposed of unopened after use and according to local guidelines. Potential infectious material should be inactivated before disposal.

Stability and Storage

- MAST ISOPLEX® *Malaria PfPv* kit is shipped at ambient temperatures. The components of the kit are freeze dried. If the kit box is damaged upon receipt, or if tubes have been compromised during shipment, contact your local distributor for assistance.
- All reagents should be stored at 2-8 °C and **used up within 4 weeks** after first opening of the vials. **Exception:** Reconstituted Lyo pellets shall be used immediately after reconstitution.
- Reconstituted IC and PC are stable for at least 2 weeks at 2-8 °C.
- Reagents should not be used after expiry!

Test Procedure

1. Specimen collection and preparation

MAST ISOPLEX® *Malaria PfPv* is intended for use with DNA extracted from human blood samples.

1. DNA Extraction

MAST ISOPLEX® *Malaria PfPv* was validated on samples after nucleic acid extraction. The assay is not validated on non-extracted specimen. Nucleic acid extraction systems should give a purity at a 260 / 280 nm ratio from 1.8 to 2.0.

2. Reconstitution of Positive and Inhibition Control

	Step	Conduction	
Positive control PC	Reconstitution	1. Bring the PC pellet to the bottom of the tube by briefly centrifuging the tube. 2. Pipette 250 µL RB-CONTROL into the tube and place on ice for 5 min to dissolve. The PC is sufficient for 50 reactions. 3. Vortex briefly to mix the reconstituted PC.	End volume 250 µL
Inhibition control IC	Reconstitution	1. Bring the PC pellet to the bottom of the tube by briefly centrifuging the tube. 2. Pipette 50 µL RB-CONTROL into the tube and place on ice for 5 min to dissolve. 3. Vortex briefly to mix the reconstituted IC.	End volume 50 µL

Reaction Preparation

Note:

- Only use reaction tubes which are recommended for your chosen amplification instrument
- Ensure reaction tubes are not scratched or cracked prior to use
- During assay set-up reagents should be kept on ice or in a cold block
- Do not mix reagents vigorously in the reaction tube
- Protect the reagents from direct sun light

A lyo-pellet (reaction pellet) is enough for 12 reactions.

1. Bring the Lyo pellet to the bottom of the tube by brief centrifugation.
2. Dissolve the pellet with 250 µL buffer.
3. Add 5 µL of the reconstituted IC control DNA.
4. Close the tube, mix gently by turning it upside down and back several times. Strong and long vortexing may have a negative effect on the polymerase.
5. Centrifuge the reaction tube briefly to collect liquids from the lid or walls, then place the master-mix tube on ice.
6. Pipette 20 µL of this master-mix into each amplification tube.

Addition of samples and controls

1. Add 5 µL of sample DNA extract into each master-mix tube, then immediately close the tube.
2. Positive Control reactions: Always include the PC provided in the kit in one tube during each run. Add 5 µL of PC to one tube instead of the patient sample.
Negative Control: Always include a Negative Control (RB-control) during each run. Add 5 µL RB-control to one tube instead of the patient sample. No template controls, or water can also be used as negative control.
3. Tightly seal the reaction tubes and spin in a mini centrifuge (holding for 8 strip tubes) to collect all liquids at the bottom of the vial.
4. Transfer reaction tubes into the amplification instrument and start the amplification reaction.

Instrument Setup

1. Define a program for isothermal amplification on the PCR or LAMP instruments.
2. Set amplification **temperature to 64 °C**.
3. Set amplification **time to 45 min**.
4. Select filters for FAM, Cy5 and ROX.
5. Allow a reading of the fluorescence signal every **30 or 60 sec**.

Control Procedure

Each reaction shall contain either an PC or IC target to prove correct assay performance.

In case of positive reactions with *P. falciparum* and / or *P. vivax* genes, the IC may show a reaction. Even without a signal of inhibition control the assay result can be valid (see table 2).

In case of a negative reaction with *P. falciparum* and / or *P. vivax* genes the IC must show a reaction.

The PC may show reactions in Cy5 and / or FAM. At least one channel must be positive. Additionally, the IC may show a positive reaction.

For correct signal appearance of IC see the Certificate of Analysis supplied with each batch.

Interpretation of Results

Figure 1: Typical amplification plots


Table 2: Interpretation criteria

Cy5 Pf	FAM Pv	ROX IC	Result	Run
+	+	+	positive for <i>P. falciparum</i> and <i>P. vivax</i>	valid
+	-	+	positive for <i>P. falciparum</i>	valid
+	+	-	positive for <i>P. falciparum</i> und <i>P. vivax</i>	valid
-	+	+	positive for <i>P. vivax</i>	valid
-	+	-	positive for <i>P. vivax</i>	valid
+	-	-	positive for <i>P. falciparum</i>	valid
-	-	+	negative	valid
-	-	-	not interpretable	invalid

A positive result is indicated by an amplification plot rising from background exponentially. In a negative test a flat line is proving that there is no amplification.

Positive:

A positive result is indicated by an amplification plot rising from background exponentially. As example see figure 1.

A positive reaction in the red (Cy5) channel indicates the presences of *P. falciparum*. A positive reaction in the green (FAM) channel indicates the presence of *P. vivax*. A positive reaction in both red and green channel indicates the presence of both *P. falciparum* and *P. vivax*.

The internal control reaction shows a positive amplification in most cases in the orange (ROX) channel. However, if the signal does not appear in the presence of a *P. falciparum* and / or *P. vivax* signal than the result is still valid.

Negative:

An IC reaction must occur in any true negative reaction, to prove that no inhibitions occurred.

The negative amplification plot of samples shows a flat line along the x-axis, a slight continuous increase of the curve must be read as negative. Sample quality may lead to a higher background causing an increase in signal performance, however, never rising to a peak during the amplification time. As an example, see in Fig. 1 the flat line in each of the channels. A slight permanent increase of the curve (ramping effect) does not mean a positive signal.

Limitations / Interferences

- The MAST ISOPLEX® *Malaria PfPv* LAMP assay is designed to detect symptomatic patients in their febrile disease state. **The assay is not designed to be used to specifically detect Malaria pathogens in asymptomatic patients and in silent carriers.**
- All assay results must be interpreted together with other clinical and prevalence information of malaria. In case of doubt other clinical and diagnostic data must be taken into consideration for a final diagnosis. A negative result may occur if the concentration of Plasmodium DNA in a sample is below the detection limit of the test. A negative test result therefore does not eliminate the possibility of an infection with Plasmodium parasites.
- The MAST ISOPLEX® *Malaria PfPv* LAMP assay is only detecting *P. falciparum* and *P. vivax*. Other Plasmodium species may be present in a patient sample depending on the area of infection.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results do not prove the presence or absence of infections with other causative agents of febrile diseases.
- Positive and negative predicted values are highly dependent on prevalence. False negative test results are more likely when prevalence of disease is high. False positive test results are more likely when prevalence is moderate to low.
- The assay has not been validated on non-extracted blood samples.

Comments on false positive, false negative and query results:

False positive results:

False positive are hard to detect, they look like true positive signal. Possible reasons could be:

- False positive results may occur in rare cases at the end of the amplification. They may be due to primer multimers which start to amplify.
- Cross contamination may lead to false-positive results, even aerosol carrying the amplification product reliably lead to false positive reactions.

False negative results:

False negative results are those, when the IC is giving a signal but the expected target regions in *P. falciparum* and *P. vivax* are not reacting. Such effects might be caused by:

- Sample volume < 5 µL
- DNA copies are below the limit of detection.
- IC (volumes significantly higher than 1 µL per reaction) compete with the target genes, delaying their reaction to become positive with the amplification time
- Incorrect amplification temperature

Performance Characteristics

MAST ISOPLEX® *Malaria PfPv* detects a mitochondrial target region of *P. falciparum* and an 18S ribosomal region of *P. vivax*.

Cross-reactivity

Potentially cross-reacting organisms has been evaluated in a scientific publication (reference 3). The cross reactivity was tested for other *Plasmodium* species as well as organisms associated with febrile symptoms.

- *Plasmodium ovale*
- *Plasmodium malariae*
- *Plasmodium knowlesi*
- Dengue virus serotypes 1 and 2
- Chikungunya virus
- *Salmonella enterica* Typhi
- *Salmonella enterica* paratyphi A

Traceability

The WHO standard *Plasmodium falciparum* 04/176 (NIBSC) can be used to calibrate the assay and check the reactivity of the assay.

Limit of detection

The assay has a limit of detection of 22.4×10^3 copies for *P. falciparum* and 300 copies for *P. vivax* per reaction.

Precision

The precision was calculated using the Ct values of two positive samples (strong and weak). The Ct value is the time at which a reaction becomes positive. These only roughly correlate with the copy number of the pathogens. For example, low pathogen concentrations can sometimes show a positive result earlier than expected. With high copy numbers, the Ct values are often closer together.

- Intra-assay precision (%):

FAM	Cy5	ROX (internal control)
0.76-3.46%	3.01-65.53%	2.59-7.95%

- Inter-assay precision (%):

FAM	Cy5	ROX (internal control)
17.16%	7.65-21.47%	18.89%

Measuring Range

The measuring range depends on the signal strength units of the respective amplification instrument and cannot be generalised.

Clinical Performance

To evaluate diagnostic performance, the "thick blood drop" was used as the gold standard.

A collective of 38 positive and 56 negative samples were analysed for *P. falciparum*.

A collective of 3 positive and 91 negative samples were analysed for *P. vivax*.

All samples were collected in an endemic area by a national central laboratory (Central Laboratory, Khartoum, Sudan).

	<i>Formula</i>	<i>P. falciparum</i>	<i>P. vivax</i>
Sensitivity	$\frac{TP}{(TP + FN)}$	95 %	>99 %
Specificity	$\frac{TN}{(TN + FP)}$	95 %	>99 %
Positive predictive value	$\frac{TP}{(TP + FP)}$	92 %	>99 %
Negative predictive value	$\frac{TN}{(TN + FN)}$	96 %	>99 %
Positive likelihood ratio	$\frac{Sens}{(1 - Spec)}$	17.68	>99
Negative likelihood ratio	$\frac{(1 - Sens)}{Spec}$	0.06	<0.01
Efficiency	$\frac{(TP + TN)}{Total}$	95 %	>99 %

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**Mast Diagnostica GmbH,**
Mast Group Ltd.