MASTAFLUOR™ FTA-ABS IgM

MASTAFLUOR™ FTA-ABS IgM

REF 630525

10 x 10 Tests

UDI-DI 4250729700071

**Gebrauchsanweisung /
Instructions for Use /
Notice d'utilisation /
Istruzioni per l'uso**

**Nur zur Verwendung durch Fachpersonal /
For professional use only /
Exclusivement pour un usage professionnel/
Da utilizzare solo da parte di personale specializzato**

CE 2797

	Deutsch	Seiten	02-06
	English	Pages	07-11
	Français	Pages	12-16
	Italiano	Pagine	17-21

MASTAFLUOR™ FTA-ABS IgM

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MASTAFLUOR™ FTA-ABS IgM

Immunofluorescence assay for the detection of IgM antibodies to *Treponema pallidum* in human serum and plasma.

Intended Purpose

Semi-quantitative (titratable) immunofluorescence assay for the detection and confirmation of IgM antibodies directed against *Treponema pallidum* in human serum and plasma as an aid to diagnosis of syphilis.

The assay is suitable for manual use and is intended for professional in-vitro diagnostic use only. All laboratory test results should be interpreted in conjunction with other clinical data. The clinical judgement and further tests have to be taken into account additionally.

Important Note for Use of these Kit Instructions

Any assay relevant changes to the kit instructions for use (IFU) will lead to a change of the version number at the bottom on the last page. All changes made will be identified on a separate sheet added to the IFU for a period of three months from the date of a change of the version. Please ensure that the latest version of the IFU is used for the assay procedure.

Principle of the Test

Diluted specimens, the negative and positive controls or calibrators are applied to the wells of microscope slides and incubated. All wells are coated with purified antigens of the very pathogen. If specific antibodies are present in the serum they will bind to the fixed pathogen antigens forming a stable antigen-antibody complex. Slides are then washed to remove any unbound material. Complexed antibodies are detected by the addition of a fluorescence labelled anti-human IgM immunoglobulin conjugate. After a further washing step to remove any unbound conjugate, slides are viewed under a fluorescence microscope.

A green fluorescence is observed if pathogen specific antibodies are present in the sample material.

Kit Contents

1. **SLIDE** Slides

Ten slides with 10 wells each coated with fixed *Treponema pallidum* cells (Nichols strain).

2. **CONTROL+** Positive control

1 x 0.5 mL human serum with *Treponema pallidum* specific antibodies (IgM), ready to use, contains < 0.1% NaN₃ and ≤ 0.005% ProClin (potentially infectious and irritating to skin, see warnings)

3. **CONTROL-** Negative Control

1 x 0.5 mL of *Treponema pallidum* IgM negative control, human, ready to use, contains < 0.1% NaN₃ and ≤ 0.005% ProClin (potentially infectious and irritating to skin, see warnings)

4. **CONJ|M** FITC Conjugate IgM

1 x 3 mL of anti-human immunoglobulin (μ-chain, rabbit) conjugated to fluorescein isothiocyanate (FITC), ready to use, contains < 0.1% NaN₃ and ≤ 0.005% ProClin (irritating to skin, see warnings)

5. **BUFFER** PBS

2 sachets with 10 g PBS powder. Dissolve 1 sachet (10 g) in 1 litre distilled or deionised water to make a solution of phosphate buffered saline (PBS) at pH 7.2 ± 0.2.

6. **MOUNTING MEDIUM** Mounting Medium

1 x 3 mL of a buffered glycerol mounting medium, contains < 0.005% Proclin (irritating to skin, see warnings), ready-to-use

7. **SORBENT** Sorbent/Ultrasonicate

1 x 4 mL culture extract of *Treponema reiteri*, ultrasonicate, liquid, ready to use, contains sodium azide (< 0.1%) and ≤ 0.005% ProClin (irritating to skin, see warnings)

8. **MASTSORB** Mastsorb

1 x 6 mL rheumatoid factor absorbent, Anti-human-IgG, goat, ready to use, contains sodium azide (< 0.1%)

Further Abbreviations

1. **RTU** ready to use

Materials Required but not Provided

1. Sterile test tubes.
2. Micropipettes and tips.
3. Staining dish or Coplin jar.
4. Moist chamber.
5. Volumetric flask.
6. Distilled water or deionised.
7. Forceps.
8. Cover slips.
9. Wash bottle.
10. Fluorescence microscope with a filter combination suitable for FITC (e.g. 490 nm excitation filter and a 510nm barrier filter).

Warnings

1. Comply with the general health and safety guidelines for working with potentially infectious materials. Wear appropriate protective clothing and use appropriate lab facilities.
2. Do not cross-contaminate reagents or interchange caps on bottles. Do not use contaminated reagents. Do not use reagents with damaged bottles or packaging. Use a separate pipette or pipette tips for each sample and reagent.
3. The control/calibrator materials of human origin provided have been tested for antibodies to HIV and for HBsAg and found to be negative. However, they should be treated as potentially infectious materials and capable of transmitting diseases. No guarantee is given that the samples are free of infections or microbial contamination.
4. The kit reagents may contain bovine components. To protect against infections with prions or zoonotic pathogens, standard occupational safety regulations (GLP) must be observed.
5. Do not use haemolysed, lipemic or icteric samples (Lit. 10, 11, 12).
6. Sodium azide and ProClin are used as preservatives. ProClin can cause skin irritation. In the event of contact with the skin, the affected areas should be washed thoroughly with soap and water. Sodium azide can lead to poisoning if ingested and can form explosive metal azides with metal compounds (copper, lead). When disposing of azide reagents down the laboratory drain, rinse with plenty of water.

Precautions

1. Read instructions carefully before conducting the assay. Do not modify the procedure without prior validation.
2. Do not use reagents beyond the expiry date.
3. Do not mix reagents between different lots, as reagents have been calibrated for each batch.
4. Use disposable plasticware wherever possible. Re-usable glassware should be washed thoroughly and rinsed free of detergents before use.

5. A new pipette tip must be used for each pipetting step.
6. Only use laboratory water of the highest quality (e.g. Aqua dest., Aqua bidest., Aqua iniectionabilia, HPLC grade).
7. Do not allow wells to dry out during the assay procedures.
8. Do not expose slides to intense sunlight or similar adverse conditions while incubating.
9. Do not expose the FITC-conjugate at any stage to strong sunlight, UV or fluorescent light. Keep in a dark place whenever possible.

Disposal

Kit reagents, samples and contaminated disposables should be disposed of in accordance with the relevant disposal guidelines and regulations for infectious materials.

Stability and Storage

MASTAFLUOR™ FTA-ABS IgM can be used until the end of expiry date as indicated on kit label. All kits components and reagents should be stored at 2–8 °C and brought to room temperature before use.

After first opening reagents are stable for 90 days. An opened slide shall be used the same day.

Reconstituted PBS powder should be stored at 2–8 °C for up to 30 days.

Samples (serum, plasma) may be stored according to general recommendation in literature. In general samples can be kept at 2–8 °C for up to 3 days prior to use or for longer storage at -20 °C. Repeated freezing and thawing of specimens should be avoided.

Sample Material

Serum / Plasma material is collected by professional personnel according to current standards / best practice guidelines. According to current understanding an infection by contaminated blood cannot be ruled out entirely. Any sample material should therefore be treated as potentially infectious.

Samples should be taken and can be stored at 2–8 °C for up to 3 days according to general recommendation in literature. Serum should be aliquoted immediately after sampling and kept at -20°C for longer storage.

The samples should not be frozen and thawed repeatedly.

After thawing samples should be briefly vortexed carefully before being used in the assay.

Lipemic, haemolytic, icteric or bacterially contaminated samples can cause false positive or false negative results (Ref. 10, 11, 12).

Test Procedure

1. Allow all materials to reach room temperature (at least 20 °C) prior to use.
2. Reconstitute one sachet of PBS powder with 1 litre of distilled water or deionised water.
3. The Sorbent and the Mastorb are ready to use.
4. Prepare samples as follows:
 - 20 µL Sorbent + 60 µL Mastorb, mix carefully.
 - Add 20 µL of undiluted sample (serum or plasma) and mix carefully.

Do not absorb the controls.

5. Incubate all samples at 37 °C for 30 min.
After incubation centrifuge the reaction for 5 min at 3000 rpm in an e.g. Eppendorf centrifuge.
6. Carefully remove the slide from its sachet. Hold the slide at the edge and do not touch the wells. Place the slides in a moist chamber.
7. Apply 20–25 µL of treated specimens and controls to respective wells on the slides according to a prepared work format. Ensure that all the wells are covered, and that serum does not escape from the wells. Direct contact of a pipette tip with the slide surface should be avoided.
8. Place the slides in the moist chamber and incubate at 37 °C for 60 min.
9. After incubation wash slides carefully with PBS from a wash bottle, taking care not to direct the jet onto the test wells. This may be done by directing the jet of PBS along the center of the slide, tilting the slide first towards wells 1–5 and then towards 6–10.
10. Immerse the slides in a staining dish or Coplin jar containing PBS and wash for 15 min at room temperature with a change of PBS after 5 min to increase washing strength.
11. Remove slides from the staining jar and drain off any excess buffer. Using a blotter, dry the area outside wells. Do not touch the well surface.
12. Immediately transfer the slide to a moist chamber and add 20–25 µL of the respective FITC-conjugate to each well. The conjugate shall cover the whole well. Do not allow the wells to dry.
13. Incubate slides at 37 °C in a covered moist chamber in the dark for 30 min.
14. After incubation wash the slides with PBS as in steps 9 and 10.
15. Apply one small drop of Mounting Medium to each well on a slide. Place a cover slip on the slide and read the results immediately. Take care not to trap air pockets under the cover slip.
Note: Remove excess Mounting Medium with a paper towel avoiding any direct movement of the cover slip. Excess Mounting Medium on a slide may result in high background fluorescence due to light scattering.
16. Examine reactions under a fluorescence microscope at a total magnification of 400x or at 800x.

Interpretation of Results

Controls have to fit the lot specific reactions stated in the QC certificate. The figures in the table below are for orientation only and do not reflect the data of the current lot.

Specimen	Expected Results
Positive Control	1 to 4 +
Negative Control	negative*

* negative: Treponemes can be visible but do not fluoresce.

For a correct interpretation the results should be compared with the controls.

If further dilution on the samples or controls are requested based on the screening titer, this can be achieved by serial dilution with PBS on screening titer.

Interpretation of Specimen Results

Positive Samples: Treponema show a green fluorescence. The fluorescence intensity depends on the concentration of the antibodies.

In an IgM positive reaction the outer shape of the Treponema may show a slight diffuse signal, which is a bit hard to focus. This can be seen as a typical IgM reaction.

The intensity can be scored for fluorescence on a scale from 1+ to 4+. For titration a 1+ intensity is counted as minimally reactive.

In general, the background reactivity in an IgM test is slightly stronger than in an IgG test.

Borderline Samples: Displaying a faint but nevertheless perceptible fluorescence (+/-) and should be repeated.

Negative Samples: In most cases no treponemes are visible. Treponemes which are stained green but are not fluorescent must be read NEGATIVE.

In some cases an inhomogeneous pattern is observed in the 1:5 screening dilution. Fluorescence is obtained in some treponemes or in parts of them, other treponemes in the same well are just stained without any fluorescence. These wells should be read as negative.

Limitations / Interferences

1. The FTA-ABS IgM test may occasionally give a false positive result under certain patient conditions or diseases e.g. pregnancy, leprosy and systemic lupus erythematosus (see references 5, 6).
2. Sera with a high anti-Borrelia antibody titer (> 1:640) may show a green fluorescence of Treponema.
3. The FTA-ABS IgM test has been demonstrated to be highly sensitive and specific, however, results should be considered along with all serological tests, clinical history and other aspects of patient management to be considered diagnostically significant.
4. Pathogenic treponemes include *T. pallidum*, the cause of syphilis, *T. pertenue*, the cause of yaws and *T. carateum*, the cause of pinta. The only treponeme of importance in Europe is *T. pallidum*. In addition, many commensal species of treponemes exist and it

is important to differentiate these from *T. pallidum* before a diagnosis of primary syphilis is given.

5. Commensal treponemes can be cultured in artificial media whereas attempts to culture pathogenic treponemes *in vitro* failed.
6. The clinical diagnosis of syphilis is confirmed in the laboratory by either demonstrating the presence of *T. pallidum* in the exudates from the lesions, or demonstrating the presence of serum antibodies against the organism. Methods used to measure antibody response to treponemal infection can be divided into two major categories:
 - a) Tests to measure antibodies against non-specific treponemal antigens i.e. cardiolipin or lipoidal antigen tests. These were formerly called 'Reagin' tests.
 - b) Tests to measure antibodies against antigens specific for pathogenic treponemes i.e. *T. pallidum* particle agglutination assay (TPPA) and the fluorescent treponemal antibody absorbance test (FTA-ABS).

Performance characteristics

Sensitivity and Specificity

The assay was validated on a panel consisting of 354 sera and plasma samples. The panel was derived from positive and negative clinical samples. MASTABLOT™ TP IgM (Mast Diagnostica, Reinfeld, Germany, REF 6653M08) and SERODIA®-TP•PA (Fujirebio Inc., Tokyo, Japan, REF 226414) were used as comparison products.

The analytical performance was calculated as follows:

	Formula	Value
Sensitivity	$\frac{TP}{(TP + FN)}$	95.73 %
Specificity	$\frac{TN}{(TN + FP)}$	91.14 %
Positive predictive value	$\frac{TP}{(TP + FP)}$	84.21 %
Negative predictive value	$\frac{TN}{(TN + FN)}$	97.74 %
Efficiency	$\frac{(TP + TN)}{total}$	92.66 %
Positive likelihood ratio	$\frac{(Sens)}{(1 - Spec)}$	10.80
Negative likelihood ratio	$\frac{(1 - Sens)}{(Spec)}$	0.05

(Abbreviations: TP: True positive, TN: True negative, FN: False negative, FP: False positive)

Trueness

Sensitivity and specificity of MASTAFLUOR™ FTA-ABS IgM is in line with performance data of the comparison products.

Cross-Reactivity

20 *Borrelia* samples were tested for cross-reactivity. Sera with a high anti-*Borrelia* antibody titer sometimes gave false positive results (see Limitations).

Traceability

The INSTAND e. V. RV 311, 2012, Pr. 61 (Serum RVTP 5) sample is an IgM positive sample and can be used to adjust the serum being used as positive control.

Assay Precision


When positive control titration (serial dilution 1:5 to 1:640) was tested 10 times from different lots according to the test procedure, all results were found to be within ± 1 level of fluorescence. Test within the same lot similarly demonstrated repeatability. Thus, assay precision was demonstrated.

High dose hook effect (prozone effect)

The MASTAFLUOR™ FTA-ABS IgM assay is an indirect immunofluorescence detection assays using an FITC-coupled antibody (2-step assay). A washing step has been added in which all unspecific or unbound antibodies are removed. For this reason, the high-dose hook effect is prevented.

Measuring range

When positive control titration (IgM: 1:2-1:256) and sample dilution (1:5) was tested according to the test procedure, following results for measuring range were documented:

Cut-off

The cut-off (1:5 dilution, minimum 1+ fluorescence intensity) was defined on a panel consisting of 354 positive and negative sera and plasma samples corresponding to the mentioned sensitivity and specificity.

Clinical Performance

The clinical performance was demonstrated by a peer-reviewed paper and by proficiency testing.

Müller et al., 2006 (reference 9) showed by using a meta-analysis of external quality control surveys by the German Infection Serology Proficiency Testing Program a mean accuracy of 88% (qualitative) / 65.8% (quantitative) for the fluorescent treponemal antibody absorption test.

For the MASTAFLUOR™ FTA-ABS IgM assay, proficiency tests (INSTAND e.V. Ringversuche, Section 311, *Treponema pallidum*) are performed 1-2 times per year. Data on the interlaboratory comparisons will be provided on request.

Availability of the summary of safety and performance (Art. 29, IVDR)

The report will be available in the EUDAMED database (as soon as the module is available).

(<https://ec.europa.eu/tools/eudamed/#/screen/home>).

The SSP report will be provided on request.

Reporting serious incidents

All serious incidents that have occurred in connection with the device must be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

References and Change History

You can find the references and the Change History at the end of the instruction for use.

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Referenzen / References / Références / Riferimenti:

1. Schöfer H, Brockmeyer NH, Hagedorn HJ, Hamouda O, Handrick W, Krause W, Marcus U et al. (2006): Syphilis. Leitlinie der Deutschen STD Gesellschaft zur Diagnostik und Therapie der Syphilis. J. Dtsch. Dermatol. Ges. 4 (2): 160-177.
2. Rath PM, Marsch WC, Brade V, Fehrenbach FJ (1994): Serological distinction between syphilis and Lyme borreliosis. Zentralbl. Bakteriol. 280 (3).
3. Hans-Jochen Hagedorn, Syphilis; in Lothar Thomas, Labor und Diagnose 8. Auflage (2012), TH-Books; Kapitel 42.14: 2006 - 2017
4. Annual epidemiological report 2014 - Antimicrobial resistance and healthcare-associated infections; Chapter: Syphilis, p. 51-53; <http://ecdc.europa.eu/>
5. Smikle MF, James OB, Prabhakar P (1990): Biological false positive serological tests for syphilis in the Jamaican population. Sexually Transmitted Infections, 66, p. 76-78
6. Wright DJM (1973): The significance of the fluorescent treponemal antibody (FTA-ABS) test in collagen disorders and leprosy Journal of Clinical Pathology, 26, p. 968-972
7. Syphilis, in Surveillance Report, Sexual Transmitted Infections in Europe, 2013, p. 25-31; www.ecdc.europa.eu
8. Bosshard PP (2013): Usefulness of IgM-specific enzyme immunoassays for serodiagnosis of syphilis: Comparative evaluation of three different assays. J. Infec., 67 (1), p. 35-42.
9. Müller I, Brade V, Hagedorn HJ, Straube E, Schörner C, Frosch M et al. (2006): Is Serological Testing a Reliable Tool in Laboratory Diagnosis of Syphilis? Meta-Analysis of Eight External Quality Control Surveys Performed by the German Infection Serology Proficiency Testing Program. Journal of clinical microbiology, 6, 1335-1341.
10. Thomas, L.: Haemolysis as influence & interference factor (2002) eJIFCC 13(4) 95-8.
11. Nikolac, N.: Lipemia: causes, interference mechanisms, detection and management (2014) Biochem 24(1): 57-67
12. Mainali, S. et al: Frequency of icteric interference in clinical chemistry laboratory tests and causes of severe icterus (2021) Pract. Lab. Med. e00259

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Change history includes changes to significant aspects of the assay / IFU.

Chapter	Description of change
Contents	Detailed description of ProClin concentration and link to warnings
Warnings	Update of the warning for ProClin

Mast Diagnostica GmbH,
Feldstraße 20,
23858 Reinfeld,
Deutschland
Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
email: mast@mast-diagnostica.de
Web: www.mast-group.com

Mast Group Ltd.
Mast House, Derby Road
Bootle, Merseyside, L20 1EA
United Kingdom
Tel: +44 (0)151 472 1444
Fax: +44 (0)151 944 1332
email: sales@mastgrp.com
Web: www.mast-group.com

Mast Diagnostic
12 rue Jean-Jacques Mention
CS91106, 80011 Amiens, CEDEX 1
France
Tél: +33 (3) 22 80 80 67
Fax: +33 (3) 22 80 99 22
email: info@mast-diagnostic.fr
Web: www.mast-group.com

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Icons are used according to DIN EN ISO 15223-1

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