eNAT[®]

Instructions for Use

English	3
Italiano	7
Español	11
Deutsch	15
Français	19
Português	23
Česky	27
Dansk	31
Norsk	35
Polski	39
Srpski	43
УКРАЇНСЬКА	173

Copan eNAT® Collection and Preservation System
Instructions for use

INTENDED USE

Copan eNAT® System is intended for the transport and preservation of clinical specimens to be analyzed by nucleic acids amplification techniques. eNAT® medium stabilizes and preserves RNA/DNA for prolonged time periods and is compatible with commercial nucleic acid extraction and amplification platforms.

SUMMARY AND PRINCIPLES

Clinical specimens stored and transported in eNAT® medium can be processed, using standard clinical laboratory operating procedures, for the detection of nucleic acids of Viruses, Bacteria, Chlamydia, Protozoa and Mycoplasma with molecular amplification assays. The primary purpose of nucleic acids amplification techniques is to screen for a wide range of infectious diseases, so nucleic acids integrity of clinical specimens during transport and storage should be preserved⁽¹⁾. eNAT® medium contains a detergent and a protein denaturant to prevent microbial proliferation, thus eNAT® is not intended to be used for culture based techniques.

PRODUCT DESCRIPTION

eNAT® is available in the product configurations indicated in the table below:

CODE	DESCRIPTION
606C	12x80mm tube containing 2 mL of eNAT® transport and preservation medium
608C	12x80mm tube containing 1 mL of eNAT® transport and preservation medium
6E021S	eNAT® Collection kit comprises: - 12x80 mm screw cap tube containing 1 ml of eNAT® transport and preservation medium - one regular applicator FLOQSwabs® - one 2 ml Pasteur pipet

REAGENTS

Guanidine thiocyanate
Tris-EDTA
HEPES
Detergent

REQUIRED MATERIALS THAT ARE NOT INCLUDED

Collection device for formats different from the kit. Appropriate materials for molecular testing according to recommended protocols as per laboratory reference manuals. Urine cup for urine collection.

STORAGE OF THE PRODUCT

This product is ready for use and no further preparation is necessary. The product should be transported and stored in its original container at 5-25°C until use. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual specimen transport vial label.

WARNINGS AND PRECAUTIONS

1. Single-use in vitro diagnostic device for professional use.
2. eNAT® medium is not for external or internal use in humans or animals.
3. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
4. The use of this product in association with diagnostic instrumentation should be previously validated by the user.
5. Do not use if the product is visibly damaged
6. eNAT® contains a cell lysing agent, so a pelleting procedure is not recommended for nucleic acid concentration.
7. Before transporting, make sure eNAT® screw cap tube is tightly closed
8. eNAT® was tested for microbial viability inactivation of Gram positive bacteria, Gram negative bacteria, yeasts and molds. However universal precautions for safe handling of biological fluids should be practiced at all times.
9. Dispose of unused reagents, waste and specimens in accordance with local regulations.
10. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.
11. Avoid contact of eNAT® medium with skin, mucous membranes. If contact does occur, immediately wash with a large amount of water.
12. eNAT® medium contains guanidine thiocyanate. Avoid direct contact between guanidine thiocyanate and sodium hypochlorite (bleach) or other highly reactive reagents such as acids and bases. These mixtures could release noxious gas.
- 13.

Danger
Contains Guanidine thiocyanate.

- H302 Harmful if swallowed
- H314 Causes severe skin burns and eye damage

- H412 Harmful to aquatic life with long lasting effects
- P264 Wash hands thoroughly after handling
- P273 Avoid release to the environment
- P280 Wear protective gloves/protective clothing/eye protection/face protection
- P301+P330+P331 IF SWALLOWED: rinse mouth. Do not induce vomiting
- P303+P361+P353+P310 IF ON SKIN (or hair): remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Immediately call a POISON CENTER or doctor/physician.
- P305+P351+P338+P310 IF IN EYES: rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing. Immediately call a POISON CENTER/doctor
- EUH032 Contact with acids liberates very toxic gas

MSDS available on request from Copan Italia S.p.A. via F. Perotti 10, 25125 Brescia Italy.

PRODUCT DETERIORATION

eNAT® should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed or (4) there are other signs of deterioration.

INSTRUCTIONS FOR USE

Proper specimen collection from the patient is extremely critical for the successful identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published standard collection manuals^(2,3). Qualified medical personal only shall collect clinical specimens using proper sampling devices.

Do not use eNAT® medium for pre-moistening or pre-wetting the sampling device prior collecting the sample or for rinsing or irrigating the sampling site.

1. Open the kit package and remove the test tube. Take the swab out of its pouch and collect the specimen from the patient (Fig. 1). In order to prevent the risk of contamination, make sure that the swab tip comes into contact only with the sampling site. **Note: valid for kit version.**
2. After collecting the specimen, unscrew and remove the cap from eNAT® tube making sure not to spill the medium.
3. Insert the specimen into the tube.

If a swab is used for the collection follow the below procedure.

Suggested procedure for Copan swab:

- **for swab with breaking point:** insert the swab into the tube until the breaking point reaches the level of the opening of the tube (Fig. 2). Bend and break the swab at the breaking point holding the tube away from the face. If needed gently bend the swab shaft up to 180 degrees angle and rotate the swab shaft to complete the breakage (Fig. 3a and Fig.3b). Discard the upper part of the swab shaft.

- **for swab without breaking point:** insert the swab into the tube and cut off the excess part of the shaft.

4. Replace the cap on the tube and close tightly (Fig. 4).
- NOTE:** For swab other than Copan swab different procedure may apply.
5. Write patient information on the tube label or apply patient identification label (Fig. 5).
6. Send the sample to the test laboratory.


Fig. 1


Fig. 2


Fig. 3a


Fig.3b


Fig 4



Fig.5

Urine specimen collection (REF 6E021S)

1. Open the eNAT® sample collection pouch. Remove the eNAT® tube and discard the pouch containing the FLOQSwabs®. Take the Pasteur pipet making sure not to touch the tip. Do not lay the pipet down on a surface.
2. Ask the patient to collect the first 20 to 30 ml of voided urine (the first part of the stream) into a urine cup (not supplied with the kit).
3. Unscrew and remove the cap from eNAT® tube making sure not to spill the medium.
4. Transfer urine from the cup into the collection tube
NOTE: In order to avoid excessive medium dilution, the volume of liquid specimen to be added to eNAT® medium should never exceed 1:3 ratio. Maximum fill volume is 6ml.
 Two suggested procedures:
 1. Transfer 2 ml of urine from the cup into the collection tube by using the Pasteur pipet provided with the eNAT® Collection Kit. The Pasteur pipet has a scale with a mark at every 0.5 ml step. Squeeze the bulb in order to aspirate 2 ml of urine. Make sure to dispose volume exceeding 2 ml back into the cup. Take special care not to introduce contamination in the medium tube.
 2. Transfer 3 ml of urine from the cup into the collection tube by using the Pasteur pipet provided with the eNAT® Collection Kit. The Pasteur pipet has a scale with a mark at every 0.5 ml step. Squeeze the bulb in order to aspirate 2 ml of urine. Dispense 2 ml urine into collection tube. Make sure to dispose volume exceeding 2 ml back into the cup. Repeat squeezing the bulb in order to aspirate 1 ml of urine. Dispense 1 ml urine into collection tube. Make sure to dispose volume exceeding 1 ml back into the cup. Take special care not to introduce contamination in the medium tube.
5. Replace the cap on the tube and close tightly.
 Mix the urine with the transport medium by vortexing the tube 5 seconds.
6. Write patient information on the tube label or apply patient identification label. Send the sample to the test laboratory.

Use in the laboratory

Processing eNAT[®] specimens for molecular testing in the laboratory.

Specimens received in the laboratory for nucleic acid detection should be processed when received in the laboratory. In case of delay, please refer to the appropriate specimen storage conditions. eNAT[®] medium preserves nucleic acids for up to 4 weeks at room temperature and at 4°C⁽¹⁾ and up to 6 months at -20°C to -80°C.

Specimens preserved in eNAT[®] medium may need to be extracted and purified before amplification, depending on the extraction method used. eNAT[®] medium has been tested with automated platforms like, for example, NucliSENS[®] easyMAG[®] (Biomereux), Abbott m2000 system (Abbott Diagnostics), QIAAsymphony (Qiagen)⁽⁴⁾, Microlab NIMBUS (Hamilton)⁽⁵⁾, xTAG[®] Technology (Luminex)⁽⁶⁾ and other manual extraction methods based on silica columns and magnetic beads. Other extraction methods may also be applicable prior validation by the user. For the complete list of the extraction and purification methods tested with eNAT, contact Copan Italia Customer Care (customercare@copangroup.com).

The following steps must be performed:

1. Wear gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations^(7,8,9,10,11).
2. When working with NAAT assays, care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential


VORTEXING METHOD:

3. Vortex eNAT[®] specimen tube for 10 s
NOTE 1: if the sample appears too mucousy, the specimen may remain mainly attached to the swab, extend the vortex time in order to break down mucus clamps and to the sample from the swab.
NOTE 2: If some foam appears in the tube after vortexing, wait a few seconds before opening it.
4. Unscrew the cap and transfer the appropriate amount of sample (e.g. 200ul-400ul or as per the protocol in use for the extraction) directly into the extraction buffer tube. When opening the eNAT[®] tube take care not to spill the medium. **NOTE:** If using automated systems, refer to the instructions for use of the instrument.
5. Continue as per extraction and amplification kits procedures.

If unable to use the vortexing method, use the following alternative protocol:

MANUAL SHAKING METHOD:

1. Hold the eNAT[®] tube from the cap making sure that it is closed tightly.
2. Shake the tube 4 times downward with rapid movements of the wrist (see picture)
NOTE: Inverting the tube up and down is not recommended. If the sample appears too mucousy, the specimen may remain mainly attached to the swab, it is recommended to extend the shaking times in order to break down mucus clamps and to easier release the sample from the swab
3. Unscrew the cap and transfer the appropriate amount of sample (e.g. 200ul-400ul or as per the protocol in use for the extraction) directly into the extraction buffer tube.
NOTE: If using automated systems, refer to the instructions for use of the instrument.
4. Continue as per extraction and amplification kits procedures.

Manual Shaking protocol has been validated using vaginal samples. The use of this protocol with other type of samples, should be previously validated by the user. Other extraction methods may also be applicable prior validation. For updated list of evaluated kits, please contact Copan Italia.

QUALITY CONTROL

ANTIMICROBIAL ACTIVITY

eNAT[®] medium is routinely tested for its antimicrobial activity against a panel of bacteria (*E. coli*, *S. aureus* and *C. albicans*). A complete microbial viability inactivation is obtained within 30 minutes for these three strains, starting from $\geq 10^3$ viable CFU/ml inoculated in 1 mL of eNAT[®] medium.

PERFORMANCE CHARACTERISTICS

NUCLEIC ACIDS PRESERVATION

eNAT[®] medium preserves nucleic acids for up to 4 weeks at room temperature and at 4°C and up to 6 months at -20°C to -80°C. Performance testing with Copan eNAT[®] was conducted using laboratory strains spiked onto a swab. Performance testing was not conducted using human specimens. Tested strains include:

ORGANISM	REF NUMBER	ANALYTE TYPE
INFLUENZA A VIRUS	ATCC [®] VR-822	RNA VIRUS
INFLUENZA B VIRUS	ATCC [®] VR 786	RNA VIRUS
CITOMEGALOVIRUS	ATCC [®] VR-977	DNA VIRUS

HERPES VIRUS TYPE I	ATCC® VR-539	DNA VIRUS
CHLAMYDIA TRACHOMATIS	ATCC® VR-880	CHLAMYDIA
CHLAMYDIA PNEUMONIAE	ATCC® VR-1360	CHLAMYDIA
NEISSERIA GONORRHOEAE	ATCC® 43069	GRAM NEGATIVE BACTERIUM
BORDETELLA PERTUSSIS	ATCC® 8467	GRAM NEGATIVE BACTERIUM
PSEUDOMONAS AERUGINOSA	ATCC® 27853	GRAM NEGATIVE BACTERIUM
MRSA	ATCC® 43300	GRAM POSITIVE BACTERIUM
STAPHYLOCOCCUS AERUS	ATCC® 6538	GRAM POSITIVE BACTERIUM
MYCOPLASMA HOMINIS	ATCC® 23114	MYCOPLASMA
MYCOPLASMA PNEUMONIAE	ATCC® 15331	MYCOPLASMA
TRICHOMONAS VAGINALIS	Biomed culture	PROTOZOAN

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

MICROBIAL VIABILITY INACTIVATION

eNAT® allows to inactivate microbial viability in a short time after inoculum.

Starting from an initial strain concentration of $\geq 10^5$ CFU (or IFU)/ ml, bacteria, yeasts and viruses are completely inactivated in ≤ 30 minutes.

Molds are completely inactivated in ≤ 1 hour.

Viability inactivation testing with Copan eNAT® was conducted using highly concentrated laboratory strains.

Strain	ATCC®	Type	Starting spiked concentration	Growth after 15'	Growth after 30'	Growth after 1 hour
<i>Streptococcus pneumoniae</i>	ATCC® 6305	Gram positive bacteria	$>10^6$ CFU/ml	No growth	No growth	No growth
<i>Streptococcus agalactiae</i>	ATCC® 12386		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Staphylococcus aureus</i>	ATCC® 6538		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Peptostreptococcus anaerobius</i>	ATCC® 27337		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Escherichia coli</i>	ATCC® 25922	Gram negative bacteria	$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Pseudomonas aeruginosa</i>	ATCC® 27853		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Shigella sonnei</i>	ATCC® 9290		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Salmonella typhimurium</i>	ATCC® 14028		$>10^7$ CFU/ml	No growth	No growth	No growth
<i>Haemophilus influenzae</i>	ATCC® 10211		$>10^6$ CFU/ml	No growth	No growth	No growth
<i>Neisseria gonorrhoeae</i>	ATCC® 43069		$>10^6$ CFU/ml	No growth	No growth	No growth
<i>Candida albicans</i>	ATCC® 10231	Fungi	$>10^5$ CFU/ml	No growth	No growth	No growth
<i>Aspergillus brasiliensis</i>	ATCC® 16404		$>10^5$ CFU/ml	140 CFU/100 uL	9 CFU/100 uL	No growth
<i>Chlamydia trachomatis</i>	ATCC® VR-880	Intracellular Gram negative bacteria	$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected
Influenza B	ATCC® VR-786	RNA virus	$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected
Influenza A	ATCC® VR-822		$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected
Herpes simplex virus type 1	ATCC® VR-539	DNA virus	$>10^5$ IFU/ml	No IFU detected	No IFU detected	No IFU detected

Legend: CFU= colony forming unit, IFU= infectious units

Results obtained will largely depend on the correct and adequate sampling as well timely transport and processing in the laboratory.

Please refer to symbol table at the end of the instructions for use

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BIBLIOGRAPHY:

- 1) Short-Term Stability of Pathogen-Specific Nucleic Acid Targets in Clinical Samples. Mohammad R. Hasan, Rusung Tan, Ghada N Al-Rawahi, Éva Thomas, Peter Tilley. J Clin Microbiol. 2012 Dec; 50(12): 4147–4150. Clinical and Laboratory Standards Institute 2005. Collection, transport, preparation, and storage of specimens for molecular methods; approved guideline. MM13-A CLSI, Wayne, PA.
- 2) World Health Organization 2019. Guidance on regulations for the transport of infectious substances 2019–2020.
- 3) Evaluation of Anatomically Designed Flocked Rectal Swabs for Molecular Detection of Enteric Pathogens in Children Admitted to Hospital with Severe Gastroenteritis in Botswana. David M. Goldfarb, Andrew P. Steenhoff, Jeffrey M. Pernica, Sylvia Chong, Kathy Luinstra, Margaret Mokomane, Loeto Mazhani, Isaac Quayle, Irene Goercke, James Mahony, Marek Smiejaa. Journal of Clinical Microbiology, p.3922-3927, Volume 52, November 2014.
- 4) Validation of Copan eNat, a Molecular Transport Medium, For the Collection and Preservation of Urine Specimens For the Detection of STI Infections with the Seegene Anyplex II STI-7 V1.1 Assay, S. Razeti, BMJ Journals, 2013.
- 5) COPAN ENAT MEDIUM STABILIZES NUCLEIC ACIDS FOR THE DETECTION OF VIRUSES WITH THE XTAG FAST RESPIRATORY VIRUS PANEL. Alice Squassina, Arnalda Giambra, Santina Castriciano, Kathy Luinstra, and Marek Smieja. 26th European Congress of Clinical Microbiology and Infectious Diseases 2016.
- 6) J. Michael Miller, Shelley A. Miller, 2017. A Guide to Specimen Management in Clinical Microbiology, Third Edition. ASM, Washington DC.
- 7) Centers for Disease Control and Prevention (CDC), 2016. Guide for Shipping Infectious Substances.
- 8) 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. October 1, 2007. Interstate Shipment of Etiologic Agents.
- 9) Centers for Disease Control and Prevention (CDC), 2009. Biosafety in Microbiological and Biomedical Laboratories 5th Edition. Ex 6
- 10) Hansen DJ. Healthcare, Laboratories and Biosafety. Vol 2., 1992. CRC Press. Ex 7.
- 11) NHS. UK Standards for Microbiology Investigations- Good practice when performing molecular amplification assays. Quality Guidance | Q 4 | Issue no: 5 | Issue date: 19.02.18.

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Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com