UriSponge™

Instructions for Use

Copan Urisponge™ - Urine Collection, Transport and Preservation System Instructions for Use

INTENDED USE

Copan Urisponge™ - Urine Collection, Transport and Preservation System is intended for the collection, transport and preservation of urine specimens from the collection site to the testing laboratory. In the laboratory, Urisponge™ specimens are processed using standard clinical laboratory operating procedures for the cultivation of uropathogenic bacteria and yeasts.

SUMMARY AND PRINCIPLES

Bacteria quantification of clean-voided midstream collected urine is widely used to evaluate urinary tract infections¹. The maintenance of the bacterial load in urine samples depends on many factors including the type and concentration of the microorganisms, duration of transport and storage temperature². The Urisponge™ device incorporates preservative substances (urine maintenance formula) onto the applicator sponge within a conical bottom, screw-cap tube, which sustains the viability of clinically important organisms during transport for up to 48 hours at 2-25 °C.

PRESERVATIVE

Copan Urisponge™ incorporates preservative substances (urine maintenance formula) onto the applicator sponge. The concentration of preservative per milliliter of urine sample collected with Urisponge™ device is estimated as follow:

Boric Acid	7.3 mg/ml
Sodium Formate	2.4 mg/ml

PRECAUTIONS

1. For In Vitro Diagnostic Use only.
2. For professional use only. Use Urisponge™ in accordance with the Package Insert.
3. Do not use if the package or tube is open, there is evidence of damage, deterioration, contamination or the expiration date has passed.
4. Do not re-sterilize. Do not re-pack.
5. Not suitable for the recovery of fastidious microorganisms, e.g. anaerobes, viruses, chlamydiae, mycoplasmas, ureaplasmas trichomonas, in culture.
6. Not suitable for collection of specimens for microscopic and macroscopic examination.
7. Use of this device in association with diagnostic kits and/or instrumentation should be validated prior to use.
8. Wear protective gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe appropriate CDC Biosafety recommendations. After use, tubes must be disposed of according to laboratory regulations for infectious waste^{7,8,9}.
9. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.

DEVICE STORAGE

Store in the original packaging at 2–25 °C until used. Do not freeze prior to use.

DEVICE DESCRIPTION / HOW SUPPLIED

Copan Urisponge™ - Collection, Transport and Preservation system is ready-to-use. Product descriptions and packaging configurations are listed in Table 1.

The .A codes are compatible with the automated device.

Figure 1: Urisponge™ device


Figure 2: Urisponge™ components

Table 1. Urisponge™ Description and Packaging Configurations

PRODUCT CODE	PRODUCT DESCRIPTION	PACK SIZE	URINE SAMPLE VOLUMES
802C	Sterile plastic sealed bag containing 100 pieces of single use devices composed of polypropylene mini size rigid plastic tube screw cap. Inside the tube the cap holds a plastic stick with two cylinders of hydrophilic polyurethane foam incorporating preservative substances.	100 devices per shelf pack 100 x 5 devices per box	Absorbed urine volume from saturated sponges will be approximately 1,0 -1,5 ml.
802CE.A	Sterile plastic sealed bag containing 50 pieces of single use devices composed of polypropylene mini size rigid plastic tube screw cap. Inside the tube the cap holds a plastic stick with two cylinders of hydrophilic polyurethane foam incorporating preservative substances.	50 devices per shelf pack 50 x 6 devices per box	Absorbed urine volume from saturated sponges will be approximately 1,0 -1,5 ml.
802CS01	Sterile single use sample collection pack containing polypropylene mini size rigid plastic tube screw cap. Inside the tube the cap holds a plastic stick with two cylinders of hydrophilic polyurethane foam incorporating preservative substances.	50 devices per shelf pack 10 x 50 devices per box	Absorbed urine volume from saturated sponges will be approximately 1,0 -1,5 ml.
803CE.A	Sterile plastic sealed bag containing 100 pieces of single use devices composed of polypropylene regular size rigid plastic tube screw cap. Inside the tube the cap holds a plastic stick with three cylinders of hydrophilic polyurethane foam incorporating preservative substances.	100 devices per shelf pack 100 x 5 devices per box	Absorbed urine volume from saturated sponges will be approximately 2,5 -3,0 ml.
800C	Sterile plastic sealed bag containing 100 pieces of single use devices composed of polyethylene regular size squeezable plastic tube screw cap. Inside the tube the cap holds a plastic stick with two cylinders of hydrophilic polyurethane foam incorporating preservative substances.	100 devices per shelf pack 100 x 5 devices per box	Absorbed urine volume from saturated sponges will be approximately 2,0 -2,5 ml.
800CS01	Sterile single use sample collection pack containing polyethylene regular size squeezable plastic tube screw cap. Inside the tube the cap holds a plastic stick with two cylinders of hydrophilic polyurethane foam incorporating preservative substances	50 devices per shelf pack 10 x 50 devices per box	Absorbed urine volume from saturated sponges will be approximately 2,0 -2,5 ml.

MATERIALS REQUIRED BUT NOT SUPPLIED

Sterile urine collection cup and appropriate materials for the isolation and culture of uropathogenic bacteria and yeast, e.g. culture media and incubation systems are not provided. Refer to laboratory reference manuals for recommended procedures for the culture and identification of uropathogenic bacteria and yeasts from clinical samples^{1,4}.


INSTRUCTIONS FOR USE

Specimen Collection and Transport with Urisponge™

Proper specimen collection from the patient is critical for successful isolation and identification of infectious organisms. Specimens obtained for assessment of urinary tract pathogens should be collected and handled following published manuals and guidelines^{1,2,3,4}.

1. Obtain a clean-catch urine sample from the midstream portion into a sterile container.
2. Open the Urisponge™ tube and, holding by the cap, dip the sponge applicator into the urine sample. Submerge the sponges for 5 seconds. The polyurethane sponges are extremely hydrophilic and spontaneously absorb urine.
3. Remove the sponge applicator from the urine sample and return it to the Urisponge™ tube. See Figure 3. Do not add urine to the tube other than the urine spontaneously absorbed by the sponge.
4. Screw the cap to securely close the container.
5. Label and transport Urisponge™ specimens in accordance with institutional, local, state, and federal requirements⁶.
6. Immediately transfer Urisponge™ specimen(s) to the laboratory, preferably within 2 hours of collection². If immediate delivery or processing is delayed Urisponge™ specimens may be refrigerated at 2-8 °C or stored at room temperature (20-25 °C) and processed within 48 hours.

Figure 3: Transfer of urine sample to Urisponge™ device.

Testing Urisponge™ Specimens for Culture in the Laboratory

1. Urisponge™ specimens should be processed immediately upon receipt and within 48 hours of collection.
2. Urisponge™ tubes may be centrifuged to extract the urine from the sponge. Centrifuge at 104 g for the mini size tube and 149 g for the regular size tube for 3 seconds in a swing-out rotor centrifuge. Urisponge™ may be re-centrifuged as needed. Alternatively, the specimen can be manually extracted from the sponge by shaking the device down three times using quick, sharp, downward wrist motions.
3. Using aseptic technique, unscrew the cap and remove the sponge applicator. Gently mix the content of the tube.
4. Follow laboratory internal Standard Operating Procedures (SOPs) for inoculating the urine sample from the tube onto culture media or refer to published microbiology manuals and guidelines^{1,4}.

LIMITATIONS

- Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection^{1,2,3,4}.
- The microbial load in urine from a given patient may be influenced by the time of collection and fluid intake. Symptomatic patients may have counts below 10⁵ microorganism/ml if specimens are collected late in the day or if diuresis is occurring¹⁰.
- The sponge preservatives will not inactivate antibiotics.

QUALITY CONTROL

Urisponge™ is verified to maintain viability after storage at 2-8 °C (refrigerated) and 20-25 °C (room temperature) for up to 48 hours using the organisms recommended by the Clinical and Laboratory Standards Institute (CLSI) M40-A2⁵.

Procedures for quality control of microbiology transport devices should be conducted using testing methods described in CLSI M40-A2⁵. If aberrant quality control results are noted, patient results should not be reported.

PERFORMANCE CHARACTERISTICS

The performance of Urisponge™ has been verified per CLSI M40-A2 for "Urine Culture Containers With Preservative" using a known concentration of each organism suspended in filtered synthetic urine. All organisms were ATCC® derived and were maintained in accordance with the manufacturer's instructions.

Microorganism recovery (viability) study

The ability of Urisponge™ to maintain viability of urinary tract microorganisms was assessed after storage at refrigerated and controlled room temperature for 24 and 48 hours. Recovery of each organism on non-selective agar media plates after being held under the challenge conditions was verified to remain within 1 log₁₀ of the initial concentration. See Table 2.

Table 2. Performance of Urisponge™ in maintaining the viability of urinary tract microorganisms

		Average CFU recovered N=9 (from three lots)			
Organism	Holding Temperature	T=0	T=24	T=48	T=48 hours Log reduction (-) or Log increase (+)
Recovery of Urine Culture Organisms in Urisponge™ per CLSI M40-A2					
<i>C. albicans</i> (ATCC® 24433)	2-6 °C	7,37 x 10 ²	4,70 x 10 ²	3,92 x 10 ²	-0,27
	20-24 °C	7,37 x 10 ²	3,79 x 10 ²	3,11 x 10 ²	-0,37
<i>E. coli</i> (ATCC® 25922)	2-6 °C	2,15 x 10 ³	1,78 x 10 ³	1,16 x 10 ³	-0,28
	20-24 °C	2,15 x 10 ³	1,22 x 10 ³	4,26 x 10 ²	-0,71
<i>E. faecalis</i> (ATCC® 29212)	2-6 °C	2,11 x 10 ³	1,64 x 10 ³	1,72 x 10 ³	-0,09
	20-24 °C	2,11 x 10 ³	1,18 x 10 ³	1,04 x 10 ³	-0,31
<i>P. aeruginosa</i> (ATCC® 27853)	2-6 °C	2,40 x 10 ³	1,96 x 10 ³	1,70 x 10 ³	-0,15
	20-24 °C	2,40 x 10 ³	1,34 x 10 ³	5,31 x 10 ²	-0,66
<i>P. mirabilis</i> (ATCC® 7002)	2-6 °C	4,33 x 10 ²	3,76 x 10 ²	3,23 x 10 ²	-0,13
	20-24 °C	4,33 x 10 ²	2,96 x 10 ²	1,47 x 10 ²	-0,47
<i>S. saprophyticus</i> (ATCC® 15305)	2-6 °C	9,67 x 10 ²	6,11 x 10 ²	5,02 x 10 ²	-0,28
	20-24 °C	9,67 x 10 ²	5,54 x 10 ²	5,11 x 10 ²	-0,28
Recovery of Additional Relevant Pathogenic Organisms of the Urinary Tract in Urisponge™					
		Average CFU recovered N=9 (from three lots)			
Organism	Holding Temperature	T=0	T=24	T=48	T=48 hours Log reduction (-) or Log increase (+)
<i>C. freundii</i> (ATCC® 8090)	2-6 °C	2,36 x 10 ³	2,06 x 10 ³	1,76 x 10 ³	-0,13
	20-24 °C	2,36 x 10 ³	1,75 x 10 ³	1,55 x 10 ³	-0,19
<i>C. glabrata</i> (ATCC® MYA-2950)	2-6 °C	1,52 x 10 ³	7,78 x 10 ²	6,21 x 10 ²	-0,39
	20-24 °C	1,52 x 10 ³	5,38 x 10 ²	3,23 x 10 ²	-0,67

<i>E. cloacae</i> (ATCC® 13047)	2-6 °C	1,01 x 10 ³	8,53 x 10 ²	6,71 x 10 ²	-0,29
	20-24 °C	1,01 x 10 ³	6,83 x 10 ²	1,40 x 10 ³	0,15
<i>M. morganii</i> (ATCC® 25829)	2-6 °C	2,27 x 10 ³	2,19 x 10 ³	1,38 x 10 ³	-0,22
	20-24 °C	2,27 x 10 ³	1,61 x 10 ³	1,40 x 10 ³	-0,21
<i>S. agalactiae</i> (ATCC® 12386)	2-6 °C	2,26 x 10 ³	1,85 x 10 ³	1,23 x 10 ³	-0,29
	20-24 °C	2,26 x 10 ³	1,15 x 10 ³	8,39 x 10 ²	-0,45

Microorganism release study (centrifugation and manual release methods)

The ability of UriSponge™ to release organisms suspended in filtered synthetic urine at known concentrations using centrifugation and manual method was compared with the concentration of the organisms in the initial spiked synthetic urine. See Table 3.

Table 3. Performance of UriSponge™ in releasing a select set of uropathogenic microorganisms using the centrifuge and manual release methods

Organism	Average CFU/mL Recovered: Spiked Urine	Average CFU/mL Recovered: Time 0 hrs Centrifugation Release	Average CFU/mL Recovered: Time 0 hrs Manual Release	Centrifugation Release vs. Spiked Urine Log reduction (-) or increase (+)	Manual Release vs. Spiked Urine Log reduction (-) or increase (+)
<i>C. albicans</i> (ATCC® 24433)	6.17 x10 ²	5.38 x10 ²	5.47 x10 ²	-0.06	-0.05
<i>E. coli</i> (ATCC® 25922)	1.39 x10 ³	1.49 x10 ³	1.46 x10 ³	0.03	0.02
<i>S. saprophyticus</i> (ATCC® 15305)	1.19 x10 ³	1.19 x10 ^{3a}	1.22 x10 ³	-0.002	0.01

^aCFU/mL for this sample was 1.186 x10³ before rounding to two decimal places.

Fill volume study

The ability of UriSponge™ to maintain the viability of organisms was evaluated in case of non-saturated sponges, verifying the lack of toxicity effect on the urine flora. See Table 4.

Table 4. Results of UriSponge™ fill volume study

Organism	Average CFU/mL Recovered: Time 0 hrs	Average CFU/mL Recovered: Time 48 hrs at 2-6 °C	T = 48 hours Log reduction (-) or increase (+)
<i>E. coli</i> (ATCC® 25922)	1.73 x10 ⁴	2.67 x10 ³	-0.81
	1.90 x10 ⁴	3.63 x10 ³	-0.72
	2.17 x10 ⁴	2.35 x10 ³	-0.97
<i>P. aeruginosa</i> (ATCC® 27853)	1.22 x10 ⁴	1.69 x10 ³	-0.86
	9.70 x10 ³	1.35 x10 ³	-0.86
	1.02 x10 ⁴	1.33 x10 ³	-0.89
<i>P. mirabilis</i> (ATCC® 7002)	1.96 x10 ⁴	2.97 x10 ³	-0.82
	2.13 x10 ⁴	2.48 x10 ³	-0.93
	1.50 x10 ⁴	2.80 x10 ³	-0.73

Please refer to symbol table at the end of the instructions for use.

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BIBLIOGRAPHY

1. Yvette S. Mccarter, Eileen M. Burd, Gerri S. Hall, And Marcus Zervos Laboratory Diagnosis of Urinary Tract Infections. Cumitech 2C, 2009. ASM, Washington DC.
2. Miller JM. A Guide to Specimen Management in Clinical Microbiology. Third Edition. American Society for Microbiology. Washington, DC. 2017.
3. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 2001. Urinalysis and Collection, Transportation, and Preservation of Urine Specimens, Approved Guideline. Second Edition. GP16-A2. Vol.21 No.19.
4. Lynne S. Garcia, Clinical Microbiology Procedures Handbook, Third Edition ASM, Washington DC 2010.
5. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 2014. Quality Control of Microbiological Transport Systems; Approved Standard. Second Edition M40-A2 Vol. 34 No. 9.
6. 542CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
7. Fleming D. Biological Safety: Principles and Practices. January 2000. ASM, Washington DC.
8. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA.
<http://www.cdc.gov/od/ohs/symp5/jyrtext.htm>.
9. Richardson JH. Biosafety in Microbiological and Biomedical Laboratories. December 1994. Diane Publishing Company.
10. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957; 100:700-714.

INDEX OF SYMBOLS

Symbol / Simbolo / Simbolo / Symbol / Symbole / Símbolos / Símbol / Символ /	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Značenje / Anlamı / Значения
	Manufacturer / Fabbricante / Fabricante / Hersteller / Fabricant / Fabricante / Proizvođač / Üretici / Виробник
	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro / Dispositivo de diagnóstico in vitro / In vitro dijagnostički uređaj / för in vitro-diagnostik / Медичний вибір для діагностики in vitro
	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE / Marcação CE / CE oznaka / CE işareti / Маркування CE
	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiação ionizante / Sterilisanu upotrebom jonizujućeg zračenja / İşinlama kullanılan yöntem / Стерилізований завдяки опроміненню
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Nemojte ponovo koristiti / Tekrar kullanmayın / Не використовувати повторно
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Kataloški broj / Katalog numarası / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Limites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Opseg temperature / Sıcaklık limitleri / Температурні обмеження
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Upotrebiti do / Son kullanma tarihi / Використати до
	Consult the operating instructions supplied with the device or available in electronic format, and which can be identified by the e-IFU indicator on the packaging label / Consultare le istruzioni per l'uso fornite con il dispositivo oppure disponibili in formato elettronico ed identificate dall'e-IFU indicator sull'etichetta imballo / Consultar las instrucciones de uso suministradas con el dispositivo o disponibles en formato electrónico e identificadas por el indicador e-IFU de la etiqueta del embalaje / Konsultieren Sie die Gebrauchsanweisung, die mit dem Gerät geliefert wird oder in elektronischem Format vorliegt und durch den e-IFU-Indikator auf dem Verpackungsetikett gekennzeichnet ist / Voir le mode d'emploi fourni avec l'appareil ou disponible au format électronique et identifiable grâce à l'indicateur e-IFU sur l'étiquette de l'emballage / Ver as instruções de utilização fornecidas com o dispositivo ou disponíveis em formato eletrônico e identificadas pelo indicador e-IFU na etiqueta da embalagem / Pogledajte uputstva za upotrebu koja su isporučena sa uređajem ili su dostupna u elektronskom formatu i označena indikatorom „e-IFU“ na etiketi pakovanja / Cihaz ile tedarik edilen veya elektronik şekilde ve ambalaj etiketinde e-IFU kodu ile tespit edilen kullanım talimatlarına bakınız / Див. друковані інструкції з використання або інструкції з використання в електронному форматі, якщо наявна позначка "eIFU Indicator"
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Serijski (lot) broj / Lot (seri) kodu / Koç serisi (naprıtı)
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Sadržaj dovoljan za testova / <n> test için yeterli içerik / Придатний для проведення <n> кількості випробувань

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562 USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com