

Tel: + 44 (0) 151 472 1444 email: sales@mast-group.com Web: www.mast-group.com

Mast Diagnostica GmbH

Feldstrasse 20 DE-23858 Reinfeld Germany

Tel: + 49 (0) 4533 2007 0 Fax: + 49 (0) 4533 2007 68 email: mast@mast-diagnostica.de Web: www.mast-group.com

Mast Diagnostic

12 rue Jean-Jacques Mention CS91106, 80011 Amiens, CEDEX 1 France

Té1: +33 (0) 3 22 80 80 67 Fax: +33 (0) 3 22 80 99 22 email: info@mast-diagnostic.fr Web: www.mast-group.com



MAST® *ID* Sorbitol Agar

IDM9/A/NCE

Intended use

Sorbitol Agar dehydrated culture medium powder is intended for use to produce a medium to aid differentiation of members of the Enterobacterales family based on their ability to ferment the specific carbohydrates incorporated into the basal medium. When prepared in accordance with the instructions for use, it produces a semi-solid medium capable of supporting the growth of Enterobacterales. The metabolic action of the organism able to ferment the carbohydrate substrate results in acidification of the medium, detected by a pH indicator dye which is interpreted as a positive result. This medium can also be used in conjunction with additional identification products to produce a phenotypic biochemical profile of the bacterial isolate in order to generate specific fingerprints or datasets that can be used, for example, to detect or rule out cross-infection or elucidate bacterial transmission.

Sorbitol Agar is intended to be used in conjunction with other phenotypic tests to aid epidemiological typing of previously isolated and identified pure cultures of members of the Enterobacterales family derived from, animal, food, environmental or human samples. It is a non-automated, qualitative device, intended to be used by professional, trained laboratory users for in vitro use and is not intended for use in the diagnosis of disease or other conditions in humans or as the basis of treatment or case management decisions.

Contents

See pack label.

Formulation*

Material:	Concentration in medium:
Peptone mixture	10.0g/litre
Sorbitol	10.0g/litre
Bromothymol Blue	0.06g/litre
Agar	20.0g/litre
Final pH:7.2 \pm 0.2	

Storage and shelf life

All dehydrated culture media containers should be kept tightly closed and stored in a dry place at 10 to 25°C until the expiry date shown on the pack label.

Precautions

For in vitro use only. Observe approved hazard precautions and aseptic techniques. To be used only by adequately trained and qualified laboratory personnel. Sterilise all biohazard waste before disposal. Refer to Product Safety Data sheet (available on request or via MAST® website). Sorbitol Agar is not intended for use in the diagnosis of disease or other conditions in humans.

Materials required but not provided

Standard microbiological supplies and equipment such as loops, swabs, applicator sticks, incinerators and incubators, etc.

Procedure

- Refer to pack label for quantities and volumes required. Prepare MAST® /D Sorbitol Agar (IDM9/A) by suspending the powder in distilled or deionised water. For sachet packs, dissolve the entire contents of the sachet in the volume shown on the label.
- Sterilise by autoclaving at 121°C (15 p.s.i.) for 15 minutes. Do not overheat any carbohydratecontaining medium.
- Mix well and pour culture plates (15 to 20ml per plate) into Petri dishes which have been labelled using the self-adhesive labels provided. Self adhesive labels are provided in each box of preweighed sachets.
- 4. Prepared medium may be used immediately or stored in plastic bags at 2 to 8°C for up to one week before
- 5. Prepare a suspension of each organism equivalent in density to a 0.5 McFarland standard. Inoculate the surface of a well dried plate using a replicating device, e.g. the SCANURIDOT Multipoint Inoculator, to deliver each inoculum onto the agar surface.
- 6. Allow the inoculum drops to dry before disturbing and incubate plates aerobically for 18 to 24 hours at 35 to 37°C.

Interpretation of results

After incubation record growth and colour development in the medium. A positive result indicating fermentation of the carbohydrate is shown by a colour change to yellow. A negative result is no change from the original blue/green.

Quality control

Check for signs of deterioration. Quality control must be performed with at least one organism to demonstrate expected performance. Do not use the product if the result with the control organism is incorrect. The list below illustrates a range of performance control strains which the end user can easily obtain.

Test Organisms		
Escherichia coli	Positive	
ATCC® 25922		
Salmonella typhimurium	Positive	
ATCC® 14028		
Shigella sonnei	Negative	
ATCC® 25931		

References

Bibliography available on request.