eSwab<sup>®</sup>

## Instructions for Use

CE 0123

## Copan Liquid Amies Elution Swab (ESwab®) Collection and Preservation System

### Instructions for use

#### INTENDED USE

Copan Liquid Amies Elution Swab (ESwab®) Collection and Transport System is intended for the collection and transport of clinical specimens containing aerobes, anaerobes, fastidious bacteria, viruses and Chlamydia from the collection site to the testing laboratory.

ESwab® medium preserves the viability of aerobes, anaerobes, fastidious bacteria from swab specimens for bacterial culture purposes and can be used for the preservation of bacterial, viral or Chlamydial antigens and nucleic acids from swab specimens.

#### SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of bacteriological infections involves the collection and safe transportation of swab samples. This can be accomplished using the Copan Liquid Amies Elution Swab (ESwab®) Collection and Transport System. Copan ESwab® incorporates a modified Liquid Amies transporting medium, which can sustain the viability of a plurality of organisms that include clinically important aerobes, anaerobes and fastidious bacteria such as *Neisseria gonorrhoeae*. ESwab® can also be used for transport and preservation of nucleic acids and bacterial, viral and Chlamydia antigens for molecular amplification techniques and rapid lateral flow immunochromatographic methods. ESwab® is free of enzymes and inhibiting effects that may interfere with molecular amplification assays and rapid lateral flow immunochromatographic methods.

The ESwab® transport medium is a maintenance medium comprising inorganic phosphate buffer, calcium and magnesium salts, and sodium chloride with a reduced environment due to the presence of sodium thioglycollate <sup>(1)</sup>.

Copan ESwab® consists of a sterile package containing two components: a pre-labeled polypropylene screw-cap tube with conical or round shaped bottom filled with 1 ml of Liquid Amies transport medium and a specimen collection swab which has a tip flocced with soft nylon fiber. Three types of applicator shafts are available: regular size flocced nylon applicator intended for the collection of samples from nose, throat, vagina, rectum, faeces or wounds; minitip size flocced nylon applicator intended for the collection of samples from small or less accessible areas such as the eye, ear, nasal passages, throat and urogenital tract; pernasal nylon flocced applicator intended for the collection of samples from the nasopharynx and pediatric sample collection.

In order to improve the sample collection efficiency and minimize patient's discomfort, two specific applicators have been additionally designed: the former for urethral sample collection and the latter for pediatric sample collection. For more details please refer to table 1. Once a swab sample is collected, it should be placed immediately into the ESwab® transport tube where it comes into contact with the transport medium. Swab specimens for bacterial culture, collected using ESwab®, should be transported directly to the laboratory, preferably within 2 hours of collection <sup>(2, 3, and 4)</sup> to maintain optimum organism viability. If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours except for *Neisseria gonorrhoeae* cultures which should be processed within 24 hours. Independent scientific studies on swab transport systems have shown that for certain bacteria viability is superior at refrigerated temperatures compared with room temperature <sup>(12 – 16)</sup>. Swab specimens for bacterial, viral and chlamydial antigens and nucleic acids investigations should be processed before five days when stored at room temperature (20 – 25°C); within 7 days if stored at 4°C and within 6 months when stored at -20°C. When the sample is processed both for bacterial culture and antigens/nucleic acids investigations, the appropriate time and temperature conditions during transport and storage must be considered, as above indicated.

#### REAGENTS

Copan ESwab® incorporates a modified Liquid Amies medium.

#### ESwab® MEDIUM FORMULATION

Sodium chloride  
Potassium chloride  
Calcium chloride  
Magnesium chloride  
Monopotassium phosphate  
Disodium phosphate  
Sodium thioglycollate  
Distilled water

#### TECHNICAL NOTE

The modified Liquid Amies Medium in ESwab® transport tubes can have a cloudy appearance. This is normal and is due to the presence of salts in the medium formulation.

The swab tip may naturally appear from slightly yellow to brownish color, this well-known phenomenon has no impact on product safety and performance.

ESwab® formula contains Sodium Thioglycollate, an important component for the performance of the product and the maintenance of organism viability. Sodium Thioglycollate has a natural sulfur-like odor. It may be possible to detect this odor momentarily when first opening the ESwab® peel pouch. This odor is a perfectly normal and completely harmless characteristic.

#### PRECAUTIONS

1. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
2. All specimens and materials used to process them should be considered potentially infectious and handled in a manner which prevents infection of laboratory personnel. Sterilize all biohazard waste including specimens, containers and media after their use. Observe other CDC Biosafety Level 2 recommendations <sup>(34, 35, 36, and 37)</sup>.
3. Directions should be read and followed carefully.

#### STORAGE

This product is ready for use and no further preparation is necessary. The product should be stored in its original container at 5 – 25°C until used. Do not overheat. Do not incubate or freeze prior to use. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual sterile collection pouch and the specimen transport tube label.

#### SPECIMEN COLLECTION, STORAGE AND TRANSPORTATION

Specimens collected for bacteriological investigations which comprise the isolation of aerobes, anaerobes and fastidious bacteria such as *Neisseria gonorrhoeae* should be collected and handled following published manuals and guidelines <sup>(2, 3, 18, 19, 20, 21, 22, and 23)</sup>. To maintain optimum organism viability, antigen and nucleic acids integrity, transport specimens collected using ESwab® directly to the laboratory, preferably within 2 hours of collection <sup>(2, 3, and 4)</sup>.

If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours except for *Neisseria gonorrhoeae* cultures which should be processed within 24 hours. Swab specimens for bacterial, viral and chlamydial antigens and nucleic acids investigations should be processed before five days when stored at room temperature (20 – 25°C); within 7 days if stored at 4°C and within 6 months when stored at -20°C. When the sample is processed both for bacterial culture and antigens/nucleic acids investigations, the appropriate time and temperature conditions during transport and storage must be considered, as above indicated. Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations <sup>(19, 22, and 23)</sup>. Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

#### MATERIALS SUPPLIED

Fifty (50) units are contained in a "Vi-Pak" shelf pack and 8 x 50 or 6 x 50 units are contained in a box. Each unit consists of a sterile package containing two components: a pre-labeled screw-cap tube with conical or round shaped bottom filled with 1 ml of Liquid Amies transport medium and a specimen collection swab which has a tip flocced with soft nylon fiber (see Fig 1). For more details on available configurations please refer to Table 1. The collection swab applicators provided with ESwab® can have a molded breakpoint in the shaft of the applicator which is highlighted with a colored indication line marked on the shaft of the applicator.

After the sample is collected from the patient, the molded breakpoint facilitates easy breakage of the swab applicator into the eSwab® tube of transport medium. eSwab® tube capture caps have an internal moulded design that is able to capture the swab shaft when it is broken off into the tube and the cap is closed. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a moulded docking receptacle in the cap (Fig.2). In the testing laboratory when the swab cap is unscrewed and removed, the swab applicator is attached to the cap. This feature allows the operator to conveniently remove the swab from the transport tube and perform microbiological analyses using the tube cap as a handle to hold the swab applicator. The capture cap feature is not applicable to pernasal and pediatric swabs (see Table 1 for capture cap feature applicability).

Fig 1 eSwab® Components


Fig 2. Capture of broken swab applicator stick by eSwab® tube cap

MATERIALS REQUIRED BUT NOT SUPPLIED

Appropriate materials for isolating and culturing aerobes, anaerobes and fastidious bacteria. Appropriate materials for bacterial, viral and Chlamydia rapid antigens and nucleic acids extractions and amplifications.

DIRECTIONS FOR USE

Copan eSwab® Collection and Transport System is available in product configurations indicated in the table below.

Table 1

Catalog No.	Copan eSwab® Product Descriptions	Pack Size	Sampling Sites*	Capture Cap Feature
480CE 490CE.A	Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, vagina, rectum, faeces and wounds	YES
480CESR 490CESR.A	Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One regular size applicator swab with flocked nylon fiber tip. Product suitable for surgical room use since it is double wrapped	50 units per "Vi-Pak" shelf pack 6x50 units per box	Nose, throat, vagina, rectum, faeces and wounds	YES
481CE 491CE.A	Sterile single use sample collection pack containing: - Orange Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One minitip applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Eye, ear, nasal passages, throat, urogenital tracts.	YES
482CE	Sterile single use sample collection pack containing: - Blue Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One pernasal applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nasopharynx and pediatric sample collection.	NO
483CE	Sterile single use sample collection pack containing: - Orange Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One urethral applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Urogenital tract.	YES
484CE	Sterile single use sample collection pack containing: - Blue Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One pediatric applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Pediatric sample collection.	NO
490CE02	Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - Two regular size applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, vagina, rectum, faeces and wounds	YES

492CE03	Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One regular size applicator swab with flocked nylon fiber tip. - One pink urethral applicator swab with flocked nylon fiber tip. - One green perinasal applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Throat, vagina, wounds, rectum, faeces or urogenital tracts or Nasopharynx and pediatric sample collection	YES, only for regular size and pink swabs
493CE02	ESwab® MRSA collection system. Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One pink regular size flock swab plus one white regular size flocked swab.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, perineum	YES
493CE03	ESwab® MRSA collection system. Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - Two pink regular size flock swabs plus one white regular size flocked swab.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, perineum	YES
4E011S.A	Sterile single use sample collection pack containing: - White Polypropylene screw-cap tube with round shape bottom filled with 1ml of Liquid Amies Medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, vagina, rectum, faeces and wounds	YES
4E014S.A	Sterile single use sample collection pack containing: - Orange Polypropylene screw-cap tube with round shape bottom filled with 1ml of Liquid Amies Medium. - One minitip applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Eye, ear, nasal passages, throat, urogenital tracts.	YES
4E033S.A	Sterile single use sample collection pack containing: - Orange Polypropylene screw-cap tube with round shape bottom filled with 1ml of Liquid Amies Medium. - One pernasal applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nasopharynx and pediatric sample collection	NO
4E039S	Sterile single use sample collection pack containing: - White Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Throat, vagina, rectum, faeces and wounds	YES
4E053S.A	Sterile single use sample collection pack containing: - White Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One regular size applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, vagina, rectum, faeces and wounds	YES
4E047S02	ESwab® MRSA collection system. Sterile single use sample collection pack containing: - White Polypropylene screw-cap tube with internal conical shape filled with 1ml of Liquid Amies Medium. - One pink regular size flock swab plus one white regular size flocked swab.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, perineum	YES
4E074S.A	Sterile single use sample collection pack containing: - Blue Polypropylene screw-cap tube with round shape bottom filled with 1ml of Liquid Amies Medium. - One pernasal applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nasopharynx and pediatric sample collection	NO
4E068S.A	Sterile single use sample collection pack containing: - Blue Polypropylene screw-cap tube with round shape bottom filled with 1ml of Liquid Amies Medium. - One pediatric applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Pediatric sample collection	NO
4E135S	Sterile single use sample collection pack containing: - Pink Polypropylene screw-cap tube with round shape bottom filled with 1ml of Liquid Amies Medium. One regular applicator swab with flocked nylon fiber tip.	50 units per "Vi-Pak" shelf pack 8x50 units per box	Nose, throat, vagina, rectum, faeces and wounds	YES

Other product codes may be available. For updates please refer to our website: [www.copangroup.com](http://www.copangroup.com)

¶ This is just a suggested table. Performance testing with Copan ESwab® was conducted using laboratory strains spiked onto a swab following the test protocols described in Clinical Laboratory Standards Institute M40-A2 Approved Standard <sup>(4)</sup>. Performance testing was not conducted using human specimens. Please refer to your internal procedures to choose the most appropriate device for the specific sampling site. Educational material related to sample collection is available on Copan website.

#### Specimen Collection

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals <sup>(2, 17, 18, 20, 21, 22)</sup>.

Do not use the ESwab® medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.

**For the ESwab® collection systems composed of one test tube with medium and only 1 flocked swab (Fig.3):**

1. Open the ESwab® sample collection pouch and remove the tube and swab.
2. Collect the sample from the patient.
3. Unscrew and remove the cap from ESwab® tube making sure not to spill the medium.
4. Insert the swab into the tube until the red marked breaking point is at the level of the tube opening.
5. Break the swab off into the tube as follows:
  - With the other hand grasp the swab shaft at the very end with the thumb and first finger.
  - Lean the part of the shaft with the breaking point against the rim of the tube.
  - Bend the swab shaft at a 180 degrees angle to break it off at the colored ink breakpoint mark. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
  - Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
6. Replace cap on the tube and secure tightly.
7. Write patient information on the tube label or apply patient identification label. Send the sample to the test laboratory.

For ESwab® MRSA collection system codes 493CE02, 493CE03 and 4E047S02:

1. Open the ESwab® sample collection pouch and remove the tube and one pink swab.
  2. Use pink swab to collect first specimen (i.e: throat, perineum, nose or any other collection site).
  3. Unscrew and remove the cap from ESwab® tube making sure not to spill the medium. Insert the swab into the tube. Dip and gently stir the swab for 5 seconds.
  4. Lift up the swab from the liquid medium and swirl the swab against the tube walls 5 times to allow release of the sample from the flocked fibre holding the tube away from your face. Remove the swab and recap.
  5. Discard pink swab in the Biohazard container.
- Repeat all previous steps (2 to 5) if your ESwab® MRSA SYSTEM contains more than one pink swab and you use the second pink swab to collect the second specimen (i.e: throat, perineum, nose or another collection site). If not, proceed to step 6.
6. **Use white swab to collect the last specimen** (i.e: throat, perineum, nose or any other collection site) and then break the swab at the molded breaking point.
  7. Break the swab off into the tube as follows (Fig.3):
    - With the other hand grasp the swab shaft at the very end with the thumb and first finger.
    - Lean the part of the shaft with the breaking point against the rim of the tube.
    - Bend the swab shaft at a 180 degrees angle to break it off at the colored ink breakpoint mark. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
    - Discard the broken handle part of the swab shaft into an approved medical waste disposal container.

For ESwab® collection system codes 492CE03:

**NOTE:** Only one swab is required for the entire procedure. **The 2 remaining swabs MUST NOT be used and MUST be discarded.**

1. Open the ESwab® sample collection pouch and remove the tube and the chosen swab (Choose among pink, green or white swab to collect the specimen according to internal sampling procedures).
2. Collect the sample from the patient.
3. Break the swab off into the tube as follows (Fig. 3):
  - With the other hand grasp the swab shaft at the very end with the thumb and first finger.
  - Lean the part of the shaft with the breaking point against the rim of the tube.
  - Bend the swab shaft at a 180 degrees angle to break it off at the colored ink breakpoint mark. If needed, gently rotate the swab shaft to complete the breakage and take away the upper part of the swab shaft.
  - Discard the broken handle part of the swab shaft into an approved medical waste disposal container.


**One swab only must be placed into the tube.**

4. Replace cap on the tube and secure tightly. Write patient information on the tube label or apply patient identification label.
5. Send the sample to the test laboratory.

#### **NOTE ABOUT CAPTURE CAP FEATURE:**


The tube cap is designed to have capturing feature when white or pink flock swab, is used and broken inside the tube. If green flock swab is used then the capturing feature is not guaranteed. Green swab must be carefully removed from the tube immediately before proceeding with sample analysis. This precaution is intended to avoid accidental falling of the swab from the cap as the broken applicator may not firmly fit into the cap.

Fig. 3 Specimen Collection

Sterile gloves and protective clothing and eyewear should be worn when collecting and handling microbiology specimens and care should be taken to avoid splashes and aerosols when breaking the swab stick into the tube of medium. During sample collection when handling the swab applicator, the operator must not touch the area below the colored breakpoint indication line; that is the area from the line to the tip of the nylon flocked swab (see Fig 4), as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.

Fig4. Collection swab showing breakpoint indication line and area for holding the applicator

**NOTE:** Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in accidental breakage of the swab shaft. Swab shafts often exhibit diameter changes to facilitate different sampling requirements. Swab shafts may also have a molded breakpoint point designed for intentional breakage of the swab into the transport tube.

#### **Plating ESwab® Specimen Cultures in the Laboratory**

ESwab® samples should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines <sup>(17, 18, 21, 24, 25)</sup>.

Culture investigations of swab specimens for the presence of aerobic bacteria, anaerobic bacteria and fastidious bacteria such as *Neisseria gonorrhoeae* routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for the inoculation of ESwab® samples using regular swabs like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum is as follows:

#### Plating ESwab® Specimen with regular swab (Fig. 5.1):

1. Vigorously shake the ESwab® tube containing the swab sample between the thumb and forefinger for 5 seconds or mix the tube using a vortex mixer for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
2. When the sample is processed with molecular assays, before plating the swab for culture, transfer an aliquot of the sample in a sterile tube.
3. Unscrew the ESwab® cap and remove the swab applicator.
4. Use the tip of the ESwab® applicator to release a drop of sample suspension onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
5. If it is necessary to culture the swab specimen onto a second culture media plate, return the ESwab® applicator to the transport medium tube for two seconds to absorb and recharge the applicator tip with transport medium/patient sample suspension then repeat Step No. 3.
6. If it is necessary to inoculate additional culture media plates, return the ESwab® applicator to the transport medium tube and recharge the swab applicator tip with the transport medium/patient sample suspension before inoculating each additional plate.

5.1

#### Plating ESwab® Specimen with all product configurations (Fig. 5.2):


1. Vigorously shake the ESwab® tube containing the swab sample between the thumb and forefinger for 5 seconds or mix the tube using a vortex mixer for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
2. When the sample is processed with molecular assays, before plating the swab for culture, transfer an aliquot of the sample in a sterile tube.
3. Unscrew the ESwab® cap and and transfer 100µl volumes of the suspension onto each culture plate using a volumetric pipettor and sterile pipet tips.

5.2

Standard laboratory techniques should be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig. 6)

#### Fig 6. Procedure for streaking ESwab® specimens on agar Petri dishes for primary isolation <sup>(33)</sup>

Seed a primary inoculum of ESwab® specimen onto the surface of an appropriate agar culture plate in the first quadrant.

Use a sterile bacteriology loop to streak the primary inoculum across the surface of the second, third and fourth quadrants of the agar culture plate.

#### Processing ESwab® with automated systems

Some ESwab® product references may be processed with automated systems. Please refer to automation manufacturer's instructions for ESwab® processing. Due to the flexibility of the shaft of pernasal and pediatric swabs, the swab will coil when placed in the tube. Use tweezers for removing the pernasal or pediatric swabs from the tube before loading it into the machine, as it may interfere with the automation. Use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.

#### Preparation of Gram Stain Smears of ESwab® Specimens

Laboratory analysis of clinical swab samples collected from certain sites on the patient can routinely include microscopic examination of stained preparations ("direct Smears") using the Gram stain procedure. This can provide valuable information to physicians who are managing patients with infectious diseases <sup>(28)</sup>. There are many instances in which a Gram stain can assist in making a diagnosis; for example, with swabs taken from the endocervix or male urethra to investigate suspected *Neisseria gonorrhoeae* infections or vaginal swabs to diagnose bacterial vaginosis <sup>(27, 28, 29, 30, 31, 39)</sup>. The Gram stain can also help to judge specimen quality and contribute to the selection of culture media especially with mixed flora <sup>(32)</sup>.

Microscope slides of patient specimens transported in Copan ESwab® transport system can be prepared for Gram stain analysis, as describe below, by sampling an aliquot of vortexed suspension of the swab <sup>(21, 32)</sup>. Sample transported in ESwab® elution medium represent an homogeneous suspension in liquid phase. It can be uniformly smeared allowing clear and easy reading.

1. Take a clean glass microscope slide, place it on a flat surface and inscribe an area using a diamond-tipped or similar glass marker to identify the location of the specimen inoculum. Note: a slide with a pre-marked 20 mm well can be used.
2. Vortex mix the ESwab® tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the Liquid Amies transport medium.
3. Unscrew the ESwab® cap and using a sterile pipet, transfer 1 – 2 drops of Liquid Amies sample suspension to the inscribed area on the glass slide. Note: about 30µl would be a suitable amount of liquid for the a pre-marked 20 mm well slide.  
In case of bloody or thicker specimens particular care should be taken to thinly spread the sample on the slide. Bacteria are difficult to detect if the sample shows many red cells and debris.
4. Allow the specimen on the slide to air dry at room temperature or place the slide in an electric slide warmer or incubator set at a temperature not exceeding 42°C.
5. Fix smears using methanol. Methanol fixation is recommended as it prevents lysis of Red Blood Cells, avoids damage to all host cells and results in a cleaner background <sup>(21, 26, 32)</sup>.
6. Follow published laboratory reference manuals and guidelines for performing the Gram stain. If commercial Gram stain reagents are used, it is important to comply with instructions in the manufacturer's product insert for performance test procedure.

For further information or guidance on the preparation of specimen slides for microscopic analysis, for information on Gram staining procedures and the interpretation and reporting of microscopic analysis, consult published laboratory reference manuals <sup>(20, 24, 25, 26, 32)</sup>.

#### Processing ESwab® specimens for molecular testing in the laboratory.

Specimens received in the laboratory for nucleic acid detection should be processed when received in the laboratory. In case of delay, please refer to the appropriate specimen storage conditions.

When working with molecular methods care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over. <sup>(42)</sup>

1. Mix ESwab® tube with a vortex mixer for 10 seconds, unscrew the cap and using it as a handle, spin the cap in between the thumb and index finger to drain most of the fluid from the tip.  
Due to the flexibility of the shaft of pernasal and pediatric swabs, the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap.

In this case, use tweezers to extract the applicator from the tube after spinning the swab tip to extract most of the fluid from the tip.

- Discard the swab, and transfer the appropriate amount of sample to an extraction tube as per laboratory SOP.
- The ESwab<sup>®</sup> has been validated with the following extraction methods: Silica-gel membrane, Magnetic beads, Organic Extraction Method, Thermal extraction. Other extraction methods may be also applicable prior validation.
- Store E-Swab specimen at -20°C when unable to extract.

#### **Processing ESwab<sup>®</sup> specimens for rapid antigen testing in the laboratory.**

- Mix ESwab<sup>®</sup> tube with a vortex mixer for 10 seconds.
  - Use the sample fluid or the swab and test according to testing specifications provided with the kit and laboratory SOP.
- Due to the flexibility of the shaft of perinasal and pediatric swabs, the capture cap feature is not applicable, as the broken applicator may not firmly fit into the cap. Use tweezers to extract the applicator from the tube.

ESwab<sup>®</sup> has been tested with some rapid lateral flow immunochromatographic kits available on the market. For a complete list of the kits tested, contact Copan Customer Care ([customer care@copangroup.com](mailto:customer care@copangroup.com)).

#### **QUALITY CONTROL**

All lot numbers of the ESwab<sup>®</sup> are tested for sterility, nuclease and inhibitors and all lot numbers of swab applicators are tested to ensure they are non-toxic to bacteria. ESwab<sup>®</sup> Liquid Amies transport medium is tested for pH stability and bio-burden using Gram stain microscopic examination to ensure acceptable levels as defined in Clinical Laboratory Standards Institute M40-A2 <sup>(4)</sup>. Each production lot of ESwab<sup>®</sup> is quality control tested before release for ability to maintain viable bacteria at both refrigerated temperatures (4 – 8°C) and room temperature (20 – 25°C) for specified time points with a panel of aerobes, anaerobes and fastidious bacteria using both Roll-Plate and Swab Elution Methods <sup>(4)</sup>. Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C) which should correspond to  $\leq 1$  log increase in growth at a specified time point. Each production lot of ESwab<sup>®</sup> is analyzed for enzymatic and inhibitory activity, which may prevent nucleic acids amplification. DNase and RNase are enzymes that degrade nucleic acids thus preventing suitable nucleic acid amplification. The presence of DNase or RNase in the transport and storage medium may result in false negative results. Testing consists in adding a known amount of DNA or RNA (Kb ladder) to the ESwab<sup>®</sup> medium and evaluating the level of DNA and RNA integrity.

Procedures for quality control of bacteriology transport devices using a quantitative Swab Elution Method and qualitative Roll-Plate Method are described in Clinical Laboratory Standards Institute M40-A2 and other publications <sup>(4, 10, 12, 14, 15, 40, 41)</sup>. If aberrant quality control results are noted, patient results should not be reported.

#### **LIMITATIONS**

- In the laboratory, wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other CDC Biosafety Level 2 recommendations <sup>(34, 35, 36, 37)</sup> when handling or analyzing patient samples.
- The use of ESwab<sup>®</sup> for the collection of samples from the urogenital tract of pregnant women has not been evaluated.
- Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection <sup>(2, 3, 17, 18, 20, 21, 24)</sup>.
- ESwab<sup>®</sup> is intended for use as a collection and transport medium for aerobes, anaerobes and fastidious bacteria such as *Neisseria gonorrhoeae*, as well as for rapidly searching for bacterial, viral and Chlamydia antigens using rapid lateral flow immunochromatographic methods, and for nucleic acids. The product is not intended to maintain the viability of virus and Chlamydia.
- ESwab<sup>®</sup> Collection and Transport System is intended to be used with the medium tubes and swabs provided in the pouch. The use of tubes of medium or swabs from any other source are not qualified for use with ESwab<sup>®</sup> and could affect the performance of the product and laboratory test results.
- ESwab<sup>®</sup> has been validated with the main extraction methods: Silica-gel membrane, Magnetic beads, Organic Extraction Method, Thermal extraction. Other extraction methods can be performed prior validation.
- After DNA extraction, an aliquot of ESwab<sup>®</sup> medium can be amplified with no purification phase. In this case we suggest a dilution of ESwab<sup>®</sup> medium in 1:5.
- The capture cap feature is not applicable to perinasal and pediatric swabs (see Table 1 for capture cap feature applicability).
- Traces of nucleic acids from non-viable microorganisms could be contained into ESwab<sup>®</sup> that may be amplified by PCR-based tests depending on the analytical sensitivity of the assay. Refer to the assay manufacturers' instructions for use and internal laboratory procedures to manage results from specimen providing low amplification (high Ct value) of the target microorganism.

#### **WARNINGS**

- This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
- Do not re-sterilize unused swabs.
- Do not re-pack.
- Not suitable to collect and transport microorganisms other than aerobes, anaerobes and fastidious bacteria.
- Not suitable for any other application than intended use.
- The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
- The compatibility of ESwab<sup>®</sup> as collection and transport device suitable for use with PCR-based tests must be qualified according to internal laboratory procedures.
- Do not use if (1) there is evidence of damage (i.e., if the swab tip or swab shaft is broken) or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed, (4) the swab package is open, or (5) there are other signs of deterioration.
- Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
- Applicator swab is qualified as Class IIa Medical Device according to European Medical Device Directive 93/42/EEC - Surgically Invasive Transient Use. Class IIa means swabs can be used for sampling body surfaces, body orifices (e.g., nose, throat and vagina and deep invasive surgical wounds).
- Do not ingest the medium.
- Directions for use must be followed carefully.
- To be handled by trained personnel only.
- Due to the design of the perinasal and pediatric shaft applicators, the swab will coil when placed in the tube. Therefore, in manual plating of ESwab<sup>®</sup> specimen, it is not recommended to remove the swab from the tube. To process the specimen, collect the liquid using a sterile pipet. If the user must remove the swab, use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.
- It must be assumed that all specimens contain infectious micro-organisms; therefore all specimens must be handled with appropriate precautions. After use, tubes and swabs must be disposed of according to laboratory regulations for infectious waste. Observe CDC Biosafety Level 2 recommendations <sup>(34, 35, 36, 37)</sup>.
- Do not use the ESwab<sup>®</sup> medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
- The tube cap should be tightly closed to guarantee that the broken swab is captured by the cap.
- The capture cap feature for minitip and urethral swabs is guaranteed only if the swab shaft is completely straight. If the swab shaft is bent, the capture cap feature could be compromised and the shaft may not be captured by the cap.
- The plating procedure onto solid agar in Petri dishes using the minitip and urethral swab as inoculation wand to transfer the sample, it is not recommended. This procedure is suitable only with regular swabs.

#### **RESULTS**

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

#### **PERFORMANCE CHARACTERISTICS**

In the routine clinical laboratory, the Roll-Plate Method is the primary means of inoculating swab transport devices onto plated media. A limitation of the Roll-Plate Method <sup>(4)</sup> for bacterial viability performance testing is that it is not a quantitative method; it is, at best, a semi-quantitative approximation. On the other hand, quantitative viability performance methods such as the Swab Elution Method <sup>(4)</sup> do not reflect the standard protocol used in most clinical laboratories. Whereas the Swab Elution Method allows a quantitative measurement of the ability of a transport system to maintain viable organisms, the Roll-Plate technique takes into consideration some mechanical variables of the direct swabbing action that exist in the clinical laboratory, and which can influence the release of the sample onto culture plates.

Because of this, both methods of performing viability studies were used to determine the performance characteristics of the Copan eSwab® Collection and Transport System. The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in Clinical Laboratory Standards Institute M40-A2<sup>(4, 10, 12, 14, 15, 40, 41)</sup>.

The test organisms utilized in this study were those specifically prescribed in M40-A2 for establishing performance claims and quality control of swab transport systems and include a representative panel of aerobes, anaerobes and fastidious bacteria. An additional group of organisms not required or specified by M40-A2 were tested in order to provide further information on the survival of specific bacteria. Bacterial viability studies were performed on the Copan eSwab® at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs, 24 hrs and 48 hrs.

At the appropriate time intervals, each swab was processed according to the Roll-Plate or Swab Elution Method. eSwab® Collection and Transport System is able to preserve DNA, RNA and antigens of bacteria, viruses and Chlamydia for five days when stored at room temperature (20 – 25°C); 7 days if stored at 4°C and up to 6 months when stored at -20°C.

eSwab® is free of contamination of DNase and RNase. DNase and RNase are enzymes that interfere with the amplification process.

Organisms evaluated were divided into three main groups (see note below):

1. Aerobes and Facultative Anaerobes:  
*Pseudomonas aeruginosa* ATCC® BAA-427, *Streptococcus pyogenes* ATCC® 19615, *Streptococcus pneumoniae* ATCC® 6305, *Haemophilus influenzae* ATCC® 10211.
2. Anaerobes:  
*Bacteroides fragilis* ATCC® 25285, *Peptostreptococcus anaerobius* ATCC® 27337, *Fusobacterium nucleatum* ATCC® 25586, *Propionibacterium acnes* ATCC® 6919, *Prevotella melaninogenica* ATCC® 25845.
3. Fastidious Bacteria:  
*Neisseria gonorrhoeae* ATCC® 43069.

Additional organisms evaluated:

*Enterococcus faecalis* (Vancomycin resistant Enterococcus VRE) ATCC® 51299, *Staphylococcus aureus* (Methicillin resistant *Staphylococcus aureus* MRSA) ATCC® 43300, *Streptococcus agalactiae* (Group B Streptococcus) ATCC® 13813, *Clostridium perfringens* ATCC® 13124, *Clostridium sporogenes* ATCC® 3584, *Fusobacterium necrophorum* ATCC® 25286, *Peptococcus magnus* ATCC® 29328.

## NOTE

For product performance claims and viability performance testing, bacteria are categorized into three groups as described in Clinical Laboratory Standards Institute M40-A2<sup>(4)</sup> according to their growth responses to atmospheric oxygen:

1. Aerobes and Facultative Anaerobes:  
Aerobic bacteria require air or free oxygen to live. Facultative anaerobes are bacteria that can survive in either the presence or absence of oxygen. Many aerobic bacteria are facultative anaerobes meaning they are able to grow and survive in the absence of oxygen. For this reason, the aerobic group includes the description facultative anaerobes.
2. Anaerobes:  
Anaerobic bacteria do not require air or free oxygen to live. This category includes obligate anaerobes that can only live in the absence of oxygen.
3. Fastidious Bacteria:  
Fastidious bacteria have complicated or exacting growth requirements and this group is represented by the bacterium *Neisseria gonorrhoeae*.

In accordance with Clinical Laboratory Standards Institute M40-A2, with the exception of *Neisseria gonorrhoeae*, viability performance is measured for each test organism at the 48 hrs time point and compared with the acceptance criteria. Viability performance is measured for *Neisseria gonorrhoeae* at the 24 hrs time point. In both the Roll-Plate and Swab Elution viability performance studies, Copan eSwab® System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) and room temperature (20 – 25°C). Acceptable recovery for the Roll-Plate Method is defined as ≥5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU. Acceptable recovery for the Swab Elution Method is defined as no more than a 3 log<sub>10</sub> (1 x 10<sup>3</sup> +/- 10%) decline in CFU between the zero-time CFU count and the CFU of the swabs after the specified holding time. Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C). For the Swab Elution Method, an overgrowth assessment is made on all bacteria species tested at the 48 hrs holding time point except for *Neisseria gonorrhoeae* which is assessed at the 24 hrs holding time point. Overgrowth assessment using the Swab Elution Method is defined as greater than 1 log<sub>10</sub> increase in CFU between the zero-time CFU count and the holding time point. For the Roll-Plate Method, an overgrowth assessment is made with a separate analysis in which swabs are dosed with 100µl containing 10<sup>2</sup> CFU of *Pseudomonas aeruginosa* culture. Overgrowth under these conditions is defined as greater than 1 log<sub>10</sub> increase in CFU between zero-time CFU and the 48 hrs holding time point. Copan eSwab® Collection and Transport System demonstrated no overgrowth in either the Swab Elution or Roll-Plate Methods based on the acceptance criteria described in Clinical Laboratory Standards Institute M40-A2.

Please refer to symbol table at the end of the instructions for use

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UA.TR.116

## BIBLIOGRAPHY

1. Amies CR. A modified formula for the preparation of Stuart's medium. *Canadian Journal of Public Health*, July 1967. Vol. 58, 296 – 300.
2. Miller JM. A Guide to Specimen Management in Clinical Microbiology. Second Edition. American Society for Microbiology. Washington, DC. 1999.
3. Miller JM, Holmes HT. Specimen collection, transport, and storage. In: *Manual of Clinical Microbiology*. 6th ed. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. Washington, DC: ASM; 1995:19-20.
4. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 2014. *Quality Control of Microbiological Transport Systems*; Approved Standard- Second Edition. M40-A2 Vol. 34 No. 9.
5. Sng E-H, Rajan VS, Teo K-L, Goh A-J. The recovery of *Neisseria gonorrhoeae* from clinical specimens: effects of different temperatures, transport times, and media. *Sex Trans Dis*. 1982; 9:74-78.
6. Sun Y, Taylor T, Williams L, Sautter RL. Comparison of bacterial viability using both the EZ brand collection and transport system with the Difco swab transport pack. Presented at: 96th ASM General Meeting. 1996; Washington DC. Abstract C35.
7. Arbique JC, Forward KR, LeBlanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. *Diagnostic Microbiology and Infectious Disease*. 2000; 36:163-168.
8. Perry JL. Effects of temperature on fastidious organism viability during swab transport. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-55.
9. Wilson DA, Tuohy MS, Procop GW, Hall GS. Effects of storage on the recovery of bacteria from three swab transport systems: BD CultureSwab®, BD Culturette and Starplex StarSwab II. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-61.
10. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies described in NCCLS document M40-P. *Quality Control of Microbiology Transport Devices*. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-40.
11. Mitchell E, Berman M, Ginocchio CC. Evaluation of two new Liquid Stuart transport systems: Platinum StarSwab II (Starplex Scientific) and BBL CultureSwab® (Becton Dickinson). 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-74.
12. Perry JL, Matthews JS. Compliance of two popular swab transport systems with performance standards detailed by the new NCCLS Proposed Standard, M40-P. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-42.
13. Robinson A, Gruver ML. Comparison of bacterial survival in two transport systems stored at room temperature and refrigerator temperatures. 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-69.
14. Human RP, Jones GA. Evaluation of 4 transport systems against a published standard. 104th General Meeting of the American Society for Microbiology. 2004; New Orleans, LA. Abstract C-161.
15. Human RP, Jones GA. Evaluation of swab transport systems against a published standard. *J Clin Pathol* 2004; 57:762-763.
16. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies for anaerobic organisms described in NCCLS document M40-P. *Quality Control of Microbiology Transport Devices*. 13th European Congress of Clinical Microbiology and Infectious Disease (ECCMID). 2003; Glasgow, UK. Abstract P-652.
17. Isenberg HD, Schoenkench FD, Von Graevenitz A. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S.J. Rubin. American Society for Microbiology, Washington, DC, 1979.
18. Koneman EW, Allen SD, Janda WM, Schreckenberger PC and Winn, Jr. WC. 1992. *Color Atlas and Textbook of Diagnostic Microbiology*. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
19. 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
20. Forbes BA, Sahm DF, Weissfeld AS. 1998. *Bailey and Scott's Diagnostic Microbiology*. 10th ed. Mosby, St. Louis, MO.
21. Isenberg HD. 2004. *Clinical Microbiology Procedures Handbook*, 2nd ed. ASM, Washington, DC.
22. Isenberg HD. 1998. *Essential Procedures for Clinical Microbiology*. Chapter 14.12, Page 787. Packaging and Shipping Infectious Substances. ASM, Washington, DC.
23. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 1994. *Procedures for Handling and Transport of Diagnostic Specimens and Etiologic Agents*; Approved Standard. H5-A3.
24. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. *Manual of Clinical Microbiology*. 7th edition. Washington, DC: ASM; 1999.
25. Summanen P, Baron EJ, Citron D, Strong C, Wexler HM, Finegold SM. (1993). *Wadsworth Anaerobic Bacteriology Manual*, 5th ed. Star Publishing Company, Belmont, CA.
26. Marler LM, Siders JA, Allen SD. *Direct Smear Atlas, A Monograph of Gram-Stained Preparations of Clinical Specimens*. Lippincott Williams and Wilkins, 2001.
27. Rotimi VO, Yakubu Z, Abudu OO, Banjo TO. Direct Gram's stain of vaginal discharge as a means of diagnosing bacterial vaginosis. *Journal of Medical Microbiology*. 1991 Vol 35, Issue 2 103-106.
28. Spiegel CA, Amsel R, Holmes KK. Diagnosis of bacterial vaginosis by direct gram stain of vaginal fluid. *J. Clin Microbiol*. 1983 Jul; 18 (1):170-177.
29. Benavides MI, Moncada X, Rodriguez B, Castillo C. Gonococcal urethritis in men: clinical experience in 1978-1988. *Rev Med Chil*. 1992 Oct; 120(10):1140-3.
30. Mayaud P, Msuya W, Todd J, Kaatano G, West B, Begkoyian G, Grosskurth H, Mabey D. Rapid assessment in Rwandan refugee camps in Tanzania. *Genitourin Med*. 1997 Feb; 73 (1):33-8.
31. Deceuninck G, Asamoah-Adu K, Khonde N, Pepin J, Frost EH, Deslandes S, Asamoah-Abu A, Bekoe V, Alary M. Improvement of clinical algorithms for the diagnosis of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* by the use of Gram-stained smears among female sex workers in Accra, Ghana. *Sex Transm Dis*. 2000 Aug; 27 (7):401-10.
32. Isenberg HD. 1998. *Essential Procedures for Clinical Microbiology*. Chapter 2.1, Page 41. Gram Stain. ASM, Washington, DC.
33. Isenberg HD. 1998. *Essential Procedures for Clinical Microbiology*. Chapter 1.1, Page 27. Collection, Transport and Manipulation of Clinical Specimens. Procedure for streaking plates for primary isolation. ASM, Washington, DC.
34. Fleming D. *Biological Safety: Principles and Practices*. January 2000. ASM, Washington DC.
35. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/fjyrtxt.htm>.
36. Richardson JH. *Biosafety in Microbiological and Biomedical Laboratories*. December 1994. Diane Publishing Company.
37. Hansen DJ. *Healthcare, Laboratories and Biosafety*. Vol 2., 1992. CRC Press.
38. Greenberg AE, Clesceri LS, and Eaton AD. 9215 heterotrophic plate count. In: *Standard Methods for the Examination of Water and Waste Water*. 18th ed. Washington, DC APHA; 1992: 9-33-9-34.
39. Washington JA. 1986. Rapid diagnosis by microscopy. *Clin. Microbiol. Newsl*. 8:135-137.
40. Van Horn KG, Rankin I. Evaluation and comparison of two Stuart's Liquid Swab transport systems tested by the NCCLS M40 method. 105th General Meeting of the American Society for Microbiology. 2005; Atlanta, Georgia. Abstract C-292.
41. Bourbeau PP, Heiter BJ. Validation of QC standard for bacteriological transport devices as specified in the NCCLS Proposed Standard M40: *Quality Control of Microbiological Transport Systems*. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-46.
42. *The Nucleic Acid Amplification Assays for the Molecular Hematopathology*; approved Guideline (NCCLS MM5-A Volume 23 No 17)
43. Collection, transport, preparation and storage of specimens for molecular methods: Proposed guideline. CLSI (MM13-P volume 25 No 9)

Symbol / Simbolo / Símbolo / Symbol / Symbole / Símbolos / Symbol / Symbol / Szimbólum / Symbol / Symbol / Symbol / Símbol / Символ	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Význam / Betyder / Jelentés / Betydning / Znaczenie / Značenje / Betydelse / Anlami / Значення
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Výrobce / Producent / Gyártó / Produsent / Wytwórca / Proizvodač / Tillverkare / Üretici / Виробник
CE 0123	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Identifikační číslo notifikované osoby / Identifikationsnummer for det bemyndigede organ / A bejelentett szervezet azonosítószáma / Identifisering av godkjenningsorganet / Numer identyfikacyjny jednostki notyfikowanej / Identifikacioni broj prijavljenog tela / Det anmälda organets identifieringskod / Onaylı kurumun tanımlama numarası / Идентифікаційний номер уповноваженого органу
STERILE R	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiación ionizante / Sterilizována pomocí ionizujícího záření / Sterilisationsméde: Bestråling / Ionizáló sugárzással sterilizált / Sterilisert ved bruk av ioniserende stråling / Sterylizacja za pomocą promieniowania jonizującego / Sterilisoano upotrebom jonizujúceho žiarenja / Steriliserad med strålning / İşinlama kullanilan yöntem / Стерилізований завдяки опроміненню
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Nepoužívejte znovu / Brug det ikke igen / Ne használja újra / Må ikke gjenbrukes / Nie używać ponownie / Nemojte ponovo koristiti / Återanvänd inte / Tekrar kullanmayın / Не використовувати повторно
REF	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Číslo katalogu / Katalognummer / Katalógusszám / Katalognr / Numer katalogowy / Kataloški broj / Katalognummer / Katalog numarası / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Limites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Teplotní limity / Temperaturgränser / Hőmérsékletkorlátozás / Temperaturgrænser / Limity temperatury / Opseg temperature / Temperaturgränser / Sicaklik limitleri / Температурні обмеження
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Použijte do / Anvende s før / Felhasználható / Må brukes innen / Zużyć do / Upotrebiti do / Sista förbrukningsdag / Son kullanma tarihi / Використати до
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Consulte las instrucciones de uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consultar as instruções de utilização / Konzultujte návod k použití / Se brugsvejledningen / Tekintse meg a használati útmutatót / Se instruksjoner for bruk / Patrz instrukcja użycia / Pročitajte uputstvo za upotrebu / Läs bruksanvisningen / Kullanma talimatlarına bakınız / Користуватись інструкціями по використанню
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacar / Odrhňte a otevřete / Riv op for at åbne / Húzza le / Riv opp for å åpne / Oderwać, aby otworzyć / Ogulite / Dra för att öppna / Açmak için yırtın / Відкривати тут
LOT	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Kód šarže (dávky) / Serienummer (parti) / Gyártási szám (kötég) / Lot-nummer (parti) / Kod partii (batch) / Serijski (lot) broj / Partinummer (lotto) / Lot (seri) kodu / Код серії (партії)
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Obsah dostatečný pro <n> testů / Innhald tilstrekkelig til <n> prøver / A(z) <n> tesztekhez elegendő tartalom / Innhold tilstrekkelig for <n> prøver / Zawartość wystarczająca do przeprowadzenia <n> testów / Sadržaj dovoljan za   testova / Innehåller tillräckligt för <n> tester / <n> test için yeterli içerik / Придатний для проведення <n> кількості випробувань
	Do not use if package is damaged / Non utilizzare in caso di confezionamento danneggiato / No utilizar en caso de paquete dañado / Bei Beschädigung der Verpackung nicht verwenden / Ne pas utiliser si l'emballage est abîmé / Não utilizar se a embalagem estiver danificada / Nepoužívejte v případě poškozeného obalu / Må ikke anvendes, hvis emballagen er beskadiget / Ne használja, ha a csomagolás sérült! / Må ikke brukes hvis emballasjen er skadet / Nie używać, jeśli opakowanie jest uszkodzone / Nemojte koristiti ako je ambalaža oštećena / Använd inte vid skadad förpackning / Ambalajı hasar görmüşse kullanmayın / Не використовувати, якщо упаковка пошкоджена

**Copan**

**Copan Italia S.p.A.**

Via F. Perotti, 10  
25125 Brescia, Italy  
Tel +39 030 2687211  
Fax +39 030 2687250

Email: [info@copangroup.com](mailto:info@copangroup.com)  
Website: [www.copangroup.com](http://www.copangroup.com)

**North American Distributor:**

Copan Diagnostics Inc.  
26055 Jefferson Avenue  
Murrieta, CA 92562, USA  
Tel: 951-696-6957  
Fax: 951-600-1832

E-mail: [customerservice@copanusa.net](mailto:customerservice@copanusa.net)  
Website: [www.copanusa.com](http://www.copanusa.com)