

IVD solutions through partnership




CHROMagar™ Salmonella Plus

For detection and isolation of *Salmonella* species including lactose positive *Salmonella*

CHROMagar
The Chromogenic Media Pioneer


● CHROMagar™ Salmonella Plus



Opaque plate: with CHROMagar™ White Opaque Supplement

Plate Reading

- *Salmonella* including *S.typhi*, *S.paratyphi* A and lactose positive *Salmonella*
→ mauve
- Coliforms
→ blue
- *E.coli*
→ colourless
- *Proteus*
→ colourless or inhibited



For detection and isolation of *Salmonella* species including lactose positive *Salmonella* in food samples

Background

Mainly due to contamination in the food chain and/or during food-production processes, *Salmonella* commonly induces enteric illness whose major symptoms are abdominal cramps,

Medium Performance

1 Meets ISO 6579-1 Requirements

by detecting lactose positive *Salmonella* in intense mauve colour.

2 Easy to read by naked eye

Another feature of this new medium resides in its nice colour contrast due to the fact that *E.coli* are colourless. *E.coli*, which are usually present in abundance in samples tested for *Salmonella* and could potentially hide suspect colonies, are no longer a concern.

3 High sensitivity and specificity

Salmonella including lactose positive *Salmonella*
99%*

*Specificity and sensitivity from scientific study: "Evaluation of a new chromogenic medium CHROMagar Salmonella Plus for the detection of *Salmonella* species including lactose positive *Salmonella*, *S.Typhi* and *S.Paratyphi*" de Beaumont C., Breuil J., Dedicova D., Tran Q. 2006. Poster presented during ECC MID meeting.

Did you know? RambaQUICK™ *Salmonella* broth enrichment can be used ahead of the inoculation step of CHROMagar™ *Salmonella* Plus culture plate, following the RambaQUICK™ *Salmonella* Method to rapidly detect *Salmonella*. Product code for this broth enrichment is 206180 bzw. 15SQ001.

Medium Description

Powder Base	Total32.8 g/L Agar15.0 Peptones and yeast extract.....8.0 Salts8.5 Chromogenic mix1.3 Storage at 15/30°C - pH: 7.5 +/- 0.2 Shelf Life 3 years
+ Supplement (included in the pack)	1st: Liquid form10 mL/L Storage at 15/30 °C Shelf Life 5 years
CHROMagar White Opaque supplement	In order to obtain a white opaque background: Powder form: 1.0 g/L Storage at 15/30 °C Shelf Life 3 years

Usual Samples	food, meat, fresh eggs, milk products, clinical samples etc.
Procedure	Inoculation after selective broth enrichment of samples according to ISO 6579. Incubation 18-24h at 37 °C. Aerobic conditions.

Scientific Publications on this product: available on www.CHROMagar.com
Please read carefully the instructions for use (IFU document) available on www.CHROMagar.com

Manufacturer:

CHROMagar
4 place du 18 juin 1940 75006 Paris - France
e-mail: CHROMagar@CHROMagar.com
www.CHROMagar.com

Distributed by:
Mast Diagnostica GmbH
Feldstraße 20
DE-23858 Reinfeld

Tel: +49 (0)4533 2007 0
Fax: +49 (0)4533 2007 68
e-mail: mast@mast-diagnostica.de
www.mastgrp.com

Ordering Information

Product	Order Code
CHROMagar™ <i>Salmonella</i> Plus dry media, 5 liter	15SA162
CHROMagar™ <i>Salmonella</i> Plus ready to use plates, 20 pcs.	201421