# LBM<sup>®</sup> SLsolution<sup>™</sup>

## Instructions for Use

English	3
Italiano	5
Español	7
Deutsch	9
Français	12
Português	14
БЪЛГАРСКИ	16
Česky	18
Dansk	21
ΕΛΛΗΝΙΚΑ	23
Eesti Keel	25
Hrvatski	28
Latviešu	30
Lietuvių K.	32
Magyar	34
Nederlands	37
Norsk	39
Polski	41
Română	43
Slovenčina	45
Slovenščina	48
Srpski	50
Suomi	52
Svenska	54
Türkçe	56

## Copan SLsolution™ - Presentation and guide to use of the product

See the explanation of the symbols at the bottom of this leaflet.

### INTENDED USE

Copan SLsolution™ is a treatment reagent for fluidification of specimens collected from the respiratory tract prior to subsequent microbiological analyses.

### SUMMARY AND PRINCIPLES

Copan SLsolution™ consists of a test tube with screw cap containing a ready-to-use solution and is indicated for fluidification of specimens collected from the respiratory tract, *Mycobacterium spp* excluded.

Copan SLsolution™ is packed so as to maintain a modified atmosphere that guarantees the stability and activity of the DTT reagent<sup>1,2</sup> until the time of use of the product.

### PRODUCT DESCRIPTION

Copan SLsolution™ is available in different configurations as indicated below and supplied in a labelled and closed test tube with screw cap filled with different volumes of reagent. The kit versions also include a sterile specimen transfer device, Pasteur pipette or sputum dipper.

The A codes are compatible with the automated device.

REF	Product description	Packaging
0E004N.A	Copan SLsolution™ kit: - 50 12x80mm PET round-bottom test tubes with screw cap containing 1ml of reagent - 50 Pasteur pipettes	50 kits per package 6 x 50 kits per box
0E003N	Copan SLsolution™ kit: - 50 12x80mm PET round-bottom test tubes with screw cap containing 1ml of reagent - 50 Pasteur pipettes	50 kits per package 6 x 50 kits per box
0U020N.A	Copan SLsolution™ kit: - 50 12x80mm PET test tubes with prehensile screw cap containing 1ml of reagent - 50 sputum dippers	50 kits per package 6 x 50 kits per box
0E006N.A	Copan SLsolution™ in bulk: - 50 12x80mm PET round-bottom test tubes with screw cap containing 1ml of reagent	50 test tubes per package 6 x 50 test tubes per box
0E005N	Copan SLsolution™ in bulk: - 50 12x80mm PET round-bottom test tubes with screw cap containing 1ml of reagent	50 test tubes per package 6 x 50 test tubes per box
0U019N	Copan SLsolution™ in bulk: - 50 12x80mm PET round-bottom test tubes with prehensile screw cap containing 1ml of reagent	50 test tubes per package 6 x 50 test tubes per box
2U063S01	Copan Sputum Dipper in bulk: - 100 specimen transfer devices for prehensile cap	100 sputum dippers per package 10 x 100 sputum dippers per box

### REAGENTS

The Copan SLsolution™ components are:

Components	Quantity (g/l)
Saline matrix	46.45
DTT (DL-dithiothreitol) <sup>1,2</sup>	2.0

### REQUIRED MATERIALS BUT NOT PROVIDED

Specimen collection device; primary specimen transfer device (included and provided only for the kit versions); microorganism culture and isolation device.

### STORAGE

This product is ready to use and requires no further preparation. The product, which has a clear/colourless appearance, can be stored at a temperature of 5-25°C until the expiration date. Do not incubate or freeze. Keep away from light sources. Do not use after the expiration date. If incorrectly stored, its effectiveness will be compromised.

### LIMITATIONS

- The conditions, times and volume of the specimen collected for culture are significant variables for obtaining reliable culture results. Follow the recommended specimen collection guidelines.
- Copan SLsolution™ is not suitable for recovery of *Mycobacterium spp*.
- The performance tests on Copan SLsolution™ were conducted using laboratory strains inoculated in the Copan SLsolution™ test tube and not human clinical specimens.
- If the specimen collected has a high viscosity, it may be necessary to vortex the test tube again in order to ensure complete fluidification. The vortex time should be decided by the operator during use by visually checking specimen fluidification.
- The sputum dipper, if provided, is not a quantitative specimen transfer device. The quantity transferred may vary depending on the manual skill of the operator and the nature of the specimen.

### WARNINGS

- Single-use in vitro diagnostic device for professional use.
- Always use Copan SLsolution™ with filling proportional to the amount of specimen available in order to always ensure a 1:1 ratio.<sup>6,7,8</sup>







3. Do not perforate the membrane in any way.
4. Do not immerse the transfer device in the reagent prior to sampling.
5. Do not aspirate the reagent from the test tube or pour it into the primary specimen.
6. Do not syringe biological material or other reagents into the test tube.
7. Always use a test tube rack for the non-self-supporting round-bottom version.
8. Do not use Copan SLsolution™ if the reagent is not clear/colourless (e.g. pink or yellow).
9. The use of this product in combination with diagnostic tests or diagnostic instrumentation should be validated by the end user before use.
10. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.
11. The device is sterilized using ionizing radiation; do not use the device if the sterile barrier system is damaged.

#### INSTRUCTIONS FOR USE

1. Take the correctly identified Copan SLsolution™ test tube and unscrew the cap.
2. Transfer the specimen to the open test tube using a suitable transfer device, possibly provided in the kit version, respecting a ratio of about 1:1<sup>6,7,8</sup>

NOTE: Only if using the sputum dipper provided with the product codes 2U063S01 and 0U020N.A, break the transfer device in the test tube:

- A. Insert the sputum dipper in the test tube making sure that the tip touches the bottom.
- B. Break the sputum dipper with a decisive movement bending it 180° parallel to the test tube.
- C. Discard the upper broken off part of the sputum dipper.
3. Screw the cap back onto the Copan SLsolution™ test tube taking care that the sputum dipper (if used) is recovered from the prehensile part of the cap.
4. Immediately vortex the test tube for 30 seconds at 2000-2500 rpm.
5. Leave the test tube at ambient temperature for at least 15 minutes and do not exceed a maximum time of 6 hours.
6. Homogenize the specimen by vortexing the test tube for at least 3 seconds at 2000-2500 rpm, checking that the specimen is actually suitable for seeding.
7. Process the specimen according to the required internal procedures.

#### DISPOSAL

Dispose of all waste according to the national laws. Use precautions for infected material when necessary.

#### PERFORMANCE CHARACTERISTICS

A test on the survival of the respiratory pathogens listed below stored in Copan SLsolution™ for up to 6 hours at ambient temperature has shown that the organisms remain vital:

STRAIN	CFU AT TIME ZERO	CFU AFTER 6 HOURS AT AMBIENT TEMPERATURE
<i>Pseudomonas aeruginosa</i> ATCC® 27853	496	167
<i>Haemophilus influenzae</i> ATCC® 10211	420	88
<i>Staphylococcus aureus</i> ATCC® 6538	485	270
<i>Streptococcus pneumoniae</i> ATCC® 6305	450	184
<i>Streptococcus pyogenes</i> ATCC® 19615	466	178
<i>Moraxella catarrhalis</i> ATCC® 25238	501	137
<i>Candida albicans</i> ATCC® 10231	475	260
<i>Klebsiella pneumoniae</i> ATCC® 700603	433	321

¥ The tests were conducted using ATCC® strains rehydrated in the laboratory and not clinical specimens or artificial matrices.

#### QUALITY CONTROL PROCEDURE

- Starting from a fresh culture plate, prepare an inoculum equal to 0.5 McFarland of the organism under test in a test tube containing PBS.
- Prepare appropriate dilutions in PBS test tubes until obtaining a suspension with a concentration from  $4 \times 10^2$  to  $6 \times 10^2$  CFU/100 µl.

- Using a micropipette, inoculate 100 µl of the chosen dilution in the Copan SLsolution™ test tube.
- Vortex the inoculated test tube for 30 seconds at 2000-2500 rpm.
- Pipette 100 µl of inoculated Copan SLsolution™ onto a culture plate with an appropriate culture medium and with the aid of a spatula, evenly plate it on the entire surface of the culture medium.
- Incubate the plate at 35°C ±2°C for 24-48 hours in atmospheric conditions suitable for microorganism growth (specimen count at time zero of treatment).
- Leave the thus inoculated Copan SLsolution™ test tubes at ambient temperature (20-25°C) for 6 hours.
- Vortex the inoculated test tube for 30 seconds at 2000-2500 rpm.
- Pipette 100 µl of inoculated Copan SLsolution™ onto a culture plate with an appropriate culture medium and with the aid of a spatula, evenly plate it on the entire surface of the culture medium.
- Incubate the plate at 35°C ±2°C for 24-48 hours in atmospheric conditions suitable for microorganism growth (specimen count at 6 hours of treatment).

**ACCEPTABILITY LIMITS:** reduction of the bacterial load not higher than 1 log<sub>10</sub> after 6 hours at ambient temperature with respect to the count obtained at time zero.

#### TABLE OF SYMBOLS

See the table of symbols at the end of the instructions for use.

#### NOTES FOR THE PROFESSIONAL USER

In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

#### REVISION HISTORY

Last Revision No.*	Release date	Changes made
02	12-2022	Revision of IFU sections (first revision in IVDR)

\* Should you need earlier revisions, contact Copan Customer Service.

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## BIBLIOGRAPHY

1. A. Efthimiadis, M.M.M. Pizzichini, E. Pizzichini, J. Dolovich, F.E. Hargreave. Induced sputum cell and fluid-phase indices of inflammation: comparison of treatment with dithiothreitol vs phosphate-buffered saline. *European Respiratory Journal* 1997; 10: 1336–1340.
2. M R Hammerschlag, L Harding, A Maccone, A L Smith and D A Goldmann Bacteriology of sputum in cystic fibrosis: evaluation of dithiothreitol as a mucolytic agent. *J. Clin. Microbiol.* June 1980vol. 11 no. 6 552-557.
3. J. Michael Miller and Shelley A. Miller, *Specimen Management in Clinical Microbiology*, 3<sup>rd</sup> edition 2017, ASM Press, Washington DC, USA ; Section III- 103.
4. *Clinical Microbiology Procedures Handbook*, Lynne S.Garcia, 3<sup>rd</sup> ed. Volume 1 , Section 2: Specimen, Collection, Transport and Acceptability T2.1-1; T2.1-2.
5. *Clinical Microbiology Procedures Handbook*, Lynne S.Garcia, 3<sup>rd</sup> ed. Volume 1 , Section 3: Aerobic bacteriology, S3.11.2 lower Respiratory Tract Cultures .
6. Methods of sputum processing for cell counts, immunocytochemistry and *in situ* hybridisation; Q. Hamid, M.M. Kelly, M. Linden, R. Louis, M.M.M. Pizzichini, E. Pizzichini, C. Ronchi, F. Van Overveld, R. Djukanović *European Respiratory Journal* 2002 20: 19s-23s. DOI: 10.1183/09031936.02.00001902.
7. Induced sputum analysis: step by step; Zsoka Weiszhar, Ildiko Horvath *Breathe* 2013 9: 300-306; DOI: 10.1183/20734735.042912.
8. Methodology for Sputum Induction and Laboratory Processing; Guiot, J., Demarche, S., Henket, M., Paulus, V., Graff, S., Schleich, F., Corhay, J. L., Louis, R., Moermans, C. *Methodology for Sputum Induction and Laboratory Processing. J. Vis. Exp.* (130), e56612, doi:10.3791/56612 (2017).

Index of Symbols / Tabella dei Simboli / Tabla de símbolos / Symboltabelle / Table des Symboles / Tabela de símbolos / Индекс на символите / Seznam symbolů / Symbolforklaring / Ευρετήριο συμβόλων / Tingmärkide kirjeldused / Kazalo simbola / Simbolu rādītājs / Simbolių rodyklė / Szimbólumok jegyzéke / Index van symbolen / Symbolindeks / Indeks symboli / Index simboluri / Zoznam symbolov / Indeks simbolov / Tabela simbola / Symbolihakemisto / Förteckning över symboler / Semboller dizini

<p>Symbol / Simbolo / Simbolo / Symbol / Symbole / Símbolos / Символ / Symbol / Symbol / Συμβολο / Tingmärk / Simbol / Simbols / Simbolis / Szimbólum / Symbool / Symbol / Symbol / Simbol / Symbol / Símbol / Simbol / Simboli / Symbol / Sembol</p>	<p>Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значение / Význam / Betydning / Σημασία / Tāhendus / Značenje / Nozīme / Reikšmė / Jelentés / Betekenis / Betydning / Znaczenie / Semnificație / Význam / Pomen / Značenje / Merkitys / Betydelse / Anlam</p>
	<p>Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Производител / Výrobce / Producent / Κατασκευαστής / Tootja / Proizvođač / Ražotājs / Gamintojas / Gyártó / Fabrikant / Produzent / Producent / Producător / Výrobca / Proizvajalec / Proizvođač / Valmistaja / Tillverkare / Üretici</p>
	<p>In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In vitro Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / Инвитро диагностично изделие / Diagnostický prostředek pro in vitro diagnostiku / Udstyr til in vitro diagnostik / Διαγνωστικό προϊόν in vitro / In vitro diagnostika seade / In vitro dijagnostički proizvod / In vitro diagnostikas ierīce / „In vitro“ diagnostikos prietaisai / In vitro diagnosztikai eszköz / In vitro diagnosztikai eszköz / Hulpmiddel voor in vitrodiagnostiek / In vitro diagnostisk enhet / Urządzenie do diagnostyki in vitro / Dispositiv de diagnostic in vitro / Diagnostická pomůcka in vitro / In vitro diagnostični pripomoček / In vitro diagnostikakalaita / Dijagnostički uređaj u epruveti / Enhet för in vitro diagnostik / In vitro tani cihaz</p>
	<p>Unique Device Identifier / Identificativo unico del dispositivo / Identificador único de dispositivo / Einmalige Produktkennung / Identifiant dispositif univoque / Identificador de dispositivo único / Уникален идентификатор на изделието / Jediněčný identifikátor zařízení / Unik enheds id / Ατοκλειστική ταυτοποίηση ιατροτεχνολογικού προϊόντος / Seadme kordumatu identifitseerimistunnus / Jedinstveni identifikator proizvoda / Unikālais ierīces identifikators / Unikalus įrenginio identifikatorius / Egyedi eszközazonosító / Unieke hulpmiddelidentificatie / Identifikator for Unique enhet / Unikatowy numer urządzenia / Identificator unic dispozitiv / Unikátny identifikátor pomůcky / Edinstveni identifikator naprave / Jedinstveni identifikator uređaja / Yksilöivä laitetunniste / Unik enhetsidentifierare / Benzersiz Cihaz Tanımlayıcı</p>
	<p>CE marking / Marchio CE / Marca CE / CE Kennzeichnung / Marquage CE/ Marcação CE / CE маркировка / Označení CE / CE mærkning / Σήμανση CE / CE mārģis / Oznaka CE / CE markējums / CE ženklas / CE jelölés / CE marking / CE merking / Oznakowanie CE / Marcaj CE / Označenje CE / Oznaka CE / Oznaka CE / CE merkintä / CE märkning / CE isareti</p>
	<p>Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiação ionizante / Стерилизован с йонизиращо лъчение / Sterilizováno za použití ionizujících záření / Sterilisationsmåde: Bestrålnng / Αποστειρωμένο με ιονισμένη ακτινοβολία / Steriliseeritud ioniseeriva kiirgusega / Sterilizirano ionizirajućim zračenjem / Sterilizacija ar apstarošanu / Sterilizuota jonizuojančiąja spinduliuote / Ionizáló sugárzással sterilizálva / Gesteriliseerd met behulp van ioniserende straling / Sterilisert ved bruk av ioniserende stråling / Sterylizacja za pomocą promieniowania jonizującego/ Sterilizat cu radiații ionizante / Sterilizované ionizujúcim žiarením / Sterilizirano z uporabo ionizirajočega sevanja / Sterilisano upotrebom jonizujućeg zračenja / Steriloitu ionisoivalla säteilyllä / Steriliserad med strålning / İşinlama kullanılan yöntem</p>
	<p>Single sterile barrier system / Sistema a barriera sterile singola / Sistema de barrera estéril único / Einfaches Sterilbarriersystem / Système de barrière stérile unique / Sistema de barreira estéril simples / Единична стерилна барьерна система / Systemtjeddne sterilni bariéry / Enkelt sterilt barriersystem / Μονό σύστημα στείρου φραγμού / Ühekordne steriiline barjäärisüsteem / Jednostruki sterilni sustav barijera / Vienas sterilas barjeras sistēma / Vienkartinė sterilaus barjero sistema / Egyszeres steril gátrendszer / Enkel steriel barriersysteem / Sterilt barriersystem til engangsbruk / System pojedynczej bariery sterylnej / Sistem cu barieră sterilă unică / Sistem s jednou sterilnou bariérou / Enojnji sterilni pregradni sistem / Sistem jednostruke sterilne barijere / Yksittäinen steriili estejärjestelmä / Enkelt sterilt barriärsystem / Tek kullannmlik sterlii ambalaj sistemi</p>
	<p>Do not re-sterilize / Non risterilizzare / No reesterilizar / Nicht resterilisieren / Ne pas restériliser / Não reesterilizar / Не стерилизуйте отново / Výrobek opakovaně nesterilizujte / Må ikke gensteriliseres / Μην επαναποστείρετε / Mitte resteriliseerida / Ne sterilizirajte ponovno / Nesterilizēt atkārtoti / Nesterilizuokite pakartotinai / Ne sterilizálja újra / Niet opnieuw steriliseren / Må ikke steriliseres på nytt / Nie sterylizować ponownie / A nu se reesteriliza / Nesterilizujte znovu / Ne sterilizirajte ponovno / Nemojte ponovo sterilisati / Ei saa steriloida uudelleen / Får ej omsteriliseras / Tekrar sterilize etmeyin</p>

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**Copan Italia S.p.A.**
Copan Diagnostics Inc.