



Mast Group Ltd.  
Mast House  
Derby Road  
Bootle  
Merseyside L20 1EA  
United Kingdom

Tel: +44 (0) 151 933 7277

[www.mast-group.com](http://www.mast-group.com)

Dear Customer

As part of our ongoing product performance review, we have now changed the removing the tamper seal on our CAMP-ID and Uri<sup>®</sup>Locate bottles.

Feedback from our customers have shown that it is difficult to determine when to stop screwing the cap on due to the tightness of the anti-tamper feature and overtightening it reverts in the anti-tamper breaking resulting in leakage.

As a result of this change the appearance of the products will change and the performance in leak prevention will improve.

Integrity of the product bottles will remain the same with only the appearance being affected.



Improved



Previous



Improved



Previous



Mast Group Ltd.  
Mast House  
Derby Road  
Bootle  
Merseyside L20 1EA  
United Kingdom

Tel: +44 (0) 151 933 7277

[www.mast-group.com](http://www.mast-group.com)

Should you wish to discuss this change further please do not hesitate to contact your Account Manager.

Yours sincerely

A handwritten signature in black ink that reads "Ashleigh Michell".

Ashleigh Michell

Product Manager - Microbiology