

## SAFETY INFORMATION **ENZY-WELL TOXOPLASMA IGM**

Edition: 3

Date: 08-Jul-2014

Page: 1/9

The kit Enzy-Well Toxoplasma IgM (REF 91041 - 91179) is an immunoenzymatic capture method for the qualitative determination of IgM-class antibodies to Toxoplasma gondii in human serum.

Each kit is composed of:

MT PLATE

**MICROPLATE** 

12 x 8 wells coated with anti-human IgM monoclonal antibodies (91041)

10 x 96 wells coated with anti-human IgM monoclonal antibodies (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

CONTROL +

POSITIVE CONTROL

1 x 1.6 mL (91041)

4 x 1.6 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: diluted human serum, containing anti-Toxoplasma gondii IgM antibodies, in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.

CONTROL CUT-OFF CONTROL CUT-OFF

1 x 2.5 mL (91041)

5 x 2.5 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: diluted human serum, containing anti-Toxoplasma gondii IgM antibodies, in phosphate buffer 0.01 mol/L with BSA 1% and sodium azide 0.09%.

Ag

**ANTIGEN** 

Freeze-dried powder x 6 vials (91041)

Freeze-dried powder x 40 vials (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: Purified Toxoplasma gondii in phosphate buffer containing mouse ascitic fluid and lactose.

CONJ

**CONJUGATE** 

1 x 18 mL (91041) 7 x 18 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: monoclonal antibodies labelled with peroxidase, in phosphate buffer solution containing phenol 0.05% and Bronidox 0.02%.

CONTROL IgM -

IgM NEGATIVE CONTROL

1 x 1.6 mL (91041)

4 x 1.6 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and / or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: human serum diluted in phosphate buffer 0.01 mol/L, with BSA 1% and sodium azide 0.09%.

WASH BUF 10x

WASHING BUFFER 10X

1 x 100 mL (91041)

5 x 100 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive  $\underline{1999/45/CE} \ and \ does \ not \ contain \ hazardous \ substances \ and/or \ with \ limit \ values \ for \ exposure \ in \ concentrations$ above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: Phosphate buffered saline (PBS), concentrated 10 times; contains Brij 0.5%.

Safety Information – Enzy-Well Toxoplasma IgM (REF: 91041 - 91179)



# SAFETY INFORMATION ENZY-WELL TOXOPLASMA IgM

Edition: 3

Date: 08-Jul-2014

Page: 2/9

SAMP DIL 2

DILUENT 2

1 x 100 mL (91041) 6 x 100 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

Contents: Proteic solution in phosphate buffer with sodium azide 0.09% containing methyl orange 0.04%.

SUBS TMB

**SUBSTRATE** 

1 x 12 mL (91041) 1 x 120 mL (91179)

The substance is not classified as dangerous according to Regulation 1272/2008 and European Directive 1999/45/CE and does not contain hazardous substances and/or with limit values for exposure in concentrations above the established limits.

For this mixture, according to EC Regulations 1907/2006 and 453/2010, is not necessary to develop a material safety data sheet.

<u>Contents</u>: Tetramethylbenzidine 0.26 mg/mL and hydrogen peroxide 0.01% stabilized in citrate buffer 0.05 mol/L (pH 3.8).

H<sub>2</sub>SO<sub>4</sub> 0.3 M

STOP SOLUTION

1 x 16 mL (91041) 1 x 120 mL (91179)

The substance is classified as **dangerous** according to Regulation 1272/2008 and European Directive 1999/45/CE.

The material safety data sheet, prepared in accordance with EC regulations 1907/2006 and 453/2010, is available below.

Contents: solution of sulphuric acid, 0.3 mol/L.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 3/9

### 1. IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND THE COMPANY/UNDERTAKING

1.1 IDENTIFICATION OF THE PRODUCT

Product name: ENZY-WELL STOP SOLUTION

Product number: PF 93602
Brand: DIESSE

1.2 IDENTIFIED USES OF THE SUBSTANCE OR MIXTURE AND USES ADVISED AGAINST

**Identified use:** Professional use as laboratory reagent.

1.3 DETAILS OF THE SUPPLIER OF THE SAFETY DATA SHEET

Company: DIESSE Diagnostica Senese SpA

<u>Registered office:</u> <u>Production plant:</u>

 Via A. Solari 19
 Via delle Rose, 10 – Loc. Tognazza

 20144 Milan, Italy
 53035 Monteriggioni (SI), Italy

 Tel: +39 02 4859121
 Tel: +39 0577 587111

 Fax: +39 02 48008530
 Fax: +39 0577 318690

1.4 EMERGENCY TELEPHONE NUMBER

**Emergency number:** Centro Antiveleni, Ospedale Niguarda Ca' Granda - Milano

Tel: +39 02 66101029

Centro Antiveleni, Azienda Ospedaliera "S.G.Battista" – Molinette di Torino - Torino

Tel: +39 011 6637637

Centro Antiveleni – U.O. Tossicologia Medica, Azienda Ospedaliera Careggi – Firenze

Tel: +39 055 4277238

Centro Antiveleni, Policlinico A. Gemelli – Università Cattolica del Sacro Cuore - Roma

Tel: +39 06 3054343

Centro Antiveleni, Azienda Ospedaliera A. Cardarelli – Napoli

Tel: +39 081 7472870

web: http://www.salute.gov.it/servizio/documenti/centri antiveleni.pdf

## 2. HAZARDS IDENTIFICATION

## 2.1 CLASSIFICATION OF THE SUBSTANCES OR MIXTURE

Classification according to Regulation (EC) No. 1272/2008:

Skin Corrosion - Category 1A

**Dangerous components:** 

H314 - Causes severe skin burns and eye damage

Classificazione ai sensi delle Direttive 67/548/CEE e 1999/45/CE:

C - Corrosive

R-phrases:

R35 - Causes severe skin burns



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 4/9

#### 2.2 LABEL ELEMENTS

Pictograms: GHS05 Signal word: Danger

## **Hazard statement:**

H314 - Causes severe skin burns and eye damage

## **Precautionary statement:**

Prevention

P260 - Do not breathe dust/fume/gas/mist/vapours/spray

P264 - Wash ... thoroughly after handling

P280 - Wear protective gloves/protective clothing/eye protection/face protection

Reaction

P301 + P330 + P331 - IF SWALLOWED: rinse mouth. Do NOT induce vomiting

P303 + P361 + P353 - IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with

water/shower

P363 - Wash contaminated clothing before reuse

P304 + P340 - IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing

P310 - Immediately call a POISON CENTER or doctor/physician

P305 + P351 + P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing

Storage

P405 - Store locked up.

#### **Contains:**

Sulphuric Acid

INDEX: 016-020-00-8

#### 2.3 OTHER HAZARDS

None

## 3. COMPOSITION/INFORMATION ON INGREDIENTS

3.1 SUBSTANCE Non applicable

3.2 MIXTURE

Classification Conc.

Substance name Directive 67/548/CEE EC Regulation 1272/2008

Sulphuric Acid

Skin Corr. 1A H314 Cas No 7664-93-9 C R35 <3%

EC No 231-639-5

Index No 016-020-00-8

## 4. FIRST AID MEASURES

#### 4.1 DESCRIPTION OF FIRST AID MEASURES

General advice: Move the person from the exposure to open air. In case of severe condition consult a

physician immediately and present this safety data sheet. Eyewash and shower for the

treatment of emergency have to be present in the workplace.

Skin contact: Take off the contaminated clothes and wash with copious amounts of water. If

irritation persist consult a physician.

Eye contact: Take off lenses (if present). Flush eyes with plenty of water for at least 15 minutes

keeping eyelids opened. Consult immediately a physician.

Inhalation: In case of severe condition consult a physician.

Clean mouth with water and drink afterwards plenty of water. Do not induce vomiting. Ingestion:

Do not administer any medicine if the person is unconscious. Consult immediately a

physician.

### 4.2 MOST IMPORTANT SYMPTOMS AND EFFECTS

No data available

#### 4.3 INDICATION OF ANY IMMEDIATE MEDICAL ATTENTION AND SPECIAL TREATMENT NEEDED

No data available



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 5/9

#### 5. FIREFIGHTING MEASURES

#### 5.1 EXTINGUISHING MEDIA

Suitable extinguishing media: Use extinguishing measures that are appropriate to local circumstances and the

surrounding environment.

Unsuitable extinguishing

media: None

#### 5.2 SPECIAL HAZARDS ARISING FROM THE SUBSTANCE OR MIXTURE

The mixture is not inflammable. However, in case of fire, irritating and/or toxic fumes can be released.

#### 5.3 ADVICE FOR FIREFIGHTERS

If necessary, wear self-contained breathing apparatus. Avoid that fire extinguishing water contaminates surface water and/or groundwater.

## 6. ACCIDENTAL RELEASE MEASURES

#### 6.1 Personal precautions, protective equipment and emergency procedures

Keep people not involved in the intervention operations away. Ensure an adequate ventilation of the affected area. Wash with plenty of water all the body parts entered in contact with the mixture immediately. Do not handle damaged containers or the leaked product before wearing appropriate protective outfit.

#### **6.2 ENVIRONMENTAL PRECAUTIONS**

Avoid the contamination of surface water and soil.

#### 6.3 METHODS AND MATERIALS FOR CONTAINMENT AND CLEANING UP

Absorb on inert material (sand/earth). Collect the spilled material with mechanical means and store it in suitable containers for disposal. Use water only to remove residuals, to avoid the danger of spill of product into the sewers.

## **6.4 REFERENCE TO OTHER SECTIONS**

For more information on disposal see section 8 and 13.

### 7. HANDLING AND STORAGE

## 7.1 PRECAUTION FOR SAFE HANDLING

Work in well-ventilated area. Avoid the contact with eyes and skin and the inhalation of vapors.

## 7.2 CONDITION FOR SAFE STORAGE, INCLUDING ANY INCOMPATIBILITIES

Store at 2-30°C, away from basic substances.

#### 7.3 SPECIFIC END USE

None

#### 8. EXPOSURE CONTROL/PERSONAL PROTECTION

#### **8.1 CONTROL PARAMETERS**

Applicable exposure limits

TWA/8h

STEL/15min

Description

Type

Description Type mg/m3 ppm mg/m3 ppm

Sulphuric Acid TLV-ACGIH 2008 0.2

**8.2 EXPOSURE CONTROLS**Work and handle according to the usual precautionary measures for handling chemicals.

Do not eat, drink or smoke while handling the product; wash hands thoroughly with soap

and water before meals and after the work shift.

**Appropriate engineering** 

**controls:** Ensure an adequate ventilation of the working area.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 6/9

#### Personal Protective Equipment:

Suggestions on the use of specific PPE are indicated below. Their choice should be made according to the use of the product and the instructions given by the supplier of the device.

Hand protections:

Chemical-resistant gloves, comply with EN 374.

Eye protections:

Safety goggles side shields

**Body protections:** 

Work outfits for complete skin protection

Respiratory protection:

Not required under normal work activities

8.3 ENVIRONMENTAL EXPOSURE CONTROLS

Do not discard residuals in the environment. Discard the empty containers as hazardous

waste.

## 9. PHYSICAL AND CHEMICAL PROPERTIES

## 9.1 INFORMATION ON BASIC PHYSICAL AND CHEMICAL PROPERTIES

No data available Odour Odour Threshold No data available < 0.5 рΗ Freezing/Melting point No data available Intital boiling and melting point range No data available No data available Flash point Evaporation rate No data available Flammability (solid, gas) Not applicable Upper/lower flammability or explosive limit Not applicable Vapour pressure No data available Vapour density No data available Relative density No data available

Water solubility
Partition coefficient: n-octano/water
Autoignition temperature
Decomposition temperature
Viscosity
No data available
Explosive properties
No data available
Oxidizing properties
No data available
No data available
No data available

#### 9.2 OTHER SAFETY INFORMATION

None

## 10. STABILITY E REACTIVITY

#### 10.1 REACTIVITY

No data available

#### 10.2 CHEMICAL STABILITY

Stable under the recommended transport, handling and storage conditions.

### 10.3 POSSIBILITY OF HAZARDOUS REACTIONS

The mixture reacts with the basic and/or oxidizing substances with formation of heat.

## 10.4 CONDITIONS TO AVOID

Avoid contact with basic and/or oxidizing substances.

#### 10.5 INCOMPATIBLE MATERIALS

Basic substances.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 7/9

#### 10.6 HAZARDOUS DECOMPOSITION PRODUCTS

When heated or in case of fire, vapors that are potentially dangerous to health (SOx) may be produced.

#### 11. TOXICOLOGICAL INFORMATION

#### 11.1 Information on Toxicological Effects

The mixture is corrosive for skin and eyes.

Acute toxicity:

No data available

Skin irritation: Corrosive for skin and mucous membranes

Eye irritation:

Severe eye damages
Subcutaneous sensitization:

No data available
Germ cell mutagenesis:

No data available
Carcinogenicity:

No data available
Reproduction toxicity:

No data available

Other information No data available

## 12. ECOLOGICAL INFORMATION

#### **12.1 TOXICITY**

Use according to good working practices, avoiding release to the environment. In case the product reach waterways or sewers or contaminate soil or vegetation, inform the competent authorities.

## 12.2 PERSISTENCE AND DEGRADABILITY

No data available

#### 12.3 BIOACCUMULATIVE POTENTIAL

No data available

## 12.4 MOBILITY IN SOIL

No data available

#### 12.5 RESULTS OF PBT AND VPVB ASSESSMENT

No data available

## 12.6 OTHER ADVERSE EFFECTS

The disposal into drains or water bodies can result in a decreasing of the pH value. A low pH value harms aquatic organisms.

## 13. DISPOSAL CONSIDERATIONS

## 13.1 WASTE TREATMENT METHODS

Serum samples and reagents used should be treated as infectious waste.

The product and its containers should be considered special waste.

Their transport and disposal should be performed by authorized specialized companies according to the law.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 8/9

### 14. TRANSPORT INFORMATION

Even if the mixture is classified as hazardous by the transport regulations (see table below), the exemptions to the classification of ADR Chapter 3.5 - Dangerous goods packed in excepted quantities – can be applied.

These goods must be transported by vehicles authorized to the transport of hazardous goods according to the requirements of the current edition of ADR and to the national regulations.

People, responsible for the loading and unloading of this dangerous good must be trained on all the risks deriving from the substance and on procedures to be taken in case of emergency.

| ADR/RID:  |   |                  |          |
|---|---|------------------|----------|
| UN N°   | 2796  | Class:           | 8        |
| Classification:                                       | C1  | Packaging group: | II       |
| Tunnel restriction code:                              | (E)   |                  |          |
| If subject to the following characteristic            | s is ADR exempt:                                      |                  |          |
| Exemption for limited quantities: 1L                  |   |                  |          |
| Exemption for expected quantities:                    |   |                  |          |
| <ul> <li>Maximum quantity for the interest</li> </ul> | ernal package: 30 r                                   | nl               |          |
| <ul> <li>Maximum quantity for the ext</li> </ul>      | ernal package: 500                                    | ) ml             |          |
| Sea transport (IMDG Code):                            |   |                  |          |
| UN N°   | 2796  | Class:           | 8        |
| Packaging group:                                      | II  | EmS:             | F-A, S-B |
| Marine polluting agent                                | no  |                  |          |
| Air transport (ICAO-IATA):                            |   |                  |          |
| Class:  | 8   | Packaging group: | II       |
| Name:   | Sulphuric Acid (not containing more than 51% of acid) |                  |          |
| Label:  | n.8   |                  | •        |

## 15. REGULATORY INFORMATION

### 15.1 SAFETY, HEALTH AND ENVIRONMENTAL REGULATIONS/LEGISLATION SPECIFIC FOR THE SUBSTANCE OR MIXTURE

Classification and labeling is done in accordance with the directives 1272/2008 (EC) and its subsequent amendments.

## Legislative Decree 81/2008:

The use of this chemical entails the obligation of "Risk Assessment" by the employer in accordance with the provisions of the Decree April 9, 2008 n. 81 and subsequent amendments. Workers exposed to this chemical agent must not undergo health checks whether the results of the risk assessment show that, depending on the type and quantity of chemical agent and the method and frequency of exposure to this agent, there is only one "Low risk to the safety and irrelevant to the health" of the workers and the measures provided for in the same Decree are sufficient to reduce the risk.

<u>European Standards:1907/2006/EC</u>: Not applicable

<u>Directive 96/82/EC (Seveso Directive)</u>: Not applicable

## 15.2 CHEMICAL SAFETY ASSESSMENT

Not performed for the product.



(According to Regulation (EC) No. 1907/2006 e 453/2010)

Edition: 3

Date: 08-Jul-2014

Page: 9/9

#### 16. OTHER INFORMATION

This product has to be used for diagnostic use only by personnel who is qualified and trained on the hazards shown in this safety sheet.

## Text of hazard statements set out in paragraph 3

H314 Causes severe skin burns and eye damage

#### Text of risk phrases set out in paragraph 3

R35 Causes severe burns.

#### **REFERENCES:**

- 1.Directive 1999/45/EC and following amendments
- 2. Directive 67/548/EEC and following amendments and technical adjustments
- 3. Regulation (EC) 1907/2006 of the European Parliament (REACH) as amended
- 4. Regulation (EC) 1272/2008 of the European Parliament (CLP) as amended
- 5. Commission Regulation (EU) 453/2010
- 6. ECHA European Chemicals Agency www.echa.europa.eu
- 7. The Merck Index.
- 8. Handling Chemical Safety
- 9. NIOSH Registry of Toxic Effects of Chemical Substances
- 10. INRS Fiche Toxicologique

DIESSE Diagnostica Senese SpA declares that the information contained in this data sheet is based on the knowledge available to us on the date of the last version. Users must verify the suitability and accuracy of the provided information according to the specific use of the product. The document has not be accepted as a guarantee of any specific product property. Because of the use of this product is not subject to our direct control, the users, under their own responsibility, must observe the laws and regulations relating to health and safety. We accept no responsibility for improper use.

This data sheet supersedes each previous edition.