

# MAST<sup>®</sup> DISCMASTER 5

Antibiotic Disc Dispenser  
Superior and accurate performance



- Guaranteed in-use disc stability of 4 weeks
- Optimal moisture protection
- Further enhanced to maintain superior performance

# MAST® DISCMASTER 5

The **NEW MAST®** DISCMASTER 5 dispenser has been developed to optimise performance and provides a reliable and robust solution for the dispensing of AST, ID and combination discs within the **MASTDISCS®** range. Due to its novel and innovative design, the **MAST®** DISCMASTER 5 dispenser guarantees ***an in-use shelf life of 4 weeks for all antimicrobials including the carbapenem group*** when stored using charged desiccant.

The **MAST®** DISCMASTER 5 is available in a convenient 6-place format, with an easy to use single-handed operation, permitting smooth integration into routine laboratory workflows and high throughput laboratories.

## Guaranteed in-use stability

All **MASTDISCS®** AST including the carbapenem group are stable for 4 weeks, ensuring accurate and reliable results and convenience for stock control processes.

## Height adjustable

Fully height adjustable skirt, accommodating variable depths of culture media.

## Audible location of cartridge position

Audible 'click' providing reassurance that cartridges are correctly seated in the dispenser.

## Optimal moisture protection

New interlocking seals on the canister, providing improved protection against moisture. Designed to be used with **MAST®** Silica Gel capsules, available from stock SILICA63.

## Positive disc placement

**MASTDISCS®** discs positively placed onto the agar surface permitting precise, centred positioning on each dispensing cycle. Stainless steel base plate and ejector pins provide robust and long lasting performance.

## Lock-off indicated empty cartridge

Automatic lock-off, alerting user to an empty cartridge.



## Ordering Information

Order Code	Product	Pack Size
MDD65	<b>MAST®</b> DISCMASTER 5 Dispenser - 6 place	1
SHD5	Single Cartridge Hand Dispenser	5
SILICA63	Silica Gel Capsule for Discmaster (MDD63 and above)	4
CANISTER65	Canister for MDD65 <b>MAST®</b> DISCMASTER Dispenser	1

ZT170 v3.0 CEA 02/20

**United Kingdom**  
Mast Group Ltd.  
Mast House  
Derby Road, Bootle  
Merseyside L20 1EA

Tel: + 44 (0) 151 933 7277  
Fax: + 44 (0) 151 944 1332  
e-mail: sales@mast-group.com

**Germany**  
Mast Diagnostica GmbH  
Feldstrasse 20  
DE-23858 Reinfeld

Tel: + 49 (0) 4533 2007 0  
Fax: + 49 (0) 4533 2007 68  
e-mail: mast@mast-diagnostica.de

**France**  
Mast Diagnostic  
12 Rue Jean-Jacques Mention  
CS 91106  
80011 Amiens CEDEX 1

Tél. + 33 (0) 322 80 80 67  
Fax + 33 (0) 322 80 99 22  
e-mail: info@mast-diagnostic.fr