


MAST[®] SELECTAVIAL & SELECTATAB

Culture Media Supplements



- Easy to use
- Extensive range
- Convenient presentations
- Complementary range of culture media available

MAST® SELECTAVIAL & SELECTATAB

MAST® SELECTATAB and **MAST® SELECTAVIALS** - culture media supplements to suit all needs.

Mast Group Ltd supplements are available for the selective isolation of a wide range of organisms including Anaerobes, Bordetella, Campylobacters, pathogenic Neisseria, Gardnerella, Yersinia, Legionella, Actinomyces, Staphylococci, Streptococci, Listeria and Clostridium difficile.

Complementary media are available in the **Mast Group Ltd** range of dehydrated media.

MAST® SELECTATAB

A range of lyophilised tablet products containing antibiotics or chemical supplements for direct addition to culture media. Just add one tablet to a specified volume of media, allow the tablet to dissolve, mix well and pour your selective culture plates.

- **Convenient tablet presentation**
- **Easy to use – no weighing or diluting**
- **No antibiotic or chemical required**
- **Strictly quality controlled**

The 200ml and 250ml size **MAST® SELECTATABS** are presented in packs of 25 tablets of one type, 500ml size tablets in packs of 10 tablets. Self-adhesive plate labels are provided in each pack.



MAST® SELECTAVIAL

A range of lyophilised supplements in vials for addition to 1 litre volumes of media following reconstitution with a recommended diluent. Particularly suited to the preparation of large volumes of selective medium.

- **Suitable for large users of selective media**
- **Economical for bulk preparation**
- **Easy to reconstitute**
- **Strictly quality controlled**

MAST® SELECTAVIALS are presented in packs of 10 vials-each vial for addition to 1 litre of medium



ZT222/V1.0/SJW/MAR19

United Kingdom
Mast Group Ltd.
Mast House
Derby Road, Bootle
Merseyside L20 1EA

Tel: + 44 (0) 151 933 7277
Fax: + 44 (0) 151 944 1332
e-mail: sales@mast-group.com

Germany
Mast Diagnostica GmbH
Feldstrasse 20
DE-23858 Reinfeld

Tel: + 49 (0) 4533 2007 0
Fax: + 49 (0) 4533 2007 68
e-mail: mast@mast-diagnostica.de

France
Mast Diagnostic
12 Rue Jean-Jacques Mention
CS 91106
80011 Amiens CEDEX 1

Tél. + 33 (0) 322 80 80 67
Fax + 33 (0) 322 80 99 22
e-mail: info@mast-diagnostic.fr