



# CRYOBANK<sup>®</sup>

Safe and Secure cryopreservation  
of bacteria and fungi

**Full  
traceability**

- Bar coded for full traceability
- Up to 9 years viable storage
- 25 beads per tube for multiple use
- Unique cryoprotectant
- Colour coded caps and beads

# CRYOBANK®

**CRYOBANK®** is based on a cryovial system comprising chemically treated beads covered with a special cryogenic preserving solution. **CRYOBANK®** offers a reliable, convenient and versatile system for storing and preserving fastidious bacteria over long periods. **CRYOBANK®** permits the easy establishment of culture collections for laboratory accreditation of research purposes. Suitable for use in all laboratories dealing with bacteria including clinical, water, brewing, veterinary, food, pharmaceutical, dairy and biotechnology.

## Features and benefits:

- Bar-coded for full traceability suitable for use with all commonly available bar code readers
- Stackable boxes
- Simple to use procedure
- Proven storage of fastidious bacteria, non-fastidious yeasts and moulds over long periods suitable for use down to -80 °C
- Supplied with picker for easy selection of cryovials
- Colour coded beads and caps for easy recognition
- External threaded tube reduce risk of contamination



## Ordering information

Order Code	Product	Content
291800	Cryobank® starter set	1 Cryobank mixed colours, 1 Cryoblock
291801	Cryobank® - red	80 vials
291802	Cryobank® - green	80 vials
291803	Cryobank® - yellow	80 vials
291804	Cryobank® - blue	80 vials
291809	Cryobank® - mixed colours	80 vials
29180DRY	Cryobank® - mixed colours dry	80 vials
291690	Cryobank® storage box (bar-coded)	1 box for 81 vials
291610	Cryoblock	1 block for 18 vials
291611	Substitute polystyrene box for Cryoblock	1 pc

MD V.1.1 / 2025-03-06

**United Kingdom**  
**Mast Group Ltd.**  
**Mast House, Derby Road, Bootle**  
**Merseyside L20 1EA**

**Tel: +44 (0)151 933 7277**  
**Fax: +44 (0)151 944 1332**  
**e-mail: sales@mastgrp.com**

**Germany**  
**Mast Diagnostica GmbH**  
**Feldstraße 20**  
**DE-23858 Reinfeld**

**Tel: +49 (0)4533 2007 0**  
**Fax: +49 (0)4533 2007 68**  
**e-mail: mast@mast-diagnostica.de**

**France**  
**Mast Diagnostic**  
**12 rue Jean-Jacques Mention**  
**CS91106, 80011 Amiens, CEDEX 1**

**Tél: +33 (0)322 80 80 67**  
**Fax: +33 (0)322 80 99 22**  
**e-mail: info@mast-diagnostic.fr**

**www.mastgrp.com**